


# RIDA<sup>®</sup> GENE Color Compensation Kit I

Art. No.: PG0001

3 reactions

 -20 °C



R-Biopharm AG, An der neuen Bergstraße 17, D-64297 Darmstadt, Germany

Tel.: +49 (0) 61 51 81 02-0 / Telefax: +49 (0) 61 51 81 02-20

## 1. Intended use

RIDA<sup>®</sup>GENE Color Compensation Kit I is intended for generating a Color Compensation File for multiplex real-time PCR experiments on the LightCycler<sup>®</sup> 480. The generated Color Compensation File can be applied to analyze multiplex real-time PCR experiments of RIDA<sup>®</sup>GENE real-time PCR kits on the LightCycler<sup>®</sup> 480.

## 2. Explanation of the test

In a multiplex real-time PCR, the wavelengths of light emitted by the reporter dyes may overlap, causing one channel to pick up signals (crosstalk) from a dye measured by another channel. This crosstalk of fluorescence signal can result in incorrect data unless a correction is made by using a Color Compensation File. Color compensation is used to subtract fluorescence crosstalk from a reporter dye into inappropriate channels outside of its dominant emission channel.

## 3. Kit components

Tab.1: Reagents provided (Reagents provided in the kit are sufficient for 3 Color Compensation experiments)

Kit Code	Reagent	Volume	Lid color
1	Blank	1x 400 µl	white
2	Dye 1	1x 400 µl	green
3	Dye 2	1x 400 µl	yellow
4	Dye 3	1x 400 µl	orange
5	Dye 4	1x 400 µl	red

#### **4. Storage instructions**

- Protect RIDA<sup>®</sup>GENE Color Compensation Kit I from light and store at -20 °C.
- RIDA<sup>®</sup>GENE Color Compensation Kit I can be used until the expiration date printed on the label. After expiry the quality guarantee is no longer valid.
- Carefully thaw RIDA<sup>®</sup>GENE Color Compensation Kit I reagents before use (e.g. in a refrigerator at 2 - 8 °C).
- During Color Compensation preparation all the reagents should be stored cold in an appropriate way (2 - 8 °C).

#### **5. Additional equipment and materials**

- LightCycler<sup>®</sup> 480I, LightCycler<sup>®</sup> 480II or LightCycler<sup>®</sup> 480z (Roche)
- Real-time PCR consumables (plate, foil)
- Pipettes (0.5 – 20 µl, 20 – 200 µl, 100 – 1000 µl)
- Filter tips

#### **6. Precautions for users**


- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.

## 7. Protocol for creating a Color Compensation File on the Light Cycler® 480

### 7.1 Preparation of the Color Compensation plate

For a color compensation experiment it is necessary to pipette 5 reactions with 20 µl of each dye and also of the background (blank) into a microwell plate (s. Fig.1).

Fig.1: Pipetting scheme for Color Compensation experiment



Thaw, mix gently and centrifuge briefly the reagents before use. Keep reagents appropriately cold during working step (2 - 8 °C).

Tab.2: Preparation of the Color Compensation


Kit Code	Reagent	Volume per reaction	Pipette 20 µl into the following wells
1	Blank	20 µl	B2, C2, D2, E2, F2
2	Dye 1	20 µl	B4, C4, D4, E4, F4
3	Dye 2	20 µl	B6, C6, D6, E6, F6
4	Dye 3	20 µl	B8, C8, D8, E8, F8
5	Dye 4	20 µl	B10, C10, D10, E10, F10

Cover the microwell plate with an optical sealing foil after pipetting the reagents. Start the real-time PCR according to LightCycler® 480 Set-up.


## 7.2 LightCycler® 480Set-up

**Note:** Login as administrator to set up the detection format in the LightCycler® 480 software.

1. After opening the LightCycler® 480 software, it is necessary to click the “**Tools**” icon to program the detection format.



2. The following window opens. Select “**Detection formats**” in the Tools window. Click the “**New**” button and set-up the detection format (see Tab.3). Save the detection format as “**R-Biopharm**”. Click the “**Close**” button to exit the Tools window.




Tab.3: Detection Channels Set-up for LightCycler® 480


Filter Combination		
LightCycler® 480I	LightCycler® 480II	LightCycler® 480z
450 / 500	465 / 510	465 / 510
523 / 568	533 / 580	540 / 580
558 / 610	533 / 610	540 / 610
615 / 670	618 / 660	610 / 670

**Note:** Set the value for Quant Factor, Melt Factor as well as the Integration Time to 1.


3. After programming the detection format, click the „New Experiment“ button.



4. The following window will open. Select detection format “R-Biopharm” and enter a reaction volume of 20 µl.



5. Program the LightCycler® 480 according to the real-time PCR profile (see Tab.4). Click the “Plus” button to program the four protocol steps.




Tab.4: LightCycler® 480 real-time PCR profile

Program	Cycles / Analysis Mode	Temperature targets			
		Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp rate [°C/s]
Initial Denat.	1 / none	95	none	00:00:30	4.4
Cycling	5 / Quantification	95	none	00:00:15	4.4
		60	single	00:00:30	2.2
TM-Analyse	1 / Color Compensation	95	none	00:00:01	4.4
		50	none	00:00:30	2.2
		70	continuous		0.14 (Acquisitions per °C = 1)
Cooling	1 / none	50	none	00:00:01	2.2


**Note:** Ensure the correct setting of the number of “Analysis Mode“ and of the “Cycles”.

6. After completion of the programming the LightCycler® experiment should look like the picture below.







7. To select the right PCR layout click the „Subset Editor“ button.



8. Click on the “Plus” icon to create a new subset and enter a unique name for the layout (e.g. CC Layout). Press the Strg-key as well as the left mouse button simultaneously and mark all wells of the microwell plate with reagents. Click on the “Apply” button to finish the subset. The final screen should look like the screenshot below.




9. Click the „**Sample Editor**“ button. Select „**Color Comp**“ from the Step 1: “**Select Workflow**” field. Then select the created subset (CC Layout) in the Step 2: **Select Samples** field. To finish the layout select in the „**Dominant Channel**“ dialog field for the reagents ( Dye 1, Dye 2, Dye 3, Dye 4) the corresponding dominant channel (see Tab.5). Select for the samples with the reagent “Blank” the dominant channel “**Water**”. It is not necessary to enter a “Sample Name”.



Tab. 5: Dominant Channel setting for reagents


Reagent	Dominant Channel		
	LightCycler® 480I	LightCycler® 480II	LightCycler® 480z
Blank	Water	Water	Water
Dye 1	450 / 500	465 / 510	465 / 510
Dye 2	523 / 568	533 / 580	540 / 580
Dye 3	558 / 610	533 / 610	540 / 610
Dye 4	615 / 670	618 / 660	610 / 670

10. Place the PCR plate with the prepared reactions into the LightCycler® 480. Click the “**Experiment**” icon and start the experiment by clicking the “**Start Run**” button in the experiment window.




### 7.3 Evaluation and creation of a Color Compensation File

After the experiment ends click the **“Analysis”** button to open the **“Create New Analysis”** dialog box and select **“Color Compensation”**. Select in the appearing dialog box the appropriate subset (e.g. CC Layout) and confirm your selection. Click the button **“Calculate”** in the analysis window to perform color compensation analysis. Click the **“Save CC Object”** button to save the Color Compensation file as **“R-Biopharm”** in the **“CCC”** folder. The created Color Compensation File can now be applied for other LightCycler® 480 experiments.



To apply Color Compensation open the multiplex assay and click the “**Analysis**” button to select the appropriate filter combination. Select “**in Database**” from the Color Compensation drop-down menu and choose the stored Color Compensation file you want to apply to the experiment. Click the “**Color Compensation (Off)**” button, and select the appropriate Color Compensation File. The “**Color Compensation (Off)**” button switches to “**Color Compensation (On)**” to confirm that Color Compensation is active. The multiplex real-time PCR experiment can now be analyzed.



**Note:** The Color Compensation File is specific for every LightCycler® 480 instrument. A new Color Compensation File has to be created if the LightCycler® 480 instrument is changed or after the optical system has been repaired.