

RIDA® GENE EAEC

REF PG2215



1. Intended use

For *in vitro* diagnostic use. RIDA[®]GENE EAEC is a multiplex real-time PCR for the direct, qualitative detection of enteroaggregative *E. coli* (EAEC) in human stool samples.^{1,2} RIDA[®]GENE EAEC multiplex real-time PCR is intended for use as an aid in diagnosis of gastroenteritis caused by enteroaggregative *E. coli*.

2. Summary and explanation of the test

Escherichia coli (*E. coli*) are gram negative, facultatively anaerobic rod bacteria which move by peritrichal flagellation and belong to the Enterobacteriaceae family. *E. coli* are part of the normal intestinal flora of humans and many farm animals and are generally nonpathogenic. Some *E. coli* strains are pathogenic to humans through the acquisition of certain virulence factors (e.g. genes for toxins).

The six known intestinal pathogenic *E. coli*: enterohämorrhagic *E. coli* (EHEC), enteropathogenic *E. coli* (EPEC), enterotoxigenic *E. coli* (ETEC), enteroinvasive *E. coli* (EIEC), enteroaggregative *E. coli* (EAEC) und diffusely adherent *E. coli* (DAEC) can be differentiated by the virulence factors.³

Enteroaggregative *E. coli* (EAEC) were first identified and described in the stool of a child from Chile in 1987.⁴ The defining feature of EAEC is its characteristic aggregative adherence (AA) phenotype. In the gold standard HEp-2 cell adherence assay EAEC adhere to the epithelial cell surface in a “stacked-brick” formation. EAEC are defined as *E. coli* that do not secrete heat-labile or heat-stable enterotoxins and adhere to HEp-2 cells in an AA pattern.⁵ Certain EAEC strains carry a high molecular weight plasmid (pAA) associated with AA, on which a number of virulence genes (e.g. *aggR*, *aggA*, *aafA*, *agg3* and *aatA*) are located.^{6,7,8} Important virulence genes for EAEC detection by PCR are the *aatA* gene (anti-aggregation protein transporter gene, referred to as CVD432 or EAEC probe) and the *aggR* gene (master regulator of the EAEC plasmid virulence genes).^{8,9,10,11} EAEC strains that carry pAA are regarded as typical EAEC while strains that lack pAA are regarded as atypical EAEC.¹² The most common clinical manifestation of EAEC Infection is watery diarrhea. Less common associated clinical symptoms are low-grade fever, nausea, vomiting, abdominal pain and the presence of fecal blood, mucus or leukocytes. The incubation ranges from 8 - 18 h. A voluntary study with an inoculum of 10¹⁰ c.f.u. of EAEC caused diarrheal illness. Due to the high required infective dose suggests a fecal-oral transmission of EAEC by food or water.

EAEC is the cause of acute and chronic (> 14 days) diarrhea among children, adults and HIV-infected persons, in both developing and industrialized countries.⁸

It is the second most common cause of travelers' diarrhea after ETEC among travelers to developing countries, such as Mexico, India and Jamaica.¹³ Outbreaks of EHEC diarrhea have been reported and linked to the consumption of contaminated food.¹⁴ EAEC has been isolated from 2 % - 68 % of patients with diarrhea and from 0 % - 15 % of controls from India, South America, Europe and the Middle East.¹³ In a meta-analysis of published studies, EAEC was a cause of acute diarrhea in a median

of 15 % of children living in developing countries and 4 % of children living in industrialized countries.¹⁵ EAEC were isolated from 2 % of pediatric patients with diarrhea in Germany compared with none of healthy controls.¹² In Swiss children EAEC was isolated in 10.2 % of children with diarrhea compared with 2.2 % of children without diarrhea. In the US EAEC was isolated 4.5 % of case patients versus 1.7 % of controls.¹⁶

3. Test principle

The RIDA[®]GENE EAEC is a multiplex real-time PCR for the direct, qualitative detection of enteroaggregative *E. coli* (EAEC).

After DNA isolation, amplification of the gene fragments specific for EAEC (*aat*, *aggR*, if present) occurs. The amplified targets are detected with hydrolysis probes, which are labeled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the presence of a target the probes hybridize to the amplicons. During the extension step the **Taq-Polymerase** breaks the reporter-quencher proximity. The reporter emits a fluorescent signal which is detected by the optical unit of a real-time PCR instrument. The fluorescence signal increases with the amount of formed amplicons. The RIDA[®]GENE EAEC assay contains an **Internal Control DNA** (ICD) as an internal control of sample preparation procedure and to determine possible PCR-inhibition.

4. Reagents provided

Tab. 1: Reagents provided (Reagents provided in the kit are sufficient for 100 determinations)

Kit Code	Reagent	Amount		Lid Color
1	Reaction Mix	2x	1050 µl	yellow
2	Taq-Polymerase	1x	80 µl	red
D	Internal Control DNA	2x	1700 µl	orange
N	No Template Control	1x	450 µl	white
P	Positive Control	1x	200 µl	blue

5. Storage instructions

- Protect all reagents from light and store at -20 °C. All reagents can be used until the expiration date. After expiry the quality guarantee is no longer valid.
- Carefully thaw reagents before using (e.g. in a refrigerator at 2 - 8 °C).
- Reagents can sustain up to 20 freeze/thaw cycles without influencing the assay performance (e.g. after the first thawing separate it in aliquots and freeze immediately).
- During PCR preparation all the reagents should be stored cold in an appropriate way (2 - 8 °C).

6. Additional necessary reagents and necessary equipment

The RIDA[®] GENE EAEC multiplex real-time PCR assay is suitable for use with following extraction platforms and real-time PCR instruments:

Tab. 2: Necessary equipment

Extraction platforms	
R-Biopharm	RIDA [®] Xtract
Promega	Maxwell [®] RSC
bioMérieux	NucliSENS [®] easyMAG [®]
Real-time PCR instruments	
Roche	LightCycler [®] 2.0, LightCycler [®] 480II
Agilent Technologies	Mx3005P
Applied Biosystems	ABI 7500
Bio-Rad	CFX96 [™]
QIAGEN	Rotor-Gene Q

Note: Only use 0.1 ml tubes on the Rotor-Gene Q (QIAGEN).

If you want to use other extraction platforms or real-time PCR instruments please contact R-Biopharm at mdx@r-biopharm.de.

- RIDA[®] GENE Color Compensation Kit IV (PG0004) for use with the LightCycler[®] 480II
- RIDA[®] GENE Color Compensation Kit II (PG0002) for use with the LightCycler[®] LC2.0
- Real-time PCR consumables (plates, tubes, foil)
- Centrifuge with a rotor for the reaction vials
- Vortexer
- Pipettes (0.5 – 20 µl, 20 – 200 µl, 100 – 1000 µl)
- Filter tips
- Powder-free disposal gloves
- PCR water (BioScience grade, nuclease-free)

7. Precautions for users

For *in-vitro* diagnostic use.

This test must only be carried out by trained laboratory personnel. The guidelines for working in medical laboratories have to be followed. The instruction manual for the test procedure has to be followed. Do not pipet samples or reagents by mouth. Avoid contact with bruised skin or mucosal membranes. During handling reagents or samples, wear appropriate safety clothing (appropriate gloves, lab coat, safety goggles) and wash your hands after finishing the test procedure. Do not smoke, eat or drink in areas where samples or reagents are being used.

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulations for disposal.

For more details see Safety Data Sheets (SDS) at www.r-biopharm.com.

8. Collection and storage of samples

8.1 Sample preparation from stool samples

For DNA isolation of human stool samples, use a commercially available DNA isolation kit (e.g. RIDA[®] Xtract (R-Biopharm)) or DNA extraction system (e.g. Maxwell[®] RSC (Promega)). Extract DNA according to the manufacturer's instructions.

We recommend to dilute the stool samples before extraction 1:3 with water. Vortex the diluted stool sample intensely and centrifuge at 1000 x g for 30 sec. Use from the supernatant the appropriate volume according to the manufacturer's instruction.

The RIDA[®]GENE EAEC assay contains an **Internal Control DNA** that detects PCR inhibition, monitors reagent integrity and confirms that nucleic acid extraction was sufficient. The **Internal Control DNA** can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as PCR inhibition control.

If the **Internal Control DNA** is used only as a PCR inhibition control, 1 µl of the **Internal Control DNA** should be added to the Master-Mix (see Tab.4).

If the **Internal Control DNA** is used as an extraction control for the sample preparation procedure **and** as PCR inhibition control, 20 µl of the **Internal Control DNA** has to be added during extraction procedure. The **Internal Control DNA** should always be added to the specimen-lysis buffer mixture

and must **not** be added directly to the specimen. We also recommend to add 1 µl of the **Internal Control DNA** to the negative control and positive control PCR Mix.

9. Test procedure

9.1 Master-Mix preparation

Calculate the total number of PCR reactions (sample and control reactions) needed. One positive control and one negative control must be included in each assay run.

We recommend to calculate an additional volume of 10 % to compensate imprecise pipetting (see Tab. 3, Tab. 4). Thaw, mix gently and centrifuge briefly the **Reaction Mix**, the **Taq-Polymerase**, the **Positive Control**, the **No Template Control** and the **Internal Control DNA** before using. Keep reagents appropriately cold during working step (2 - 8 °C).

Tab. 3: Calculation and pipetting example for 10 reactions of the Master-Mix (ICD as extraction and PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
	Total	20 µl	220 µl

Mix the components of the Master-Mix gently and briefly spin down.

Tab. 4: Calculation and pipetting example for 10 reactions of the Master-Mix (ICD only as PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
D	Internal Control DNA	1.0 µl	11 µl
	Total	21.0 µl	231.0 µl

Mix the components of the Master-Mix gently and briefly spin down.

9.2 Preparation of the PCR-Mix

Pipette 20 µl of the Master-Mix in each reaction vial (tube or plate).

Negative control: Add 5 µl **No Template Control** to the pre-pipetted Master-Mix.

Note: If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the negative control.

Sample: Add 5 µl DNA extract to the pre-pipetted Master-Mix.

Positive control: Add 5 µl **Positive Control** to the pre-pipetted Master-Mix.

Note: If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the positive control.

Cover tubes or plate. Spin down and place in the real-time PCR instrument. The PCR reaction should be started according to the PCR instrument set-up see (Tab. 5, Tab. 6, Tab.7, Tab. 8).

9.3 PCR Instrument set-up

9.3.1 DNA real-time PCR profile

Tab. 5: DNA real-time PCR profile for LightCycler® series and Rotor-Gene Q

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

Tab. 6: DNA real-time PCR profile for Mx3005P, ABI7500, CFX96™

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

9.3.2 Universal real-time PCR profile

Note: The universal real-time PCR profile should only be used for DNA assays when combining RIDA® GENE DNA and RNA real-time PCR assays in one run.

Tab. 7: Universal real-time PCR profile for LightCycler® series

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

Tab. 8: Universal real-time PCR profile for Mx3005P, ABI7500, CFX96™ and Rotor-Gene Q

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

9.4 Detection channel set-up

Tab. 9: Selection of appropriate detection channels

Real-time PCR instrument	Detection	Detection channel	Note
Roche LightCycler® 2.0	EAEC	530	RIDA® GENE Color Compensation Kit II (PG0002) is required
	ICD	560	
Roche LightCycler® 480II	EAEC	465/510	RIDA® GENE Color Compensation Kit IV (PG0004) is required
	ICD	533/580	
Agilent Techn. Mx3005P	EAEC	FAM	Check that reference dye is none
	ICD	HEX	
ABI 7500	EAEC	FAM	Check that passive reference option ROX is none
	ICD	VIC	
Bio-Rad CFX96™	EAEC	FAM	-
	ICD	VIC	
Qiagen Rotor-Gene Q	EAEC	Green	The gain settings have to be set to 5 according to the default settings
	ICD	Yellow	

10. Quality control

The analysis of the samples is done by the software of the used real-time PCR instrument according to the manufacturer's instructions. Positive control and negative control have to show correct results (see Table 10, Fig. 1) in order to determine a valid run.

The **Positive Control** has a concentration of 10^3 copies/ μ l for stx1, stx2 and eae. In each PCR run it is used in a total amount of 5×10^3 copies.

Tab. 10: For a valid run, the following conditions must be met:

Sample	Assay result	ICD Ct	Target Ct
Positive control	Positive	NA ^{*1}	See Quality Assurance Certificate
Negative control	Negative	Ct > 20	0

**1 No Ct value is required for the ICD to make a positive call for the positive control.*

If the positive control is not positive within the specified Ct range but the negative control is valid, prepare all new reactions including the controls.

If the negative control is not negative but the positive control is valid, prepare all new reactions including the controls.

If the required criteria are not met, following items have to be checked before repeating the test:

- Expiry of the used reagents
- Functionality of the used instrumentation
- Correct performance of the test procedure




Fig.1: Correct run of the positive and negative control (EAEC) on the LightCycler[®] 480II

11. Result interpretation

The result interpretation is done according to Table 11.

Tab.11: Sample interpretation

Virulence factor genes		
aat/aggR	ICD	Result
positive	positive/negative	EAEC detected
negative	positive	Target genes not detected
negative	negative	Invalid

EAEC is detected, if the sample DNA and the Internal Control DNA show an amplification signal in the detection system.

EAEC is also detected, if the sample DNA shows an amplification signal but none for the Internal Control DNA in the detection system. The detection of the internal amplification control is not necessary because high concentrations of the amplicon can cause a weak or absent signal of the Internal Control DNA.

EAEC is not detected, if the sample DNA shows no amplification signal, but an amplification signal for the Internal Control DNA in the detection system. An inhibition of the PCR reaction can be excluded by the detection of the Internal Control DNA.

A sample is invalid, if the sample DNA and Internal Control DNA show no amplification signal in the detection system. The sample contains a PCR inhibitor or a failure occurred in the extraction procedure. The extracted sample needs to be further diluted with PCR water (1:10) and re-amplified, or the isolation and purification of the sample has to be improved.

12. Limitations of the method

1. The result of molecular analysis should not lead to the diagnosis, but always be considered in the context of medical history and symptoms of the patient.
2. This assay is only validated for stool samples.
3. Inappropriate specimen collection, transport, storage and processing or a pathogen load in the specimen below the analytical sensitivity can result in false negative results.
4. The presence of PCR inhibitors may cause invalid results.
5. Mutations or polymorphisms in primer or probe binding regions may affect detection of new variants resulting in a false negative result with the RIDA[®]GENE EAEC assay.
6. As with all PCR based *in vitro* diagnostic tests, extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
7. A positive test result does not necessarily indicate the presence of viable organisms. However, a positive result is indicative for the presence of target genes (aatA, aggR).

13. Performance characteristics

13.1 Analytical sensitivity

The RIDA[®] GENE EAEC multiplex real-time PCR has a detection limit of ≥ 10 DNA copies per reaction.

The following figure2 shows a dilution series of aat/aggR (10^5 - 10^1 DNA copies per μ l) on the LightCycler[®] 480II.




Fig. 2: Dilution series EAEC (10^5 - 10^1 DNA copies/ μ l) on the LightCycler[®] 480II

The detection limit of the whole procedure depends on the sample matrix, DNA extraction and DNA concentration.

13.2 Analytical specificity

The RIDA® GENE EAEC multiplex real-time PCR is specific for enteroaggregative *E. coli*. No cross-reaction could be detected for the following species (see Tab.12):

Tab. 12: Cross-reactivity testing










Adenovirus 1, human, strain Adenoid 71	-	<i>Campylobacter lari</i> subsp. <i>lari</i>	-	<i>Cryptosporidium parvum</i>	-	Norovirus GG II	-
Adenovirus 7, human, strain Gomen	-	<i>Campylobacter upsaliensis</i>	-	<i>E. coli</i> (O157:H7)	-	<i>Proteus vulgaris</i>	-
Adenovirus 40, human, strain Dugan	-	<i>Candida albicans</i>	-	<i>E. coli</i> (O26:H-)	-	<i>Pseudomonas aeruginosa</i>	-
Adenovirus 41, human, strain Tak	-	<i>Citrobacter freundii</i>	-	<i>E. coli</i> (O6)	-	Rotavirus	-
<i>Aeromonas hydrophila</i>	-	<i>Clostridium bifermentans</i>	-	<i>Entamoeba histolytica</i>	-	<i>Salmonella enteritidis</i>	-
<i>Arcobacter butzleri</i>	-	<i>Clostridium difficile</i>	-	<i>Enterobacter cloacae</i>	-	<i>Salmonella typhimurium</i>	-
Astrovirus	-	<i>Clostridium novyi</i>	-	<i>Enterococcus faecalis</i>	-	<i>Serratia liquefaciens</i>	-
<i>Bacillus cereus</i>	-	<i>Clostridium perfringens</i>	-	<i>Giardia intestinalis</i> Portland 1	-	<i>Shigella flexneri</i>	-
<i>Bacteroides fragilis</i>	-	<i>Clostridium septicum</i>	-	<i>Giardia intestinalis</i> WB Clone C6	-	<i>Staphylococcus aureus</i>	-
<i>Campylobacter coli</i>	-	<i>Clostridium sordellii</i>	-	<i>Giardia lamblia</i>	-	<i>Staphylococcus epidermidis</i>	-
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	-	<i>Clostridium sporogenes</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Vibrio parahaemolyticus</i>	-
<i>Campylobacter jejuni</i>	-	<i>Cryptosporidium muris</i>	-	Norovirus GG I	-	<i>Yersinia enterocolitica</i>	-

14. Version history

Version number	Chapter and designation
2013-01-30	Release version
2018-09-03	General revision
2018-09-03	4. Reagents provided 5. Storage instructions 6. Additional necessary reagents and necessary equipment 9. Test procedure 10. Quality control 11. Result interpretation 13. Performance characteristics 14. Version history 15. Explanation of symbols

15. Explanation of symbols

General symbols

	For in vitro diagnostic use
	Consult instructions for use
	Lot number
	Expiry
	Store at
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

Testspecific symbols

16. Literature

1. Müller D *et al.* Identification of unconventional intestinal pathogenic *Escherichia coli* isolates expressing intermediate virulence factor profiles by using a novel single-step multiplex PCR. *Appl Environ Microbiol* 2007, 73 (10): 3380-3390.
2. Cordeira F *et al.* Evaluation of a Multiplex PCR for Identification of Enteroaggregative *Escherichia coli*. *J Clin Microbiol* 2008, 46 (2): 828-829.
3. Kaper JM *et al.* PATHOGENIC *ESCHERICHIA COLI*. *Nature Reviews Microbiology* 2004, 2:123-140.
4. Nataro JP. Enteroaggregative *Escherichia coli* pathogenesis. *Curr Opin Gastroenterol* 2005, 21: 4-8.
5. Nataro JP and Kaper JB. Diarrheagenic *Escherichia coli*. *Clinical Microbiology Reviews* 1998, 11(1): 142-201.
6. Law D and Chart H. Enteroaggregative *Escherichia coli*. *J Appl Microbiol* 1998, 84: 685-697.
7. Vial PA *et al.* Characterization of enteroadherent-aggregative *Escherichia coli*, a putative agent of diarrheal disease. *J Infect Dis* 1988, 158: 70-79.
8. Huang DB *et al.* A review of an emerging enteric pathogen: Enteroaggregative *Escherichia coli*. *J Med Microbiol* 2006, 55: 1303-1311.
9. Baudry B *et al.* A sensitive and specific DNA probe to identify enteroaggregative *Escherichia coli*, a recently discovered diarrheal pathogen. *J Infect Dis* 1990, 161: 1249-1251.
10. Schmidt H *et al.* Development of PCR for screening of enteroaggregative *Escherichia coli*. *J Clin Microbiol* 1995, 33: 701-705.
11. Nishi J *et al.* The Export of Coat Protein from Enteroaggregative *Escherichia coli* by a Specific ATP-binding Cassette Transporter System. *J Biol Chem* 2003, 278 (46): 45680-45689.
12. Weintraub A. Enteroaggregative *Escherichia coli*: epidemiology, virulence and detection. *J Med Microbiol* 2007, 56: 4-8.
13. Adachi JA *et al.* Enteroaggregative *Escherichia coli* as a major etiologic agent in traveler's diarrhea in 3 regions of the world. *Clin Infect Dis* 2001, 32: 1706-1709.
14. Nataro JP *et al.* Enteroaggregative *Escherichia coli*. *Emerg Infect Dis* 1998, 4(2): 251-261.
15. Huang DB *et al.* Enteroaggregative *Escherichia coli* is a cause of acute diarrheal illness: a meta-analysis. *Clin Infect Dis* 2006, 43: 556-563.
16. Cennimo DJ *et al.* Enteroaggregative *Escherichia coli*: A Review of Trends, Diagnosis, and Treatment. *Infect Med* 2007, 24: 100- 110.