# Instructions for use / Testanleitung

## Pancreatic Elastase ELISA

ELISA for the determination of human pancreatic elastase  
ELISA zur Bestimmung der humanen pankreatischen Elastase

R-Biopharm Order Code No.: G09040 Version SK15

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---

**Intended Use**

The Pancreatic Elastase ELISA from BIOSERV Diagnostics is an enzyme-linked immunosorbent assay for the quantitative determination of human pancreatic elastase in faeces.

---

**Clinical Significance**

Pancreatic Elastase is a proteolytic enzyme exclusively produced in pancreas. The enzyme stability is remarkably high during the intestinal passage and is even accumulated in a six-fold concentration in stool, compared to the concentration in the duodenal juice. The determination of enzyme concentration in faeces reflects the exocrine secretory capacity of the pancreas.

The polyclonal antibodies used in this assay are specifically directed against defined peptide sequences of the human pancreatic elastase molecule. These sequences are species- and organ specific.

---

**Fields of Application**

The Pancreatic Elastase ELISA from BIOSERV Diagnostics is a part in clinical practice to diagnose exocrine pancreatic insufficiency, which may be caused by:

- Chronic pancreatitis
- Cystic fibrosis
- Diabetes mellitus
- Cholelithiasis
- Hereditary Pancreatitis
- Chronical- inflammatory bowel disease
- Pancreatic carcinoma
- Autoimmunologically caused pancreatitis
- Shwachman-Diamond-Syndrome
- Autoimmune- and connective tissue diseases
- Zollinger-Ellison-Syndrome

---

**Principles of the Assay Method**

The Pancreatic Elastase ELISA (**E**nzyme **L**inked **I**mmuno**S**orbent **A**ssay) from BIOSERV Diagnostics is a solid-phase enzyme-linked immunosorbent assay based on a double-sandwich technique for the quantitative determination of human pancreatic elastase in faeces. Two polyclonal antibodies are recognizing several different epitopes on defined species- and organ-specific human pancreatic elastase peptide sequences. BIOSERV Diagnostics is using polyclonal antibodies which recognize several different epitopes in parallel to reach a higher diagnostic sensitivity and specificity (see bibliography).

The ELISA microplate is coated with antibodies directed against human pancreatic elastase binding the pancreatic elastase contained in the patient samples or in the standards, respectively. In the following step the second antibody, labelled with biotin, binds to the immobilized pancreatic elastase. To visualize the bound pancreatic elastase, the biotin binds in the following step to streptavidin-labelled horseradish-peroxidase. The peroxidase then oxidizes the substrate TMB (3,3',5,5'-tetramethylbenzidine). The reaction will be stopped by addition of 0.25 mol/l H<sub>2</sub>SO<sub>4</sub>. The developed dye (oxidized TMB) can be measured photometrically at 450 nm. A reference measurement using a reference filter at ≥ 550 nm wavelength is recommended but not mandatory.

---

**Reagents** (sufficient for 96 determinations)

1. Reference standard set (ready for use) - per vial	0.7 ml
- Standard SK15 ( 15 µg elastase/g – colourless screw cap with pink insert)	
- Standard 1 ( 50 µg elastase/g – colourless screw cap)	
- Standard 2 ( 100 µg elastase/g – white screw cap)	
- Standard 3 ( 200 µg elastase/g – yellow screw cap)	
- Standard 4 ( 500 µg elastase/g – blue screw cap)	
2. Control 1 (equivalent to 70 - 110 µg elastase/g – brown screw cap)	0.7 ml
3. Control 2 (equivalent to 160 - 240 µg elastase/g – green screw cap)	0.7 ml
4. Washing solution (10x)	2 x 100 ml
5. Elastase-antibody biotinylated (201x, red screw cap)	0.12 ml
6. Streptavidin conjugate (ready for use)	8 ml
7. Substrate buffer solution (solution of TMB, ready for use)	13 ml
8. Stop solution (0.25 mol/l H <sub>2</sub> SO <sub>4</sub> , ready for use)	13 ml
9. Microtiter strips + holder	96 wells
10. Dilution bottle for elastase-antibody biotinylated	1 x

---

**Materials Required but not Included**

1. Microplate reader with 450 nm filter, optionally with a reference filter  $\geq 550$  nm.
2. Microliter pipettes with disposable tips: 5 µl, 50 µl, 100 µl and 1000 µl.
3. Tubes for the dilution of the samples.
4. Ultrapure water (recommendation: according to ISO 3696 grade 1).
5. Absorbent paper.
6. Accuracy weighing scales for weighing the stool samples or BIOSERV Stool Sample Preparation Set (Cat.-No. GZ3008).
7. Incubator respectively incubation chamber for automates.
8. Extraction buffer (Cat.-No. G09038P01) or BIOSERV Stool Sample Preparation Set (Cat.-No. GZ3008)

---

**Warnings and Precautions**

1. This kit is intended for *in vitro* use only.
2. Avoid contact with the stop solution, it may cause skin irritations and burns.
3. Do not pipette reagents by mouth.
4. Please use only calibrated pipettes and instruments.
5. Please regard all samples as potentially infectious and handle them with utmost care.
6. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation where this exists.

---

**Instructions for Reagent Preparation**

1. The components of this kit are intended for use as an integral unit and should not be interchanged with the components of other kits.
2. All reagents and specimens have to be brought to room temperature (18 – 25°C / 64 – 77 °F) before use.
3. All reagents have to be mixed without foaming.
4. Once the test procedure has been started, all steps should be continued without interruption.
5. Pipette all reagents and samples onto the bottom of the wells. Mixing or shaking after pipetting is not required.
6. Use new disposable tips for each specimen.
7. Before starting the assay, all reagents to be used should be prepared and ready for immediate use, all needed strips should be secured in the holder etc. This will ensure equal time periods for each pipetting step without interruption.
8. For optimal results it is important to wash the wells thoroughly after incubation and to remove even the last water drops by hitting the plate on absorbent paper or cloth.
9. It is recommended to conduct all tests in double determination in order to minimize the

consequences of pipetting or handling errors.

---

### Storage Instructions and Shelf Life Information

1. Store the reagents at 2 – 8 °C (36 – 46 °F). **Do not freeze the kit or kit components!**
2. The reagents which are not yet opened remain stable until the expiration date of the kit.  
Attention please: Once the kit is in use, the remaining shelf life is much shorter than mentioned at the kit label. All ingredients have to be closed accurately immediately after use, stored at 2 – 8 °C temperature and should then be used up **within the next 4 weeks**. Otherwise the results could be differing.
3. The diluted washing solution is stable for 4 weeks at refrigerator temperatures (2 – 8 °C / 36 – 46 °F).
4. The diluted extraction buffer can be stored for 4 weeks at refrigerator temperatures (2 – 8 °C / 36 – 46 °F, **not included!**).
5. Put caps back on the vials immediately after use.
6. Store the microtiter strips in a dry bag with desiccants. The remaining strips must be stored in the tightly resealed bag together with the desiccants. Under these storage conditions, they are stable at least for 4 weeks after opening of the sealed bag.

---

### Sample Material

Human Faeces

---

### Specimen Collection and Preparation

1. Use fresh faeces samples. Handle all samples with utmost care and treat them as potentially infectious material.
2. Samples may be stored at different temperatures for the following time-spans:  
Environmental temperature up to 40 °C (104 °F): up to five days  
Refrigerator temperature 2 – 8 °C (36 – 46 °F): up to one week  
Household freezer temperature ca. -18 °C (-0.4 °F): up to one year
3. Procedure for watery stool:  
In case of a watery stool sample, try to collect another sample with a more solid consistency. If this is not possible, the sample can be pipetted. Since the dilution of the watery stool sample may cause false pathological (false low) results, the consistency must be noted. In case of a pathological result the test needs to be repeated with a formed stool sample.

#### Precautions for users:

The test kit does not contain any substances of human origin. However, the stool samples must be treated as potentially infectious and handled in accordance with national safety regulations. Please adhere to the applicable regulations for accident prevention, environment protection and handling of hazardous substances. Please wear suitable protective clothing and disposable gloves in any case.

#### Preparation of the faeces samples

1. Preparation of the Extraction buffer (10 x concentrated, **not included!**):
  - Dilute the extraction buffer 1:10 with ultrapure water (e.g. 50 ml + 450 ml). The diluted solution can be stored for 4 weeks at refrigerator temperatures (2 – 8 °C / 36 – 46 °F).
  - Use a 12 ml tube with a spatula (or with an inoculation loop) to weigh 30 – 100 mg faeces on a scales with a sensitivity of 1 mg. Add 1 ml diluted extraction buffer per 10 mg faeces (e.g. 70 mg faeces in 7 ml buffer, 83 mg faeces in 8.3 ml buffer).
  - Homogenise the samples thoroughly by means of a vortex mixer (2 minutes). After sedimentation of solid constituents (at least 15 – 30 minutes), the supernatant can be used for the determination of the Elastase after dilution. It is also possible to let the sedimentation take place overnight at 2 – 8 °C (36 – 46 °F).
2. Usage of the BIOSERV Stool Sample Preparation Set.

- Optionally to manual weighing of faeces samples BIOSERV Diagnostics offers the completely ready-to-use Stool Sample Preparation Set (Cat.-No. GZ3008) for convenient, reliable and hygienic preparation of faeces samples. Each device is already filled with extraction buffer.

**Attention:** If in extraction buffer or washing solution salt crystals have been falling out, please dissolve them by shaking carefully and by growing warm in a water bath at 30 – 35 °C before use.

For later examinations the supernatants can be frozen at -18 °C for a time period of four weeks. **Attention:** It is absolutely necessary to separate clearly the supernatants from the sediment before storage.

---

## Assay Procedure

1. Preparation of the washing solution (10 x concentrated): The concentrated washing solution (100 ml) must be diluted with 900 ml ultrapure water (according to ISO 3696 grade 1). **Attention:** Use ultrapure water only for preparation of washing buffer.
2. Warm all reagents to room temperature and mix thoroughly before use.
3. Fix the required number of coated wells or strips in the strip holder.
4. Dilute the supernatants of extracted faeces 1:201 (1+200) with washing solution (e.g. 25 µl supernatant + 5 ml washing solution). The dilutions are now ready for use. Please use a new disposable tip in each case.
5. Dispense 50 µl of the blank (= zero standard, = washing solution), elastase standards, controls, each extracted and diluted faeces samples with new disposable tips into the respective wells.
6. Incubate for 60 minutes at 37 °C / 98 °F.
7. Briskly shake out the contents of the strips and then rinse the wells 3 times with 200 µl washing buffer each.
8. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
9. Dilute the biotinylated anti-elastase antibody 1:201 (red screw cap) with washing solution. Use the included empty dilution bottle. **Attention: The dilution bottle is intended for one-time use only!** (1 part biotinylated antibody + 200 parts washing solution). Examples:
  - 5 µl biotinylated anti-elastase antibody in 1 ml washing solution (for 2 strips)
  - 10 µl biotinylated anti-elastase antibody in 2 ml washing solution (for 4 strips)
  - 15 µl biotinylated anti-elastase antibody in 3 ml washing solution (for 6 strips)
  - 25 µl biotinylated anti-elastase antibody in 5 ml washing solution (for 10 strips)Dispense 50 µl into each well.  
**Attention:** The biotinylated antibody solution has to be prepared freshly before every run of the test.
10. Incubate for 30 minutes at 37 °C / 98 °F.
11. Briskly shake out the contents of the strips and then rinse the wells 3 times with washing solution, each time using 200 µl.
12. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
13. Dispense 50 µl of streptavidin peroxidase conjugate (ready for use) into each well.
14. Incubate for 30 minutes at 37 °C / 98 °F.
15. Briskly shake out the contents of the strips and then rinse the wells 3 times with washing solution, each time using 200 µl.
16. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
17. Dispense 100 µl of substrate solution into each well.
18. Incubate for 20 minutes at 37 °C / 98 °F. **Attention:** Set time at start of pipetting the first substrate sample. Incubate in the dark.

19. Stop the enzymatic reaction by adding 100 µl stop solution to each well. **Attention:** Apply the stop solution then in the same chronological order and using the same time intervals as when adding the substrate.

Read the absorbance of each well at 450 nm with a microplate reader. It is recommended that the wells be read within 10 minutes after stopping the enzymatic reaction. Any 96 well microplate reader capable of determining the absorbance at 450 nm may be used. A reference measurement using a wavelength  $\geq 550$  nm is recommended, but not absolutely necessary.

### Pipetting Scheme for the Pancreatic Elastase ELISA SK15 from BIOSERV Diagnostics

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK	BLANK	P	1	P	9	P	17	P	25	P	33
B	SK	15	P	2	P	10	P	18	P	26	P	34
C	S	1	P	3	P	11	P	19	P	27	P	35
D	S	2	P	4	P	12	P	20	P	28	P	36
E	S	3	P	5	P	13	P	21	P	29	P	37
F	S	4	P	6	P	14	P	22	P	30	P	38
G	C	1	P	7	P	15	P	23	P	31	P	39
H	C	2	P	8	P	16	P	24	P	32	P	40

In this pipetting scheme the recommended positions for the blank (zero standard), standards (SK15 – S4), controls (C1 and C2) and for the patient samples (P1 – P40) are shown as double determinations.

### Calculation of the Results

1. The manufacturer recommends for test analysis the linear regression respectively a four parameter fit.
2. Calculate the average absorbance values for each set of reference standards, controls and patient samples.
3. The average absorbance (Y) of each reference standard is plotted against its corresponding concentration (X).
4. Use the average absorbance of each patient sample to determine the corresponding elastase value by simple interpolation from this standard curve.
5. The value of control 1 should be within a range of 70 µg/g and 110 µg/g and control 2 should be within a range of 160 µg/g and 240 µg/g. If either control is outside its corresponding range the test should be repeated. In this case please check all incubation steps and times.

### Limitations of the Assay

- At temperatures higher than 40 °C (104 °F) the samples should be transported cooled or refrigerated.
- Watery faeces from patients with diarrhoea may lead to falsely low readings because of a dilution effect. For further advice regarding treatment of watery stool, please refer to point 3 under paragraph "Specimen Collection and Preparation".
- If the optical density of the blanks (= zero standard, = washing solution) is higher than OD 0.15 the test should be repeated.
- It could be possible that a cross-reaction with porcine proteins may occur. Therefore it would be recommended to interrupt an enzyme substitution therapy. But there is no proven knowledge for what concerns the necessity to interrupt a substitution therapy before conducting this assay.

## Expected Values

- |  |                                |
|--|--------------------------------|
| - severe exocrine pancreatic insufficiency   | < 100 µg elastase/g faeces     |
| - moderate exocrine pancreatic insufficiency | 100 – 200 µg elastase/g faeces |
| - normal exocrine pancreatic function        | > 200 µg elastase/g faeces     |

## Assay Performance Characteristics

### 1. Diagnostic specificity: 95 %

Samples from 609 healthy individuals were investigated. 78 of these subjects were healthy blood donors, for the remaining 531 patients a pancreatic disease had been excluded by other diagnostic methods (ultrasonography, ERCP).

### 2. Diagnostic sensitivity:

Severe chronic pancreatitis (below 100 µg pancreatic elastase per g stool) **94 %**

Samples from 46 patients suffering from a severe chronic pancreatitis, diagnosed by ultrasonography and ERCP, were investigated.

Mild chronic pancreatitis (between 100 and 200 µg pancreatic elastase per g stool) **63 %**

Samples from 56 patients suffering from a mild chronic pancreatitis, as diagnosed by ultrasonography and ERCP, were investigated.

### 3. Diagnostic sensitivity for pancreatic carcinoma: 61 %

Samples from 51 patients suffering from clinically diagnosed pancreatic cancer were investigated.

### 4. Diagnostic sensitivity for cystic fibrosis: 100 %

Samples from 36 patients suffering from a clinically diagnosed cystic fibrosis were investigated.

### 5. Intra-Assay Reproducibility

The intra-assay reproducibility was determined by measuring six clinically defined faeces samples (covering the whole measuring range of the kit) in 36 replicates each. The testing was conducted in three different lots. (N = 36 x 3)

	Elastase concentration (mean value in µg/g)	Standard deviation (SD in µg/g)	Coefficient of variation (CV in %)
Sample 1	45	3,80	5,96 (3,88 – 7,10)
Sample 2	105	8,73	4,93 (3,82 – 6,03)
Sample 3	250	15,11	3,51 (2,43 – 4,85)
Sample 4	182	12,2	3,32 (2,63 – 4,07)
Sample 5	327	12,5	2,77 (2,45 – 3,14)
Sample 6	385	23,4	2,92 (2,21 – 4,01)

## 6. Inter-Assay Reproducibility

The inter-assay reproducibility was determined by measuring four clinically defined faeces samples (covering the whole measuring range of the kit) on ten individual days. The measurements were conducted in two separate runs each and as double replicate. (N = 10 x 2 x 2)

	Elastase concentration (mean value in µg/g)	Standard deviation (SD in µg/g)	Coefficient of variation (CV in %)
Sample 1	55	5,17	10,21
Sample 2	264	16,67	6,62
Sample 3	193	8,79	8,26
Sample 4	328	17,66	5,62

## 7. Inter-Lot Reproducibility

The inter-lot reproducibility was determined by measuring six clinically defined faeces samples (covering the whole measuring range of the kit) in three different lots (N = 3). The measurements were conducted in two separate runs each as double replicate.

	Elastase concentration (mean value in µg/g)	Standard deviation (SD in µg/g)	Coefficient of variation (CV in %)
Sample 1	92	3,85	4,17
Sample 2	248	8,03	3,24
Sample 3	102	5,30	5,19
Sample 4	187	7,68	4,10
Sample 5	316	14,32	4,53
Sample 6	396	15,56	3,93

---

## **Bibliography Regarding the Pancreatic Elastase ELISA from BIOSERV Diagnostics**

Weiss FU, Budde C, Lerch MM (2016): Specificity of a Polyclonal Fecal Elastase ELISA for CELA3. *PLoS One*. 2016 Jul 26;11(7):e0159363.

Erickson JA, Aldeen WE, Grenache DG, Ashwood ER (2008): Evaluation of a fecal pancreatic elastase-1 enzyme-linked immunosorbent assay: Assessment versus an established assay and implication in classifying pancreatic function. *Clinica Chimica Acta*, Vol 397, issues 1-2, November 2008, pages 87-91.

Qualia CM, Villalona JF, Ren C, Rossi TM (2007): Specificity of the Polyclonal Antibody Test System For Human Elastase in Stool. *Journal of Pediatric Gastroenterology and Nutrition*, Vol 45:E26 # 66, 2007

Weiss U, Ruthenbueger M, Hammer E, Voelker U, Lerch MM (2006): Assessment of Isoform Specificity of a Polyclonal Elastase ELISA. *Journal of Pediatric Gastroenterology and Nutrition* 43:E32 # 58, 2006.

Abdel Rahman A, Abdul Wahab A, Abdel Rahman MO (2006): Faecal elastase-1 concentration in cystic fibrosis patients with CFTR I1234V mutation. *Acta Paediatrica*, 95: 1066-1069.

Hahn JU, Bochnig S, Kerner W, Koenig H, Sporleder B, Lankisch PG, Maisonneuve P, Lowenfels AB (2005): A New Fecal Elastase 1 Test Using Polyclonal Antibodies for the Detection of Exocrine Pancreatic Insufficiency. *Pancreas* Vol 30, No 2, pages 189 – 190.

Miendje Y, Maisin D, Sipewa MJ, Deprez P, Buts JP, De Nayer P, and Philippe M (2004): Polyclonal versus monoclonal ELISA for the determination of faecal elastase-1: Diagnostic Value in Cystic Fibrosis and Chronic Pancreatic Insufficiency. *Clinical Laboratory*, Vol. 50, Nr. 7+8, pages 419 – 414

Keim V, Teich N, and Moessner J (2003): Clinical Value of a New Fecal Elastase Test for Detection of Chronic Pancreatitis. *Clinical Laboratory* Vol. 5 +6, pages 209 – 215.

Garcia-Bueno CA, Rossi TM, Lee KW, Yuwono, MT, Robinson A, Tjota A (2002): Quantification of fecal elastase-1 using either polyclonal or monoclonal antibodies. *Gastroenterology* Vol. 122 (4), page A-510.

Deprez PH, Natale M Del, Deji MYV, Pauwels S, Philippe M (2002): Comparative Evaluation of 13C-Mixed Triglyceride Breath Test and Faecal Elastase 1 Tests in the Assessment of Exocrine Pancreatic Insufficiency. *Gastroenterology*, Vol 122, no 4, page A-510.

Tomášová H, Zemková D, Bartosová J, Skalická, V, Koloušková S, Nevoral J, Macek MJ and Vávrová V (2002): Proper interpretation of Elastase 1 activities in Stool: experience of the Prague CF centre. *Proceedings of the 25th European Cystic Fibrosis Conference (Volume ISBN 88-323-2622-1, CD ISBN 88-323-2623-X)*; page 111 – 114

Keim V, Teich N, and Moessner J. (2000): Value of polyclonal elastase ELISA for diagnosis of chronic pancreatitis. *Pancreas*, Vol. 21, Number 4, November 2000, page 453.

**Also available from BIOSERV Diagnostics:**

Please contact us for more information: [info@bioserv-diagnostics.com](mailto:info@bioserv-diagnostics.com)

<b>Gastroenterology</b>	<b>Order Code</b>	<b>Determinations per Kit</b>
Human Pancreatic Elastase ELISA (pancreatic insufficiency)	BS - 86 - 01	96
Accessories Kit: Stool Sample Preparation Kit with 45 Tubes	BS - 00 - 03	
Pankrin™ ELISA (acute pancreatitis)	BS - 86 - 02	96
<b>Infertility</b>		
Anti-Spermatozoa Antibody Latex Agglutination Test	BS - 10 - 10	50
Anti-Spermatozoa Antibody ELISA for Serum	BS - 10 - 20	96
Anti-Spermatozoa Antibody ELISA for Seminal Plasma	BS - 10 - 21	96
Anti-Spermatozoa Antibody Haemagglutination Test	BS - 10 - 30	40
Anti-Spermatozoa Antibody ELISA - Ig Classifying for Serum	BS - 10 - 50	96
Anti-Zona Pellucida Antibody Latex Agglutination Test	BS - 20 - 10	50
Anti-Zona Pellucida Antibody ELISA	BS - 20 - 20	96
Anti-Zona Pellucida Antibody Haemagglutination Test	BS - 20 - 30	40
Anti-Zona Pellucida Antibody ELISA - Ig Classifying	BS - 20 - 50	96
Anti-Ovary Antibody Latex Agglutination Test	BS - 40 - 10	50
Anti-Ovary Antibody ELISA	BS - 40 - 20	96
Anti-Ovary Antibody Haemagglutination Test	BS - 40 - 30	40
Anti-Ovary Antibody ELISA - Ig Classifying	BS - 40 - 50	96
<b>Monitoring of Risk Pregnancies</b>		
IGF-BP1 ELISA (PP12)	BS - 30 - 10	96
Glycodelin ELISA (PP14)	BS - 30 - 20	96

**BIOSERV Diagnostics GmbH**

Doberaner Str. 151  
18057 ROSTOCK / GERMANY

Phone: +49 (0) 381 3758 2091

Fax: +49 (0) 381 3758 2099

E-mail: [info@bioserv-diagnostics.com](mailto:info@bioserv-diagnostics.com)

Web: [www.bioserv-diagnostics.com](http://www.bioserv-diagnostics.com)