# Pancreatic Elastase ELISA

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# Pancreatic Elastase ELISA

Enzyme immunoassay for the quantitative detection of pancreatic elastase in stool samples of human origin

**REF** G09038  96  
**IVD** In-vitro diagnostic medical device **CE**

 **Seramun Diagnostica GmbH** • Sprenhagener Str. 1 • 15754 Heidesee • Germany •  
T +49 33767 791-10 • [info@seramun.com](mailto:info@seramun.com) • [www.seramun.com](http://www.seramun.com)

 <b>IVD</b> In-vitro diagnostic medical device	 <b>UDI</b> Unique device identifier	 Manufacturer
 Country of manufacture and date of manufacture	 <b>REF</b> Article number	 <b>SN</b> Serial number
 Keep away from sunlight	 Humidity limitation	 <b>LOT</b> Batch code
 Consult instructions for use	 Temperature limit	 Do not reuse
 Sufficient for <i>n</i> tests	 Biohazard	 Use-by date
		 Attention

## Intended Use

Pancreatic Elastase ELISA is an IVD test for the quantitative determination of pancreatic elastase in stool samples of human origin by a laboratory professional user.

The test can be used in combination with accessories Stool Preparation Kit from Bioserv Diagnostics GmbH (GZ3008) or Stool Preparation Kit from Seramun Diagnostica GmbH (ACS-001). Alternatively, it can be used with Extraction Buffer (10x) from Bioserv Diagnostics GmbH (G09038P01) or Extraction buffer (10x) from Seramun Diagnostica GmbH (ACS-002).

It is intended to aid in the diagnosis of exocrine pancreatic insufficiency in samples of patients with symptoms of an exocrine pancreatic insufficiency.

The test must not be used with specimen materials other than stool samples of human origin, for screening, monitoring, diagnosis, prediction, prognosis, as companion diagnostic, in the near patient setting and by lay persons.

## Principle of the Test

Pancreatic Elastase ELISA (Enzyme Linked ImmunoSorbent Assay) is a solid phase enzyme immunoassay based on the double sandwich technique for the quantitative determination of pancreatic elastase in human stool samples.

The wells of the microtiter plate are coated with polyclonal antibodies, which immobilize human pancreatic elastase molecules from patient samples by antibody-antigen binding in the first incubation step. In the second incubation step, biotin-labelled polyclonal antibodies are added that bind to the immobilized human pancreatic elastase. In the third incubation step, the streptavidin-peroxidase (HRP) conjugate binds to the biotin. To remove unbound components, a washing step is performed after each incubation step.

In the following enzymatic reaction step, the peroxidase of the streptavidin-peroxidase conjugate converts the colorless substrate solution to a blue reaction product. The reaction is stopped by the addition of stop solution, resulting in a color change from blue to yellow. The optical density (OD) of the reaction product measured at 450 nm and 620 nm reference filter, respectively, is directly proportional to the concentration of the specifically bound pancreatic elastase.

## Test Components (Delivery Scope)

			<b>For 96 wells</b>
1	<b>WELLS</b>	<b>Microtiter plate</b> coated with polyclonal anti-pancreatic-elastase-antibodies (rabbit)	12 single breakable 8-well strips, vacuum-sealed with desiccant
2	<b>WASHBUF (10x)</b>	<b>Wash buffer (10x)</b> Phosphate-based buffer	2 x 100 mL concentrate, for 2 x 1000 mL solution, colorless white cap
3	<b>STD</b>   1 - 4	<b>Standards 1 – 4</b> STD 1                      50 µg/g STD 2                      100 µg/g STD 3                      200 µg/g STD 4                      500 µg/g	0.7 mL each, ready to use, colorless cap white cap yellow cap blue cap
4	<b>CONTROL</b>   1	<b>Control 1</b> 70 - 110 µg/g	0.7 mL, ready to use, colorless brown cap
5	<b>CONTROL</b>   2	<b>Control 2</b> 160 - 240 µg/g	0.7 mL, ready to use, colorless, green cap
6	<b>CONJ BIOTIN (201x)</b>	<b>Elastase-antibody biotinylated (201x)</b> biotinylated, polyclonal anti-pancreatic-elastase-antibody (rabbit)	0.12 mL, concentrate, colorless for 24 mL solution red cap

7	<b>CONJ STREPT</b>	<b>Streptavidin-HRP conjugate</b>	8.0 mL, ready to use, colored reddish, purple cap
8	<b>SUBSTR</b>	<b>Substrate</b> < 0.1 % 3,3',5,5'-tetramethyl-benzidine	13 mL, ready to use, colorless, blue cap
9	<b>STOP</b>	<b>Stop solution</b> 0.25 M sulphuric acid	13 mL, ready to use, colorless, yellow cap
10	<b>DIL BOTTLE</b>	<b>Empty bottle</b> for dilution of elastase-antibody biotinylated	empty, red cap
11		<b>Instructions for Use</b>	1 piece
12		<b>Certificate of Analysis</b>	1 piece

## Additional Materials and Aids Required for the Test Procedure

Adjustable single-channel-micropipette • 8-channel pipette or multi-pipette with pipette tips • 8-channel wash comb with vacuum pump and waste bottle or microplate washer • microplate reader with 450 nm measuring filter and  $\geq 620$  nm reference filter • deionized water • measuring cylinder • tubes for sample preparation • precision scale for sample weighing • Stool Preparation Kit (GZ3008 / ACS-001) or Extraction buffer (10x) (G09038P01 / ACS-002) • incubator

## Important Information

**This device is for *in-vitro* diagnostic use only.** Follow the instructions carefully. The kit may be performed by laboratory professionals only.

The instructions for use must be strictly followed. Do not use reagents from damaged packages or bottles. The shelf life specified must be observed. Do not mix components with reagents from other manufacturers.

**Mixing of test kit components of different lots is not allowed.**

All serious incidents occurring in relation with Pancreatic Elastase ELISA must be reported to the manufacturer and the competent authority of the EU member state in which user and/or patient are located.

### Information on Assay Procedure

All reagents should be stored at 2...8 °C. Bring all test components to room temperature before use. For larger sample series, pipetting reagents from liquid reservoirs using a multichannel pipette is recommended to avoid time delays. Follow the pipetting scheme and time schedules of the protocol. Wash solution should be allowed a minimum reaction time of 5 s in the wells per wash cycle. Remove wash buffer residues by thoroughly aspirating or tapping out the cavities! The aspiration and washing steps can be performed manually or with the help of a microplate washer or waterjet pump. Protect substrate from light!

### Safety Instructions

Reagents must not be swallowed. Contact with skin or mucous membranes should be avoided. Some reagents may contain biocides as preservative. Handle all components and patient samples as if potentially hazardous and infectious.

Additional information may be taken from the Material Safety Data Sheet.

Product contains the following hazardous component/-s:

WELLS	-	Contains material of animal origin.
WASHBUF (10x)	-	-
STD 1 - 4	-	Contains material of animal origin.
CONTROL 1	-	Contains material of animal origin.
CONTROL 2	-	Contains material of animal origin.
CONJ BIOTIN (201x)	-	Contains material of animal origin.
CONJ STREPT	-	Contains material of animal origin.
SUBSTR	EUH210	Safety data sheet available on request.
STOP	Hazard components H290	Sulphuric acid 2.5 % May be corrosive to metals.

### Limitations of the Procedure

At temperatures higher than 40 °C samples should be transported refrigerated or frozen (see Sample Shelf Life and Storage). Aqueous stool samples may result in falsely low results due to the dilution effect. In case of watery stool, another sample with more solid consistency should be obtained. If this is not possible, the stool may be pipetted (e.g. 10 µL watery stool + 1 mL Extraction buffer (1x)). Since there is a risk of false pathological (falsely low) values due to the dilution of the stool, the consistency should be noted and in case of a pathological finding, the test must be repeated with a solid stool sample.

It is possible that substitution with porcine pancreatic proteins leads to cross-reactions. Substitution therapy might therefore have to be interrupted. There is no firm evidence regarding the need to discontinue substitution therapy before performing the assay.

## Sample Treatment

### Sample Collection

Collect stool sample in suitable sampling container.

### Sample Shelf Life and Storage

Stool samples can be stored at different temperatures for the following time periods.

<u>Native human stool samples</u>		<u>Extracted sample supernatants</u>	
storage	Duration	storage	duration
2...8 °C	7 days	2...8 °C	3 days
20...22 °C	5 days	-18...-27 °C	1 month
40 °C	2 days		

Stool samples and extracted sample supernatant may be frozen and thawed up to two times. It is recommended to aliquot samples to avoid frequent freeze and thaw cycles.

### Sample Preparation

For sample weighing:

1. Preparation of Extraction buffer (10x) by 1 : 10 dilution (= 1+9) with deionized water (e.g. 50 mL Extraction buffer (10x) + 450 mL deionized water)
2. Mixing of stool samples of human origin 1 : 101 with Extraction buffer (1x) (e.g. 10 mg stool + 1 mL Extraction buffer (1x)).
3. Homogenization of the samples by means of a rotary shaker and subsequent sedimentation of solid constituents.
4. Removal of the supernatants after 15 - 30 min sedimentation at room temperature or alternatively at 2...8 °C overnight.

Using the Stool Preparation Kit:

Optionally, Stool Preparation Kit (GZ3008/ACS-001) can be used for weighing stool samples, according to its instructions for use.

Final dilution of sample:

Sample supernatants will be diluted 1 : 201 (1 + 200) with wash buffer (1x) (e.g. 25 µL supernatant + 5 mL wash buffer (1x)). The prepared sample dilutions are now ready to use for testing.

**Note:** If salt crystals have precipitated in the extraction or wash buffer, dissolve them before use by gently shaking and heating in a water bath at 30 - 35 °C. For analysis, the collected supernatants can be stored at -18 °C for up to 4 weeks. It is important that supernatants are separated from the sediment during storage.

## Reagent Treatment

### Reagent Shelf Life and Storage

The complete test kit with sealed reagent bottles and microtiter strips is stable until the printed expiration date when stored at 2...8 °C. All opened test kit components are stable for up to 3 months if stored properly at 2...8 °C. The diluted extraction buffer can be used at 2...8 °C up to 1 month. The diluted wash buffer can be used at 2...8 °C for up to 1 month. It is recommended to prepare a fresh ready to use dilution of the biotin conjugate prior to use. Further investigations have shown that storage at 2...8 °C is possible for up to 1 month.

### Reagent Preparation

Microtiter plate with breakable 8-well strips vacuum sealed with desiccant. Allow packaging to reach room temperature before opening. Protect unused wells from moisture and store refrigerated with desiccant in the original bag carefully resealed.

Dilute wash buffer (10x) 1 : 10 with deionized water.

Example: 10 mL wash buffer (10x) + 90 mL deionized water.

Dilute concentrated biotin conjugate (201x) with wash buffer (1x) 1 : 201. Use the empty bottle with red cap included in the kit for preparation of the ready to use biotin conjugate.

Note: Empty bottle is for single use only.

Examples: (1 part biotin conjugate (201x) + 200 parts wash buffer (1x)).

– 5 µL biotin conjugate (201x) in 1 mL wash buffer (1x) (for 2 microtitration strips)

– 25 µL biotin conjugate (201x) in 5 mL wash buffer (1x) (for 10 microtitration strips)

## Assay Procedure

1. Allow test reagents and required number of wells to reach room temperature (RT). Shake reagents gently before use. Avoid foaming.
2. Pipette 50 µL Blank (=zero standard, diluted wash buffer)  
50 µL 

STD	1 - 4
-----	-------

 Standards 1-4  
50 µL 

CONTROL	1
---------	---

 Control 1  
50 µL 

CONTROL	2
---------	---

 Control 2  
50 µL extracted and diluted stool specimen each.
3. Incubate for 60 min at 37 °C.
4. Decant, then wash each well 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
5. Pipette 50 µL 

CONJ BIOTIN
-------------

 diluted biotin conjugate per cavity.
6. Incubate for 30 min at 37 °C.
7. Decant, then wash each well 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
8. Pipette 50 µL 

CONJ STREPT
-------------

 Streptavidin-HRP conjugate per cavity.
9. Incubate for 30 min at 37 °C.
10. Decant, then wash each well 3 x with 300 µL diluted wash buffer. Tap dry onto absorbent paper if necessary.
11. Pipette 100 µL 

SUBSTR
--------

 Substrate per cavity.
12. Incubate for 20 min at 37 °C **protected from light**.
13. Pipette 100 µL 

STOP
------

 Stop solution per cavity, mix gently.
14. Read OD at 450 nm for measuring and  $\geq 620$  nm as reference with a microplate reader within 30 min following reaction stop.

## Evaluation of Results

### Quantitative Evaluation

A calibrator curve is generated by plotting the measured absorbances of standards 1-4 (y-axis) against the corresponding antigen concentrations (x-axis).

The 4-parameter regression model is recommended to create the standard curve.

By referring to the standard curve, the measured absorbances of the diluted samples are converted to µg pancreatic elastase/g stool (antigen concentration in µg/g stool). The indicated concentrations of the calibrators already take into account the standard dilution factor of the stool samples of 1 : 201.

The test run is valid if:

- the mean OD of the Blank (=zero standard, diluted wash buffer) is  $< 0.15$
- the mean concentration of Control 1 is between 70 µg/g and 110 µg/g
- the mean concentration of Control 2 is between 160 µg/g and 240 µg/g

If the above-mentioned quality criteria are not met, test should be repeated strictly following the test procedure (incubation times and temperatures, sample and wash buffer dilution, wash steps, etc.). In case of repeated failure of the quality criteria contact the manufacturer.

## Interpretation of Results

Normal exocrine pancreatic function	> 200 µg elastase/g stool
Mild/moderate exocrine pancreatic insufficiency	100 – 200 µg elastase/g stool
Severe exocrine pancreatic insufficiency	< 100 µg elastase/g stool

## Performance Characteristics

### Precision

To determine the precision (Intra-Assay-, Inter-Assay- and Lot-to-Lot-coefficient of variation), 8 samples were determined on 20 days by two operators in duplicates within 3 batches of Pancreatic Elastase ELISA.

Sample	Intra-Assay-coefficient of variation		Inter-Assay-coefficient of variation		Lot-to-Lot-coefficient of variation	
	$\bar{x}$ µg elastase/g stool	CV (%)	$\bar{x}$ µg elastase/g stool	CV (%)	$\bar{x}$ µg elastase/g stool	CV (%)
1	66.6	6.4	66.6	10.2	68.8	11.2
2	91.9	2.7	91.9	7.0	93.6	8.8
3	105.9	4.2	105.9	7.6	105.9	9.3
4	200.8	3.3	200.8	9.0	201.8	9.4
5	207.5	2.8	207.5	6.6	208.3	7.1
6	319.5	3.4	319.5	7.3	319.1	8.0
7	438.3	2.6	438.3	4.8	433.6	5.9
8	506.5	2.9	506.5	6.2	499.6	8.2

### Linearity

Using one batch of Pancreatic Elastase ELISA, a linear range of 0.3 to 231.6 µg elastase/g stool was determined with a coefficient of determination of 0.997.

### Limits of Detection (LoB, LoD, LoQ, High-Dose-Hook-Effect)

The Limit of Blank (LoB) of Pancreatic Elastase ELISA has been determined by 80-fold determination of diluted wash buffer at 0.6 µg elastase/g stool. The Limit of Detection (LoD) of Pancreatic Elastase ELISA has been determined by 80-fold determination of a diluted stool sample within the concentration range of LoB to 4-fold LoB at 4.0 µg elastase/g stool. The Limit of Quantification (LoQ) of Pancreatic Elastase ELISA has been determined by examination of 5 stool dilutions in a previously empirically calculated range of concentration at 61.4 µg elastase/g stool.

Due to successive test processing a high-dose hook effect is excluded, which has been proven experimentally.

### Sensitivity and Specificity

Sensitivity and specificity of Pancreatic Elastase ELISA have been determined in comparison to a commercially available ELISA.

For the determination of the sensitivity a sample collective of n = 102 clinically characterized stool samples were examined. Results were compared in the following four-cross-tabulation. The collective is composed of patients diagnosed with pancreatic cancer, diabetes mellitus, chronic pancreatitis and cystic fibrosis.

n = 102 clinical stool samples	Comparative ELISA > 200 µg elastase/g stool	Comparative ELISA < 200 µg elastase/g stool
Pancreatic Elastase ELISA > 200 µg elastase/g stool	47	2
Pancreatic Elastase ELISA < 200 µg elastase/g stool	5	48

Sensitivity: 96.0%

For the determination of the specificity a cohort of n = 200 healthy stool sample donors were examined and the results compared in the following four-cross-tabulation.

n = 200 healthy stool sample donors	Comparative ELISA > 200 µg elastase/g stool	Comparative ELISA < 200 µg elastase/g stool
Pancreatic Elastase ELISA > 200 µg elastase/g stool	195	2
Pancreatic Elastase ELISA < 200 µg elastase/g stool	1	2

Specificity: 99.5%

### Recovery

To determine the recovery of concentration measurements for Pancreatic Elastase ELISA, a weakly reactive sample (~100 µg elastase/g stool) was mixed with a highly reactive sample (> 300 µg elastase/g stool) in the ratio of 1 part to 9 parts and the deviation between calculated elastase concentration and measured elastase concentration was determined. The deviation of calculated to measured elastase concentration was between 9.0 and 16.1 %.

### Cross reactivity

The following microorganisms were diluted in wash buffer (1x) at  $\geq 10^8$  colony-forming units (CFU) per mL, alternatively clinical isolates diluted 1 : 20.000 were tested in the Pancreatic Elastase ELISA. No cross reactivity has been shown.

<i>Aeromonas hydrophila</i>	ATCC 7966	<i>Proteus vulgaris</i>	ATCC 8427
<i>Bacillus cereus</i>	ATCC 11778	<i>Pseudomonas aeruginosa</i>	ATCC 10145
<i>Bacillus subtilis</i>	ATCC 6633	<i>Pseudomonas fluorescens</i>	ATCC 13525
<i>Bacteroides fragilis</i>	ATCC 25285	<i>Pseudomonas putida</i>	ATCC 49128
<i>Campylobacter coli</i>	ATCC 33559	<i>Salmonella anatum</i>	ATCC 9270
<i>Campylobacter fetus</i>	ATCC 27374	<i>Salmonella enterica</i> Serovar <i>thyphimurium</i>	ATCC 14028
<i>Campylobacter jejunii</i>	ATCC 32291	<i>Salmonella enterica</i> ssp. <i>galolyticus/enteritidis</i>	ATCC 13076
<i>Campylobacter lari</i>	ATCC 35221	<i>Salmonella infantis</i>	ATCC 51741
<i>Campylobacter upsaliensis</i>	ATCC 43954	<i>Salmonella paratyphi A</i>	ATCC 11511
<i>Candida albicans</i>	ATCC 10231	<i>Salmonella paratyphi B</i>	ATCC 8759
<i>Citrobacter freundii</i>	ATCC 8090	<i>Salmonella paratyphi C</i>	Nr. 2 Pasteur
<i>Clostridium difficile</i>	VPI 10463	<i>Shigella flexneri</i>	ATCC 12022
<i>Clostridium sordelli</i>	ATCC 9714	<i>Shigella sonnei</i>	ATCC 25931
<i>Enterobacter aerogenes</i>	ATCC 13048	<i>Staphylococcus aureus</i>	ATCC 25923
<i>Enterobacter cloacae</i>	ATCC 13047	<i>Staphylococcus epidermidis</i>	ATCC 12228
<i>Enterococcus faecalis</i>	ATCC 29212	<i>Streptococcus agalactiae</i>	ATCC 13813
<i>Escherichia coli</i>	ATCC 25922	<i>Vibrio cholerae</i>	clinical isolate
<i>Klebsiella pneumonia</i>	ATCC 13883	<i>Vibrio parahaemolyticus</i>	ATCC 17802
<i>Lactococcus lactis</i>	DSM 20481	<i>Yersinia enterocolitica</i> 0:3	clinical isolate
<i>Morganella morganii</i>	ATCC 25830	<i>Yersinia enterocolitica</i> 0:9	clinical isolate
<i>Peptostreptococcus anaerobius</i>	ATCC 27337	<i>Yersinia enterocolitica</i> RKI 0803733	clinical isolate
<i>Proteus mirabilis</i>	ATCC 29906	<i>Yersinia enterocolitica</i> Y11	clinical isolate

## Interference

None of the following substances added to analyte-containing samples and weakly reactive samples showed a significant impact on the test result: The mentioned concentrations refer to 30 mg faeces.

(-)-Scopolamine n-Butyl bromide (0.5 %, Buscopan®), Barium sulfate (5 %), Biotin (0.001%), Bismuth(III) subsalicylate (0.5 %, Pepto-Bismol), Cyclamat (5 %), Diclofenac (0.5 %), hemoglobin human (5 %), blood human (5 %), Hylak® N (5 %), Iberogast® (5 %), Immodium® akut duo (0.03/1.9 %), Loperamid hydrochloride (5 %, Loperamid-CT akut), Metronidazole (0.5 %), Mucin (5 %), Nexium® (0.25 %), Nifuroxazide (0.5 %, Pentofuryl®), palmitic acid (20 %), pancreas powder from pig (5%, Pangrol® 25.000), pancreas powder from pig (2.5%, Pankreatan® 36.000), pancreas powder from pig (2.5%, Panzytrat® 40.000), Pankreatin (5%, Kreon® 35.000), Perenterol forte (0.5 %), Rennie® (20 %), Simigel® (1 %), stearic acid (20 %), Vancomycin (0.5 %).

## Clinical Performance

A total of 15 clinical studies have been performed on Pancreatic Elastase ELISA in the European area and other areas between 2003 and 2021. The clinical performance of these studies has been clearly demonstrated [1-15].

Pancreatic Elastase ELISA was used to test a selection of clinical samples within an external clinical performance study, from which the following performance parameters were calculated retrospectively.

<b>Diagnostic sensitivity</b>	60 % Chronic pancreatitis 79 % Cystic fibrosis 41 % Diabetes mellitus 61 % Pancreatic cancer <b>61% Total</b>
<b>Diagnostic specificity</b>	<b>88 % Total</b>
<b>Positive predictive value</b>	59.6 % Chronic pancreatitis 79.2 % Cystic fibrosis 41.2 % Diabetes mellitus 61.1 % Pancreatic cancer <b>91.9% Total</b>
<b>Negative predictive value</b>	31.3 % Chronic pancreatitis 9.8 % Cystic fibrosis 17.9 % Diabetes mellitus 13.2 % Pancreatic cancer <b>51.7% Total</b>
<b>Likelihood-ratio</b>	5.2 = (LR+), 0.5 = (LR-) Chronic pancreatitis 6.9 = (LR+), 0.2 = (LR-) Cystic fibrosis 3.6 = (LR+), 0.7 = (LR-) Diabetes mellitus 5.3 = (LR+), 0.4 = (LR-) Pancreatic cancer <b>6 = (LR+), 0.4 = (LR-) Total</b>

From a comparative study using clinically characterized samples of Pancreatic Elastase ELISA with comparative ELISA, the following results are obtained:

<b>Diagnostic sensitivity</b>	92.6 %
<b>Diagnostic specificity</b>	97.6 %
<b>Positive predictive value</b>	89.3 %
<b>Negative predictive value</b>	98.4 %
<b>Likelihood-ratio</b>	38.27 (LR+) 0.08 (LR-)


## Application

### Automated Workflow

Performing Pancreatic Elastase ELISA on fully automated microplate processors (e.g. DS2® or DSX) may cause elevated absorbances in comparison to the manual procedure due to individual differences concerning wash procedures and general technical specifications of the equipment.


### Correlation: manual – automated workflow

A panel of 160 stool specimens was processed manually and automatically in parallel in 3 batches (DSX®, Dynex Technologies). The average correlation coefficient was calculated at  $r = 0.986$ .

## Equivalence Study Stool Preparation Systems

A panel of 120 stool specimens was investigated in parallel, by sample preparation with Stool Preparation Kit (ACS-001) and RIDA®TUBE Calprotectin (R-Biopharm, GZ3016) in one batch of the Pancreatic Elastase ELISA (G09040). The average correlation was calculated at  $r = 0.974$ .

### Sample Preparation (optional)

Using the RIDA®TUBE Calprotectin (R-Biopharm, GZ3016):

As an alternative to stool sample weighing the stool preparation kit RIDA®TUBE Calprotectin from R-Biopharm can be used.

Final dilution of sample:

Sample supernatants have to be diluted 1 : 81 (1 + 80) with wash buffer (1x) (e.g. 25 µL supernatant + 2 mL wash buffer (1x)). The prepared sample dilutions are ready to use for testing.

### Sample Shelf Life and Storage for RIDA®TUBE Calprotectin (R-Biopharm, GZ3016)

Stool samples in extraction tubes or extracted sample supernatants can be stored at different temperatures for the following time periods.

<u>Stool samples in extraction tubes</u>		<u>Extracted sample supernatants</u>	
storage	duration	storage	duration
2...8 °C	3 days	2...8 °C	overnight
23 °C	<b>not possible</b>	23 °C	<b>not possible</b>
-20 °C	<b>not possible</b>	-20 °C	<b>not possible</b>

**Note:** It is recommended to perform sample extraction, sample dilution with RIDA®TUBE Calprotectin (R-Biopharm, GZ3016) and sample testing in Pancreatic Elastase ELISA on the same day. In addition, it is possible to store the sample extracts generated with the RIDA®TUBE Calprotectin (R-Biopharm, GZ3016) in the extraction tube for up to 3 days at 2...8 °C before testing. Separated sample supernatants may be stored overnight at 2...8°C before testing. It is not possible to keep sample extracts, sample supernatants or sample dilutions for a longer period than specified or to store them at higher or lower temperatures than 2...8 °C before testing.

**Note:** Further information can be found in the instructions for use of the RIDA®TUBE Calprotectin (R-Biopharm, GZ3016).

## History of Changes

Version	Section	Modifications
2022-10_v01	Entire document	Document creation

## References

Citations belonging to the section Clinical Performance.

1. Pietzner, M., et al., Exocrine Pancreatic Function Modulates Plasma Metabolites Through Changes in Gut Microbiota Composition. *J Clin Endocrinol Metab*, 2021. **106**(5): p. e2290-e2298.
2. Boga, S., et al., Liver and pancreas: 'Castor and Pollux' regarding the relationship between hepatic steatosis and pancreas exocrine insufficiency. *Pancreatology*, 2020. **20**(5): p. 880-886.
3. Aksoz, Z., et al., The easy way of evaluating exocrine pancreatic insufficiency in type 2 diabetes : listen to the patients' complaints and look in their eyes! *Acta Gastroenterol Belg*, 2020. **83**(3): p. 407-412.
4. Akay, S., B. Sirin, and B. Unsal, Fecal Elastase Levels Predict Honeycombing in Pancreas Detected with Endoscopic Ultrasound. *Can J Gastroenterol Hepatol*, 2018. **2018**: p. 4625247.
5. Ellemunter, H., et al., Fecal Calprotectin in Cystic Fibrosis and Its Relation to Disease Parameters: A Longitudinal Analysis for 12 Years. *J Pediatr Gastroenterol Nutr*, 2017. **65**(4): p. 438-442.
6. Aoufi Rabih, S., et al., Exocrine pancreatic insufficiency and chronic pancreatitis in chronic alcoholic liver disease: coincidence or shared toxicity? *Pancreas*, 2014. **43**(5): p. 730-4.
7. Fernandez-Lorenzo, A.E., et al., V232D mutation in patients with cystic fibrosis: Not so rare, not so mild. *Medicine (Baltimore)*, 2018. **97**(28): p. e11397.
8. Demoulin, N., et al., Enteric hyperoxaluria in chronic pancreatitis. *Medicine (Baltimore)*, 2017. **96**(19): p. e6758.
9. Shin, Y.C., et al., Effects of pancreatectomy on nutritional state, pancreatic function, and quality of life over 5 years of follow up. *J Hepatobiliary Pancreat Sci*, 2020.
10. Srivastava, A., et al., Prevalence and predictive factors of undernutrition and low bone mineral density in children with chronic pancreatitis. *Pancreatology*, 2021. **21**(1): p. 74-80.
11. Singh, N., et al., Antioxidants for Pancreatic Functions in Chronic Pancreatitis: A Double-blind Randomized Placebo-controlled Pilot Study. *J Clin Gastroenterol*, 2020. **54**(3): p. 284-293.
12. Chowdhury, S.D., et al., Pancreatic exocrine insufficiency: Comparing fecal elastase 1 with 72-h stool for fecal fat estimation. *Indian J Gastroenterol*, 2016. **35**(6): p. 441-444.
13. O'Sullivan, B.P., et al., Evolution of pancreatic function during the first year in infants with cystic fibrosis. *J Pediatr*, 2013. **162**(4): p. 808-812 e1.
14. Malluta, E.F., et al., Pancreatic endosonographic findings and clinical correlation in Crohn's disease. *Clinics (Sao Paulo)*, 2019. **74**: p. e853.
15. Bartels, R.H., et al., Pancreatic Enzyme Replacement Therapy in Children with Severe Acute Malnutrition: A Randomized Controlled Trial. *J Pediatr*, 2017. **190**: p. 85-92 e2.

Common citations.

1. Basturk, A., Y. Curek, R. Felek, G. Celmeli and R. Artan (2021). "Exocrine pancreas functions in children with type 1 diabetes mellitus." *Arab J Gastroenterol* **22**(3): 236-239.
2. Beyer, G., Hoffmeister, A., Michl, P., Gress, T. M., Huber, W., Algül, H., ... & Mayerle, J. (2022). S3-Leitlinie Pankreatitis—Leitlinie der Deutschen Gesellschaft für Gastroenterologie, Verdauungs-und Stoffwechselkrankheiten (DGVS)—September 2021—AWMF Registernummer 021-003. *Zeitschrift für Gastroenterologie*, 60(03), 419-521.
3. Dominici, R. and C. Franzini (2002). "Fecal elastase-1 as a test for pancreatic function: a review." *Clin Chem Lab Med* **40**(4): 325-332.
4. Gonzales, A. C., S. M. Vieira, R. L. Maurer, F. A. Silva and T. R. Silveira (2011). "Use of monoclonal faecal elastase-1 concentration for pancreatic status assessment in cystic fibrosis patients." *J Pediatr (Rio J)* **87**(2): 157-162.
5. Gullo, L., L. Graziano, S. Babbini, A. Battistini, R. Lazzari and R. Pezzilli (1997). "Faecal elastase 1 in children with cystic fibrosis." *Eur J Pediatr* **156**(10): 770-772.
6. Keim, V., Teich, N., & Moessner, J. (2003). Clinical value of a new fecal elastase test for detection of chronic pancreatitis. *Clinical laboratory*, 49(5-6), 209-215.
7. Löser, C., A. Möllgaard and U. R. Fölsch (1996). "Faecal elastase 1: a novel, highly sensitive, and specific tubeless pancreatic function test." *Gut* **39**(4): 580-586.
8. Naruse, S., H. Ishiguro, S. B. Ko, T. Yoshikawa, T. Yamamoto, A. Yamamoto, S. Futakuchi, H. Goto, Y. Saito and S. Takahashi (2006). "Fecal pancreatic elastase: a reproducible marker for severe exocrine pancreatic insufficiency." *J Gastroenterol* **41**(9): 901-908.
9. Rothenbacher, D., M. Löw, P. D. Hardt, H. U. Klör, H. Ziegler and H. Brenner (2005). "Prevalence and determinants of exocrine pancreatic insufficiency among older adults: results of a population-based study." *Scand J Gastroenterol* **40**(6): 697-704.