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---

**Intended Use**

The sample preparation set enables you to easily and reliably prepare stool samples for the Pancreatic Elastase ELISA developed and produced by BIOSERV Diagnostics.

---

**Instructions for Reagent Preparation**

The sample preparation set is already filled with extraction buffer. Bring the sample preparation set to room temperature before use. In case of a precipitation of salt crystals in the extraction buffer these precipitates should be disintegrated by warming in a water bath at 30 - 35 °C (86 – 95 °F) and shaking for 5 - 10 min. A variability of the filling volume of 5 % (1,425 to 1,575 ml) could arise due to technological reasons. The measuring results in the Pancreatic Elastase ELISA are not considerably influenced thereof.

---

**Warnings and Precautions**

Sample material collected from patients should principally be considered as potentially infectious. Handling and disposal should be in accordance with the procedures defined by the national biohazard safety guidelines or regulations, where these exist. Samples collected from risk patients should be clearly marked and should be handled under special security precautions. In particular suitable protective laboratory clothing and gloves should be used.

---

**Storage Instructions and Shelf Life Information**

The extraction buffer included in the sample preparation set may be stored at a temperature of 2 - 8 °C (36 °F – 46 °F) until the expiry date.

---

**Information with Reference to the Handling of the Sample Material**

1. The stool samples to be examined can be stored before and after the examination at different temperatures for the following time-spans:

|  |                 |
|--|-----------------|
| Environmental temperature up to 40 °C (104 °F)     | up to five days |
| Refrigerator temperature 2 – 8 °C (36 – 46 °F)     | up to one week  |
| Household freezer temperature ca. -18 °C (-0.4 °F) | up to one year  |
2. If stool samples which are extracted with the sample preparation set must be stored for later examinations, it is necessary in every case to detach the supernatant from the stool sediment (by means of pipetting). The supernatant can be stored under the following conditions without influence on the elastase concentration:

|  |                   |
|--|-------------------|
| Refrigerator temperature 2 – 8 °C (36 – 46 °F) | two to three days |
| Household freezer temperature -18 °C (-0.4 °F) | up to one month   |

---

**Procedure of the Sample Preparation**

1. The sample preparation set is already filled with extraction buffer and ready for use.
2. The sample preparation set is to be marked with the patients data and to be put into a test tube holder.
3. Unscrew the yellow cap and take the integrated stick totally out of the lower part of the set.
4. Stick the lower end of the stick into the stool sample. The notches at the end of the stick must be totally filled with stool.
5. Put the stick back into the sample preparation set and screw down tightly. All surplus stool is wiped off by the cone in the middle part of the sample preparation set. This ensures the extraction of a defined stool quantity.
6. Mix the sample with the extraction buffer vigorously for at least 2 min, for example by using a Vortex™.
7. Please note: All stool contained in the notches at the end of the stick must be brought into suspension.
8. Leave the sample preparation set for 15 min for extraction at room temperature. It is also possible to let the extraction take place over night at 2 – 8 °C.
9. Repeat the mixing procedure (see point 6).
10. Leave the sample preparation set again for 15 min for further extraction and for sedimentation of the suspended stool particles.
11. Unscrew carefully without dispersing the sediment the blue middle part of the sample preparation system without unscrewing the yellow upper part. Discard the yellow upper part and the blue middle part.
12. Take out the volume of the supernatant you need for carrying out the test.

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