

## RIDA® GENE *Pneumocystis jirovecii*

**REF** PG1905



## 1. Intended use

For *in vitro* diagnostic use. RIDA<sup>®</sup>GENE Pneumocystis jirovecii is a multiplex real-time PCR for the direct qualitative and quantitative detection of *Pneumocystis jirovecii* from human bronchoalveolar lavage fluid (BAL).<sup>1,2</sup>

The RIDA<sup>®</sup>GENE Pneumocystis jirovecii multiplex real-time PCR is intended to use as an aid in diagnosis for respiratory infections caused by *Pneumocystis jirovecii*.

## 2. Summary and explanation of the test

*Pneumocystis jirovecii* (former *P. carinii*) belongs to the family of *Pneumocystidaceae* and may lead to an interstitial pneumonia. Opportunistic infections are a major problem in immunocompromised patients, for example HIV/AIDS patients, chemotherapy-treated patients and patients receiving an organ transplant. *Pneumocystis jirovecii* causes respiratory infections and is the most common opportunistic illness in HIV-infected people. *Pneumocystis jirovecii* does not cause any harm in healthy people and is widely spread among the normal population. However, immunocompromised people infected with *Pneumocystis jirovecii* develop pneumonia with symptoms including dry cough, shortness of breath, tachypnea and fever.<sup>3</sup> Although HAART therapy decreased the *Pneumocystis jirovecii* incidence by 3.4 % per year after 1996, it is estimated that still 9 % among hospitalized HIV/AIDS patients and 1 % among solid organ transplant recipients are infected.<sup>4</sup> According to the Center for Disease Control (CDC), *Pneumocystis jirovecii* causes 100 % mortality in patients without treatment and the mortality rate in immunocompromised patients is between 5 % - 40 % in treated patients.<sup>4</sup> The mortality from *Pneumocystis jirovecii* in HIV-uninfected patients can be as high as 40 %.<sup>5</sup> Until now, detection of *Pneumocystis jirovecii* was done by immunofluorescence staining. However, due to its low sensitivity this is now substituted by PCR.<sup>2</sup>

## 3. Test principle

RIDA<sup>®</sup>GENE Pneumocystis jirovecii is a multiplex real-time PCR for the direct, qualitative and quantitative detection of *Pneumocystis jirovecii* from human bronchoalveolar lavage fluid (BAL).

After DNA isolation, amplification of the gene fragment (if present) specific for *Pneumocystis jirovecii* (mt LSU; large subunit) occurs. The amplified targets are detected with hydrolysis probes, which are labeled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the presence of a target the probes hybridize to the amplicons. During the extension step, the **Taq-polymerase** breaks the reporter-quencher proximity. The reporter emits a fluorescent signal, which is detected by the optical unit of a real-time PCR instrument. The fluorescence signal increases with the amount of formed amplicons. With the standards, **Standard A**, **Standards B** and **Standards C**, included in the kit, it

is possible to quantify the results. The RIDA® GENE Pneumocystis jirovecii multiplex real-time PCR kit contains an **Internal Control DNA** (ICD) that detects PCR inhibition, monitors reagent integrity and confirms that nucleic acid extraction was sufficient.

#### 4. Reagents provided

**Tab. 1: Reagents provided (Reagents provided in the kit are sufficient for 100 determinations)**

Kit Code	Reagent	Amount		Lid Color
1	Reaction Mix	2x	1050 µl	yellow
2	Taq-Polymerase	1x	80 µl	red
D	Internal Control DNA	2x	1700 µl	orange
N	No Template Control	1x	450 µl	white
P	Positive Control	1x	200 µl	blue
10 <sup>^1</sup>	Standard A	1x	100 µl	dark blue
10 <sup>^3</sup>	Standard B	1x	100 µl	dark blue
10 <sup>^5</sup>	Standard C	1x	100 µl	dark blue

#### 5. Storage instructions

- Protect all reagents from light and store at -20 °C. All reagents can be used until the expiration date. After expiry the quality guarantee is no longer valid.
- Carefully thaw reagents before using (e.g. in a refrigerator at 2 - 8 °C).
- Reagents can sustain up to 20 freeze/thaw cycles without influencing the assay performance (e.g. after the first thawing separate it in aliquots and freeze immediately).
- During PCR preparation all the reagents should be stored cold in an appropriate way (2 - 8 °C).

## 6. Additional necessary reagents and necessary equipment

The RIDA<sup>®</sup>GENE Pneumocystis multiplex real-time PCR assay is suitable for use with following extraction platforms and real-time PCR instruments:

**Tab. 2:** Necessary equipment

Extraction platforms	
R-Biopharm	RIDA <sup>®</sup> Xtract
Promega	Maxwell <sup>®</sup> RSC
Real-time PCR instruments	
Roche	LightCycler <sup>®</sup> 2.0, LightCycler <sup>®</sup> 480II
Agilent Technologies	Mx3005P
Applied Biosystems	ABI 7500
Bio-Rad	CFX96 <sup>™</sup>
QIAGEN	Rotor-Gene Q

**Note: Only use 0.1 ml tubes on the Rotor-Gene Q (QIAGEN).**

If you want to use other extraction platforms or real-time PCR instruments please contact R-Biopharm at [mdx@r-biopharm.de](mailto:mdx@r-biopharm.de).

- RIDA<sup>®</sup>GENE Color Compensation Kit II (PG0002) for use with the LightCycler<sup>®</sup> 2.0
- RIDA<sup>®</sup>GENE Color Compensation Kit IV (PG0004) for use with the LightCycler<sup>®</sup> 480II
- Real-time PCR consumables (plates, tubes, foil)
- Centrifuge with a rotor for the reaction vials
- Vortexer
- Pipettes (0.5 - 20 µl, 20 - 200 µl, 100 - 1000 µl)
- Filter tips
- Powder-free disposal gloves
- PCR water (nuclease-free)

## 7. Precautions for users

For *in-vitro* diagnostic use.

This test must only be carried out by trained laboratory personnel. The guidelines for working in medical laboratories have to be followed. The instruction manual for the test procedure has to be followed. Do not pipet samples or reagents by mouth. Avoid contact with bruised skin or mucosal membranes. During handling reagents or samples, wear appropriate safety clothing (appropriate gloves, lab coat, safety goggles) and wash your hands after finishing the test procedure. Do not smoke, eat or drink in areas where samples or reagents are being used.

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.

- Samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after the expiration date.

All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulations for disposal.

For more details see Safety Data Sheets (SDS) at [www.r-biopharm.com](http://www.r-biopharm.com).

## 8. Collection and storage of samples

### 8.1 Sample preparation from bronchoalveolar lavage (BAL)

For DNA isolation from bronchoalveolar lavage (BAL), use a commercially available DNA isolation kit (e.g. RIDA<sup>®</sup> Xtract (R-Biopharm)) or DNA extraction system (e.g. Maxwell<sup>®</sup> RSC (Promega)). Extract DNA according to the manufacturer's instructions.

The RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* assay contains an **Internal Control DNA** that detects PCR inhibition, monitors reagent integrity and confirms that nucleic acid extraction was sufficient. The **Internal Control DNA** can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as PCR inhibition control.

If the **Internal Control DNA** is used only as a PCR inhibition control, 1 µl of the **Internal Control DNA** should be added to the Master-Mix (see Tab.4).

If the **Internal Control DNA** is used as an extraction control for the sample preparation procedure **and** as PCR inhibition control, 20 µl of the **Internal Control DNA** has to be added during extraction procedure. The **Internal Control DNA** should always be added to the specimen-lysis buffer mixture and must **not** be added directly to the specimen. We also recommend to add 1 µl of the **Internal Control DNA** to the negative control and positive control PCR Mix.

## 9. Test procedure

### 9.1 Master-Mix preparation

Calculate the total number of PCR reactions (sample and control reactions) needed. One positive control and one negative control must be included in each assay run.

We recommend to calculate an additional volume of 10 % to compensate imprecise pipetting (see Tab. 3, Tab. 4). Thaw, mix gently and centrifuge briefly the **Reaction Mix**, the **Taq-Polymerase**, the **Positive Control**, the **No Template Control**, the **Internal Control DNA** and **Standard A**, **Standard B** and **Standard C** before using. Keep reagents appropriately cold during working step (2 - 8 °C).

**Tab. 3:** Calculation and pipetting example for 10 reactions of the Master-Mix (ICD as extraction and PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
	<b>Total</b>	<b>20 µl</b>	<b>220 µl</b>

Mix the components of the Master-Mix gently and briefly spin down.

**Tab. 4:** Calculation and pipetting example for 10 reactions of the Master-Mix (ICD only as PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
D	Internal Control DNA	1.0 µl	11 µl
	<b>Total</b>	<b>21.0 µl</b>	<b>231.0 µl</b>

Mix the components of the Master-Mix gently and briefly spin down.

## 9.2 Preparation of the PCR-Mix

Pipette 20 µl of the Master-Mix in each reaction vial (tube or plate).

**Negative control:** Add 5 µl **No Template Control** to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the negative control.

**Sample:** Add 5 µl DNA extract to the pre-pipetted Master-Mix.

**Positive control:** Add 5 µl **Positive Control** to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the positive control.

**Standard (A, B, C):** Add 5 µl **Standard** (A, B, C) to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the standards.

**Note:** Using the following cyclers requires to include a standard curve in each run: ABI 7500 (Applied Biosystems) and CFX96™ (Bio-Rad).

For all other cyclers, only one sample of the standard curve (**Standard B**) has to be included in the experimental set-up as calibrator for each new real-time PCR run. Here, the application of a standard curve is only required to run once per lot number.

Cover tubes or plate. Spin down and place in the real-time PCR instrument. The PCR reaction should be started according to the PCR instrument set-up (see Tab. 5, Tab. 6, Tab. 7, Tab. 8).

### 9.3 PCR instrument set-up

#### 9.3.1 DNA real-time PCR profile

**Tab. 5:** DNA real-time PCR profile for LightCycler® series and Rotor-Gene Q

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
PCR Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Tab. 6:** DNA real-time PCR profile for Mx3005P, ABI7500, CFX96™

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
PCR Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Note:** The total copy number per reaction of **Standard A**, **Standard B** and **Standard C** has to be typed in into the setup file of the software program of the respective real-time PCR cyclers. A total volume of 5 µl DNA is used resulting in following concentrations:

**Standard A:**  $5 \times 10^1$  copies/reaction

**Standard B:**  $5 \times 10^3$  copies/reaction

**Standard C:**  $5 \times 10^5$  copies/reaction

**Note:** The standard curve can be saved on the real-time PCR cycler for each parameter. Apart from the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler, the standard curve is only required to run once per lot number. Using the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler requires to include a standard curve in each run. For all other cyclers, only one sample of the standard curve (**Standard B**) has to be included in the experimental set-up as calibrator for each new real-time PCR run. The general recommendation is, however, to test always all three standards with every run.



### 9.3.2 Universal real-time PCR profile

**Note:** The universal real-time PCR profile should only be used for DNA assays when combining RIDA® GENE DNA and RNA real-time PCR assays in one run.

**Tab. 7:** Universal real-time PCR profile for LightCycler® series

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Tab. 8:** Universal real-time PCR profile for Mx3005P, ABI7500, CFX96™ and Rotor-Gene Q

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Note:** The total copy number per reaction of **Standard A**, **Standard B** and **Standard C** has to be typed in into the setup file of the software program of the respective real-time PCR cycler. A total volume of 5 µl DNA is used resulting in following concentrations:

**Standard A:** 5 x 10<sup>1</sup> copies/reaction

**Standard B:** 5 x 10<sup>3</sup> copies/reaction

**Standard C:** 5 x 10<sup>5</sup> copies/reaction

**Note:** The standard curve can be saved on the real-time PCR cycler for each parameter. Apart from the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler, the standard curve is only required to run once per lot number. Using the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler requires to include a standard curve in each

run. For all other cyclers, only one sample of the standard curve (Standard B) has to be included in the experimental set-up as calibrator for each new real-time PCR run. The general recommendation is, however, to test always all three standards with every run.

#### 9.4 Detection channel set-up

Tab. 9: Selection of appropriate detection channels

Real-time PCR Instrument	Detection	Detection Channel	Note
Roche LightCycler® 2.0	<i>Pneumocystis jirovecii</i>	530	RIDA® GENE Color Compensation Kit II (PG0002) is required
	ICD	560	
Roche LightCycler® 480II	<i>Pneumocystis jirovecii</i>	465/510	RIDA® GENE Color Compensation Kit IV (PG0004) is required
	ICD	533/580	
Agilent Techn. Mx3005P	<i>Pneumocystis jirovecii</i>	FAM	Check that reference dye is none
	ICD	HEX	
ABI 7500	<i>Pneumocystis jirovecii</i>	FAM	Check that passive reference option ROX is none
	ICD	VIC	
Bio-Rad CFX96™	<i>Pneumocystis jirovecii</i>	FAM	-
	ICD	VIC	
Qiagen Rotor-Gene Q	<i>Pneumocystis jirovecii</i>	Green	The gain settings have to be set to 5, according to the default settings
	ICD	Yellow	

## 10. Quality control

The analysis of the samples is done by the software of the used real-time PCR instrument according to the manufacturer's instructions. Positive control and negative control have to show correct results (see Table 10, Fig. 1) in order to determine a valid run.

The **Positive Control** has a concentration of  $10^3$  copies/ $\mu$ l. In each PCR run it is used in a total amount of  $5 \times 10^3$  copies.

**Tab. 10:** For a valid run, the following conditions must be met:

Sample	Assay result	ICD Ct	Target Ct
Positive control	Positive	NA <sup>*1</sup>	See Quality Assurance Certificate
Negative control	Negative	Ct > 20	0


*\*1 No Ct value is required for the ICD to make a positive call for the positive control.*

If the positive control is not positive within the specified Ct range but the negative control is valid, prepare all new reactions including the controls.


If the negative control is not negative but the positive control is valid, prepare all new reactions including the controls.

If the required criteria are not met, following items have to be checked before repeating the test:

- Expiry of the used reagents
- Functionality of the used instrumentation
- Correct performance of the test procedure



**Fig. 1:** Correct run of the positive control (red) and negative control (green) (*Pneumocystis jirovecii*) on the LightCycler® 480II



**Fig. 2:** Standard series *Pneumocystis jirovecii* with **Standard A** (10<sup>1</sup> DNA copies per µl), **Standard B** (10<sup>3</sup> DNA copies per µl) and **Standard C** (10<sup>5</sup> DNA copies per µl) on the LightCycler® 480II

## 10.1 Validity of quantitative detection

For the validity of a quantitative diagnostic test run, all control conditions of a valid qualitative diagnostic test run must be met. Furthermore, for accurate quantification results a valid standard curve has to be generated. For a valid quantitative diagnostic test run, the following control parameter values of the standard curve should be achieved.

	Control parameter	Valid value
Roche LightCycler® 2.0	Efficiency	1,9 – 2,1
Roche LightCycler® 480II	Efficiency	1,9 – 2,1
	Slope	-3,1 – -3,6
Agilent Techn. Mx3005P	Rsq	> 0,98
	Slope	-3,1 – -3,6
ABI 7500	R <sup>2</sup>	> 0,98
	Slope	-3,1 – -3,6
Bio-Rad CFX96™	R <sup>2</sup>	> 0,98
	Slope	-3,1 – -3,6
Qiagen Rotor-Gene Q	R <sup>2</sup>	> 0,98
	M	-3,1 – -3,6

## 11. Result interpretation

The result interpretation is done according to Table 11.

**Tab. 11:** Sample interpretation

Target genes		
<i>Pneumocystis jirovecii</i>	ICD	Ergebnis
positive	positive/negative	<i>Pneumocystis jirovecii</i> detected
negative	positive	Target genes not detected
negative	negative	Invalid

*Pneumocystis jirovecii* is detected, if the sample DNA and the Internal Control DNA show an amplification signal in the detection system.

*Pneumocystis jirovecii* is also detected, if the sample DNA shows an amplification signal but none for the Internal Control DNA in the detection system. The detection of the internal amplification control is not necessary because high concentrations of the amplicon can cause a weak or absent signal of the Internal Control DNA.

*Pneumocystis jirovecii* is not detected, if the sample DNA shows no amplification signal, but an amplification signal for the Internal Control DNA in the detection system. An inhibition of the PCR reaction can be excluded by the detection of the Internal Control DNA.

A sample is invalid, if the sample DNA and Internal Control DNA show no amplification signal in the detection system. The sample contains a PCR inhibitor. The extracted sample needs to be further diluted with PCR water (1:10) and re-amplified, or the isolation and purification of the sample has to be improved.

### 11.1 Quantification of samples

To quantify *Pneumocystis jirovecii* positive samples, a standard curve with the Standard A, Standard B and Standard C has to be performed separately. The standard curve measurement has to be saved separately. However, the same standard curve measurement can be used in all runs with products from the same lot number by importing the saved experiment.

**Note: This is not valid for the following cyclers: ABI 7500 (Applied Biosystems) and CFX96™ (Bio-Rad). Here, a standard curve has to be measured with each run.**

**For all other cyclers, one sample of the standard curve (Standard B) has to be included in the experimental set-up as calibrator for each new real-time PCR run.**

To quantify *Pneumocystis jirovecii* positive samples, all standard samples (A, B and C), the positive control and the negative control as well as the unknown samples to

be quantified, have to be selected and analyzed according to the instructions of the cyclor manufacturer. Correct quantification results are only reliable if Ct-values of the *Pneumocystis jirovecii* specific target gene (mt LSU; large subunit) can be detected within the standard Ct-range.

With the quantitative RIDA<sup>®</sup>GENE Pneumocystis jirovecii multiplex real-time PCR the amount of DNA in copies/reaction of the parameter is calculated. The conversion in copies/ml is done with a correction factor K and takes into account the dilutions of the extraction procedure (dependent on the extraction kit used) and the PCR set-up.

The conversion of the result of the quantitative RIDA<sup>®</sup>GENE Pneumocystis jirovecii multiplex real-time PCR in copies/ml is calculated with following formula:

$$C \text{ [copies/ml]} = c \text{ [copies/reaction]} \times K$$

- C [copies/ml] - concentration of sample in copies/ml sample
- c [copies/reaction] - DNA concentration in PCR reaction (result of quantitative PCR)
- K - correction factor

For the calculation of the correction factor, following information has to be considered:

- Sample dilution
- Starting volume of sample for DNA extraction
- DNA extract from total eluate used for PCR reaction

**Tab. 12:** Example of calculation of correction factor K using the Maxwell<sup>®</sup> RSC (Promega)

Description	Factor
300 µl sample put into extraction*, eluted in 60 µl final volume	No factor
5 µl DNA extract put into PCR (total eluate 60 µl =1/12)	X 12
300 µl sample scaled up to 1 ml*	X 3. $\bar{3}$
Correction factor K for <i>Pneumocystis jirovecii</i>	40

\* Result is based on 1 ml BAL starting material

**Note:** For further information on quantification of the samples please contact [pcr@r-biopharm.de](mailto:pcr@r-biopharm.de).

## 12. Limitations of the method

1. The result of molecular analysis should not lead to the diagnosis, but always be considered in the context of medical history and symptoms of the patient.
2. This assay is only validated for human bronchoalveolar lavage fluid (BAL).
3. Inappropriate specimen collection, transport, storage and processing or a pathogen load in the specimen below the analytical sensitivity can result in false negative results.
4. The presence of PCR inhibitors may cause invalid results.
5. Mutations or polymorphisms in primer or probe binding regions may affect detection of new variants resulting in a false negative result with the RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* assay.
6. As with all PCR based *in vitro* diagnostic tests, extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
7. A positive test result does not necessarily indicate the presence of viable organisms. However, a positive result is indicative for the presence of target genes ((mt LSU; large subunit).

## 13. Performance characteristics

### 13.1 Clinical performance

In a retrospective clinical validation study we analyzed 203 extracted specimens (BAL) with the RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* assay and an in-house real-time PCR assay in a laboratory in Germany.

**Tab. 13:** Correlation of the *Pneumocystis jirovecii* results with the RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* multiplex real-time PCR and reference in-house real-time PCR.

		In-house real-time PCR		Total	Comments
		Positive	Negative		
RIDA <sup>®</sup> GENE <i>Pneumocystis jirovecii</i>	Positive	28	2 <sup>b)</sup>	30	Pos. agreement: 91.8 %
	Negative	3 <sup>a)</sup>	170	173	Neg. agreement: 98.6 %
	Total	31	172	203	

a) Three (3) samples are below the limit of detection (LOD) of the RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* assay with a *C<sub>p</sub>* value > 35 in the reference in-house real-time PCR assay.


b) Two (2) samples are below the limit of detection (LOD) of the reference in-house real-time PCR assay with a *C<sub>p</sub>* value > 33 in the RIDA<sup>®</sup>GENE *Pneumocystis jirovecii* assay



## 13.2 Analytical sensitivity

The RIDA<sup>®</sup> GENE *Pneumocystis jirovecii* multiplex real-time PCR has a detection limit of  $\geq 10$  DNA copies per reaction for *Pneumocystis jirovecii*.

The following figure 2 shows a dilution series of *Pneumocystis jirovecii* ( $10^5 - 10^1$  DNA copies per  $\mu\text{l}$ ) on the LightCycler<sup>®</sup> 480II.



**Fig. 2:** Dilution series *Pneumocystis jirovecii* ( $10^5 - 10^1$  DNA copies per  $\mu\text{l}$ ) on the LightCycler<sup>®</sup> 480II

The detection limit of the whole procedure depends on the sample matrix, DNA extraction and DNA concentration.

### 13.3 Analytical specificity

The RIDA<sup>®</sup> GENE *Pneumocystis jirovecii* multiplex real-time PCR is specific for *Pneumocystis jirovecii*. No cross-reaction could be detected for the following species (see Tab. 14, \* detected with sequence alignment):

**Tab. 14:** Cross-reactivity testing










<i>Acinetobacter baumannii</i>	-	<i>Corynebacterium diphtheriae</i>	-	<i>Klebsiella pneumoniae</i>	-	<i>Rhizomucor pusillus</i>	-
<i>Alternaria alternata</i>	-	Coxsackie B4, human	-	<i>Lactobacillus plantarum</i>	-	<i>Saccharomyces cerevisiae</i>	-
Adenovirus 1, human, strain Adenoid 71	-	Cytomegalovirus, human	-	<i>Legionella pneumophila</i> subsp. <i>pneumophila</i>	-	<i>Scedosporium apiospermum</i>	-
Adenovirus 7, human, strain Gomen	-	<i>Doratomyces microsporus</i>	-	<i>Moraxella catarrhalis</i>	-	<i>Scedosporium prolificans</i>	-
<i>Aspergillus fumigatus</i> *	-	Epstein-Barr-Virus, strain B95-8	-	<i>Mycoplasma pneumoniae</i>	-	<i>Sporothrix schenckii</i>	-
<i>Aspergillus terreus</i>	-	<i>Fusarium solani</i>	-	<i>Neisseria meningitidis</i>	-	<i>Staphylococcus aureus</i>	-
<i>Bordetella parapertussis</i>	-	<i>Haemophilus influenzae</i> Rd	-	Parainfluenza virus 1, human strain C35	-	<i>Staphylococcus epidermidis</i>	-
<i>Bordetella pertussis</i>	-	<i>Helicobacter felis</i>	-	Parainfluenza virus 2, human, strain Greer	-	<i>Staphylococcus haemolyticus</i>	-
<i>Candida glabrata</i>	-	Herpes simplex virus 1, strain McIntyre	-	Parainfluenza virus serotype 3	-	<i>Staphylococcus hominis</i> subsp. <i>novobiosepticus</i> R22	-
<i>Candida parapsilosis</i>	-	Herpes simplex virus 2, strain MS	-	Parainfluenza virus 4b, human, strain CH19503	-	<i>Streptococcus pneumoniae</i>	-
<i>Candida tropicalis</i>	-	Human <i>Metapneumovirus</i>	-	Respiratory syncytial virus, human, strain 9320	-	<i>Streptococcus pyogenes</i>	-
<i>Cladosporium</i> spp	-	Influenza virus, infectious A/PR/8/34	-	Respiratory syncytial virus, human, strain Long	-	<i>Streptococcus salivarius</i>	-
Coronavirus 229E, human	-	<i>Klebsiella oxytoca</i>	-	Rhinovirus, genogroup A, human	-	Varicella Zoster Virus (Type B)	-

## 14. Version history

Version number	Chapter and designation
2019-07-22	4. Reagents provided 6. Additional necessary reagents and necessary equipment 9.2 Preparation of the PCR-Mix 9.3 PCR instrument set-up 10. Quality control

## 15. Explanation of symbols

### General symbols

	For in vitro diagnostic use
	Consult instructions for use
	Lot number
	Expiry
	Store at
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

### Testspecific symbols

Not applicable

## 16. Literature

1. Linssen CF *et al.* Inter-laboratory comparison of three different real-time PCR assays for the detection of *Pneumocystis jirovecii* in bronchoalveolar lavage fluid samples 2006, 55: 1229-1235.
2. Tia T *et al.* A highly sensitive novel PCR assay for the detection of *Pneumocystis jirovecii* DNA in bronchoalveolar lavage specimens from immunocompromised patients 2012, 18: 598-603.
3. Borde JP *et al.* Aktuelle Diagnostik und Therapie der *Pneumocystis-jirovecii*-Pneumonie. Dtsch Med Wochenschr 2011, 136: 1426-1430.
4. Centers for Disease Control and Prevention. *Pneumocystis pneumonia* Statistics 2012.
5. Krajicek BJ *et al.* *Pneumocystis pneumonia*: current concepts in pathogenesis, diagnosis, and treatment. Clin Chest Med 2009, 30: 265-278.