RIDASCREEN® Spec. IgG Foodscreen

REF A8020 RIDASCREEN® Spec. IgG Foodscreen Reagents

REF	A0301	RIDASCREEN® Foodscreen Allergens Customized Plate
	A8010	RIDASCREEN® Spec. IgG Foodscreen Reagents Bulk
	A8101	RIDASCREEN® Spec. IgG Foodscreen Customized Plate
	A8121-05/-15	RIDASCREEN® Spec. IgG Foodscreen Plate 1
	A8221-05/-15	RIDASCREEN® Spec. IgG Foodscreen Plate 2
	A8321-05/-15	RIDASCREEN® Spec. IgG Foodscreen Plate 3
	A8621-05/-15	RIDASCREEN® Spec. IgG Foodscreen Plate 6
	A8821-15	RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 1
	A8822-15	RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 2
	A8823	RIDASCREEN® Spec. IgG Foodscreen Plate ES
	A8824	RIDASCREEN® Spec. IgG Foodscreen Plate GS
	A8825	RIDASCREEN® Spec. IgG Foodscreen Plate GB
	A8826	RIDASCREEN® Spec. IgG Foodscreen Plate GP1
	A8827	RIDASCREEN® Spec. IgG Foodscreen Plate GP2
	A8828	RIDASCREEN® Spec. IgG Foodscreen Plate GK
	A8829-05/-15	RIDASCREEN® Spec. IgG Foodscreen Plate VEG
	A8921	RIDASCREEN® Spec. IgG Foodscreen Prescreen 1
	A8931	RIDASCREEN® Spec. IgG Foodscreen Prescreen Customized

Deutsch	3
English.....	21
Čeština	38
Dansk	55
Español.....	72
Ελληνικά	90
Français.....	109
Italiano	127
Nederlands	144
Norsk	162
Polski	179
Português	197
Suomi	215
Svenska.....	232

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Tab. 1: Inhalt von RIDASCREEN® Foodscreen Allergens Customized Plate (A0301)

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Tab. 4: Inhalt von RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101)

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








16. Literatur

Bentz S *et al.* (2010) Clinical relevance of IgG antibodies against food antigens in Crohn's disease: a double-blind cross-over diet intervention study. *Digestion* 81(4):252-64.

Alpay K *et al.* (2010) Diet restriction in migraine, based on IgG against foods: a clinical double-blind, randomised, cross-over trial. *Cephalalgia*. 30(7):829-37.

Aydinlar EI *et al.* (2013) IgG-based elimination diet in migraine plus irritable bowel syndrome. *Headache*. 53(3):514-25.

Uzunismail H *et al.* (2012) The effects of provocation by foods with raised IgG antibodies and additives on the course of Crohn's disease: a pilot study. *Turk J Gastroenterol*. 23:19-27.

G. E. Mullin *et al.* (2010) Testing for Food Reactions: The Good, the Bad, and the Ugly. *Nutr Clin Pract* 25(2):192-198.

A. Aljada *et al.* (2004) Increase in intranuclear nuclear factor B and decrease in inhibitor B in mononuclear cells after a mixed meal: evidence for a proinflammatory effect. *Am J Clin Nutr* 79:682–90.

A. R. Gaby. (1998) The Role of Hidden Food Allergy/Intolerance in Chronic Disease. *Alternative Medicine Review*, 3(2):90-100.

RIDASCREEN® Spec. IgG Foodscreen

REF A0301, A8010, A8020, A8101, A8121-05/-15, A8221-05/-15, A8321-05/-15, A8621-05/-15, A8821-15, A8822-15, A8823, A8824, A8825, A8826, A8827, A8828, A8829-05/-15, A8921, A8931

1. Intended use

For *in vitro* diagnostic use. The RIDASCREEN® Spec. IgG Foodscreen test is a quantitative or qualitative enzyme immunoassay (EIA) for the determination of specific IgG antibodies against food allergens. The test is designed for manual processing and is a combination of reagents and the associated microtiter plate variants. Only the Foodscreen product components shall be combined with each other. The test is suited for the use of human serum and human capillary blood samples. The test is a diagnostic aid in cases of suspected IgG-mediated food hypersensitivities.

2. Summary and explanation of the test

Food intolerance can be identified by the detection of specific IgG antibodies. Food intolerances can be the result of increased permeability of the intestinal wall for a specific food. Immune complexes consisting of IgG antibodies and food can form resulting from food or its ingredients penetrating through the intestinal wall into the surrounding tissue, and these complexes can trigger inflammatory responses that can lead to different diseases.

The test is an EIA. All microtiter plates listed under 4 "Reagents provided" were validated in combination with R-Biopharm AG's A8020 reagent kit. The RIDASCREEN® Spec. IgG Foodscreen reagents may not be used in combination with microtiter plates from other manufacturers, and the RIDASCREEN® Spec. IgG Foodscreen microtiter plates may not be used with reagents from other manufacturers.

3. Test principle

Food extracts that can cause an intolerance reaction are adsorbed on the inner surfaces of the wells of a microtiter plate. Human patient samples (serum or eluate from capillary blood), standard sera, and control sera are pipetted into the appropriate cavities and incubated at room temperature. During incubation, specific IgG antibodies bind to the corresponding adsorbed food antigens. Unbound material is removed by washing. An anti-human IgG antibody conjugated with alkaline phosphatase is then added. During another incubation phase, this antibody conjugate binds to the human IgG antibodies of the patient sample. Any unbound conjugate is removed by washing. Lastly, after being added, the substrate is dephosphorylated by alkaline phosphatase and turns into a yellow-colored product. The intensity of the yellow coloration is proportional to the amount of antigen-specific antibodies in the serum and can be measured by photometry at 405 nm minus the reference wavelength measurement at 620/630 nm. The standard curve is used to quantify the measurement results. (Important: Plates A0301, A8921, and A8931 allow only for a qualitative evaluation.)

4. Reagents provided

Table 1: Contents of RIDASCREEN® Foodscreen Allergens Customized Plate (A0301)

Plate	15 x 1	Microtiter plate, customer-specific, coated with various food antigens (see enclosed allergen matrix), without standard curve, qualitative
-------	--------	--

Table 2: Contents of RIDASCREEN® Spec. IgG Foodscreen Reagents Bulk (A8010)

Wash buffer salt Tween	5 x	Wash buffer salt for 5 liters of wash buffer, 10 mM PBS, 0.05% TWEEN20
Standard 1	30 ml	Standard serum 1, diluted human serum; concentration: 2.5 µg IgG/ml; in stabilized protein solution
Standard 2	30 ml	Standard serum 2, diluted human serum; concentration: 10.0 µg IgG/ml; in stabilized protein solution
Standard 3	30 ml	Standard serum 3, diluted human serum; concentration: 40.0 µg IgG/ml; in stabilized protein solution
Standard 4	30 ml	Standard serum 4, diluted human serum; concentration: 200.0 µg IgG/ml; in stabilized protein solution
Control high	30 ml	Control serum, diluted human serum; see certificate for concentration; in stabilized protein solution
Conjugate IgG	1 l	Anti-human IgG conjugate (sheep); ready to use; alkaline-phosphatase-conjugated antibody in stabilized protein solution
AllergySub	1 l	Substrate solution, ready to use, contains pNPP
AllergyStop R	1 l	Stop reagent, ready to use

Table 3: Contents of RIDASCREEN® Spec. IgG Foodscreen Reagents (A8020)

Sample Buffer IgG	120 ml	Sample dilution buffer, ready to use
Wash buffer salt Tween	6 x	Wash buffer salt for 1 liters of wash buffer, 10 mM PBS, 0.05% TWEEN20
Standard 1	2 ml	Standard serum 1, diluted human serum; concentration: 2.5 µg IgG/ml; in stabilized protein solution
Standard 2	2 ml	Standard serum 2, diluted human serum; concentration: 10.0 µg IgG/ml; in stabilized protein solution
Standard 3	2 ml	Standard serum 3, diluted human serum; concentration: 40.0 µg IgG/ml; in stabilized protein solution
Standard 4	2 ml	Standard serum 4, diluted human serum; concentration: 200.0 µg IgG/ml; in stabilized protein solution
Control high	2 ml	Control serum, diluted human serum; see certificate for concentration; in stabilized protein solution
Conjugate	80 ml	Anti-human IgG conjugate (sheep); ready to use; alkaline-phosphatase-conjugated antibody in stabilized protein solution
AllergySub	80 ml	Substrate solution, ready to use, contains pNPP
AllergyStop R	80 ml	Stop reagent, ready to use

Table 4: Contents of RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101)

Plate	15 x 1	Microtiter plate, customer-specific, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	---

Table 5: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate 1 (A8121-05/15)

Plate 1	1 x 5, 3 x 5	Microtiter plate 1, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
---------	-----------------	--

Table 6: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate 2 (A8221-05/15)

Plate 2	1 x 5, 3 x 5	Microtiter plate 2, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
---------	-----------------	--

Table 7: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate 3 (A8321-05/15)

Plate 3	1 x 5, 3 x 5	Microtiter plate 3, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
---------	-----------------	--

Table 8: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate 6 (A8621-05/15)

Plate 6	1 x 5, 3 x 5	Microtiter plate 6, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
---------	-----------------	--

Table 9: Contents of RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 1 (A8821-15)

Plate	15 x 1	Microtiter plate panel 1, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	--

Table 10: Contents of RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 2 (A8822-15)

Plate	15 x 1	Microtiter plate panel 2, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	--

Table 11: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate ES (A8823)

Plate	15 x 1	Microtiter plate ES, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	---

Table 12: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate GS (A8824)

Plate	15 x 1	Microtiter plate GS, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	---

Table 13: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate GB (A8825)

Plate	15 x 1	Microtiter plate GB, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	---

Table 14: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate GP1 (A8826)

Plate	15 x 1	Microtiter plate GP1, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	--

Table 15: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate GP2 (A8827)

Plate	15 x 1	Microtiter plate GP2, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	--

Table 16: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate GK (A8828)

Plate	15 x 1	Microtiter plate GK, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-------	--------	---

Table 17: Contents of RIDASCREEN® Spec. IgG Foodscreen Plate VEG (A8829-05/15)

Plate VEG	1 x 5, 3 x 5	Microtiter plate VEG, coated with various food antigens (see enclosed allergen matrix), including standard curve, quantitative
-----------	-----------------	--

Table 18: Contents of RIDASCREEN® Spec. IgG Foodscreen Prescreen 1 (A8921)

Plate	15 x 1	Microtiter plate Prescreen 1, coated with various food antigens (see enclosed allergen matrix), including cut-off, qualitative
Cut-Off	5 ml	Cut-off control, diluted human serum; concentration: 10 µg IgG/ml; in stabilized protein solution

Table 19: Contents of RIDASCREEN® Spec. IgG Foodscreen Prescreen Customized Plate (A8931)

Plate	15 x 1	Microtiter plate Prescreen, customer-specific, coated with various food antigens (see enclosed allergen matrix), including cut-off, qualitative
Cut-Off	5 ml	Cut-off control, diluted human serum; in stabilized protein solution

Hazardous materials are indicated according to labeling obligations. For further details, see the safety data sheets (SDSs) at www.r-biopharm.com.

5. Storage instructions

The test kit must be stored at 2 to 8°C and can be used until the expiration date printed on the label. The diluted wash buffer **Wash buffer salt Tween** has a maximum shelf life of 4 weeks when stored at 2 to 8°C. Microbial contamination must be prevented. After the expiration date, the quality guarantee is no longer valid. Contamination of the substrate solution **AllergySub** with the conjugate **Conjugate** must be avoided as this will cause the substrate **AllergySub** to discolor. Likewise, direct exposure of the substrate **AllergySub** to light must be avoided to prevent degradation or discoloration by hydrolysis. Do not use the substrate **AllergySub** if it has developed an intense yellow coloration. See 10.3 “Quality control – indication of instability or expiration of reagents.”

5.1 Storage of microtiter plates

Since the microtiter plates [Plate] are susceptible to air humidity, the aluminum bags should be securely resealed after the required number of plates [Plate] has been removed. The desiccant bag needs to stay in the aluminum bag throughout the storage period. The microtiter plates [Plate] can be stored in the correctly resealed aluminum bag for the specified shelf life at 2 to 8°C.

6. Reagents required but not provided

6.1 Reagents provided

- Distilled or deionized water

6.2 Equipment

- Photometer for microtiter plates (with 405-nm and 620/630-nm filter), measuring range 0 to 3.5 OD
- Microtiter plate washer for 8 to 96 wells
- 12-channel pipette
- Multi-step pipette
- Horizontal shaker (for the elution of capillary blood samples)

7. Warnings and precautions for the users

- For *in vitro* diagnostic use only
- This test must be carried out only by trained laboratory personnel. The guidelines for working in medical laboratories must be followed. Always adhere strictly to the user instructions for carrying out this test.
- Do not pipette samples or reagents by mouth; avoid contact with injured skin and mucous membranes. Wear disposable gloves when handling samples and wash your hands after the test.
- Do not smoke, eat, or drink in areas where samples or test reagents are being processed.
- The standard sera in the kit were tested for HIV and HCV antibodies, HBsAg, and syphilis as per CFR 21.640 and were negative. Nevertheless, the patient samples and all materials that you touch should be treated as potentially infectious. Follow the national safety regulations for handling and disposing of potentially infectious material.
- Standard sera [Standard], conjugate [Conjugate], and wash buffer concentrate [Wash buffer salt Tween] contain sodium azide as a preservative. This substance must not be allowed to come into contact with skin or mucous membranes. Contact with lead or copper pipes may cause the formation of explosive metal azides.

- Do not allow the stop reagent or substrate solution to come into contact with skin, eyes, or clothing. In case of skin contact, wash with soap and water immediately. In case of eye contact, immediately flush under running water for 15 minutes with the eyelids open. Consult an eye specialist. If swallowed, drink plenty of water, avoid vomiting, and seek immediate medical attention.
- All components of the kit must be appropriately disposed of by the user after use (see safety datasheets for the components in the kit).
- All reagents and materials that come into contact with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121°C for at least one hour.

8. Collection and storage of samples

RIDASCREEN® Spec. IgG Foodscreen was developed for testing human serum and capillary blood.

During venous blood collection, the serum should be separated from the clot as quickly as possible after complete coagulation to prevent hemolysis. The serum samples must be stored cold (2 to 8°C) or frozen (–20°C) until they are tested. Avoid repeated freezing and thawing of the serum as well as bacterial contamination. The use of heat-inactivated, lipemic, hemolytic, icteric, or turbid sera may lead to false results.

When capillary blood samples are used in the RIDASCREEN® Spec. IgG Foodscreen test, the right equipment must be used (A8025/A8025-IMU/A8025-PIM

RIDASCREEN® Foodscreen Blood Collection Kit, A8025-BCC RIDASCREEN®

Foodscreen Blood Collection Card). Only this equipment was validated for this test.


The dried capillary blood samples on the Blood Collection Card (component of A8025 and A8025-BCC) are stable for a maximum of 3 weeks when stored in a dry location at room temperature (20 to 25°C). Do not store capillary blood samples in the refrigerator.

Table 20: Stability of sample types

Undiluted serum		Diluted serum	Dried capillary blood sample	Eluate of capillary blood sample
2 to 8°C	–20°C	20 to 25°C	20 to 25°C	20 to 25°C
1 month	> 3 months	Max. 6 hours	3 weeks	Max. 6 hours

Important note for processing capillary blood samples:

The circles on the blood card need to be completely filled with blood. Also check the back of the card to ensure the circles are soaked through.

9. Test procedure

9.1 General information

All reagents, patient sera, and the microtiter plate [Plate] must be brought to room temperature before use (important). Mix the reagents well immediately before use. After use, the reagent kit must be stored at 2 to 8°C.

The coated microtiter plates [Plate] may not be used more than once. Do not use reagents or the microtiter plates [Plate] if the packaging is damaged or the containers are not tightly sealed.

The coated microtiter plates [Plate] can be used with the reagents of different kit lots.

Any deviation from the specified incubation times and temperatures will cause the standard curve to shift in comparison with the certificate. Significant differences in the values of the standard curve can lead to invalid test results. The customer needs to modify (and validate) the incubation time and/or temperature, for example, for fully automated processing of the test.

Do not conduct the test in direct sunlight. It is recommended to cover the microtiter plates [Plate] during the incubation phases.

9.2 Preparing the wash buffer

Dissolve the contents of one package of wash buffer salt [Wash buffer salt] [Tween] in 1 liter of distilled water. Use only the amount needed for each assay. Diluted wash buffer [Wash buffer salt Tween] has only a limited maximum shelf life of 4 weeks (2 to 8°C).

9.3 Preparing the microtiter plates

The coated microtiter plates [Plate] cannot be used directly. Before pipetting takes place, the plates [Plate] need to be labeled with patient name or number.

9.4 Sample dilution

9.4.1 Sample dilution; serum samples

Dilute serum sample with sample dilution buffer [Sample Buffer IgG] 1:100.

Table 21: Recommended sample dilution preparations

Up to 50 determinations		Up to 100 determinations		Up to 200 determinations		Up to 300 determinations	
Serum	Buffer	Serum	Buffer	Serum	Buffer	Serum	Buffer
40 µl	4.0 ml	70 µl	7.0 ml	150 µl	15.0 ml	200 µl	20.0 ml

9.4.2 Sample elution; capillary blood samples

Shake the dried capillary blood sample with the sample dilution buffer [Sample Buffer IgG] in an appropriate vial for 1 (to 14) hour(s) at 20 to 25°C on a horizontal shaker. For each assay, punch out the corresponding perforated circles [Card]. Disposable gloves should be worn to punch out the dried, blood-filled circles [Card] using your finger. You can also use an appropriate sterile object, such as the end of a pipette tip or a laboratory spatula, to punch out the circles [Card].

Table 22: Recommended sample elution preparations

Up to 50 determinations		Up to 100 determinations		Up to 200 determinations		Up to 300 determinations	
Sample	Buffer	Sample	Buffer	Sample	Buffer	Sample	Buffer
2 circles	4.0 ml	4 circles	8.0 ml	6 circles	12.0 ml	8 circles	16.0 ml

The eluate of the capillary blood sample already corresponds to the diluted patient sample and is pipetted onto the microtiter plate [Plate] according to the instructions.

9.5 First incubation (sample incubation)

Pipette 50 µl of the standard sera [Standard 1], [Standard 2], [Standard 3], [Standard 4], the control [Control high], or [Cut-Off] and the diluted patient serum or eluate of the capillary blood sample into the wells of the microtiter plate [Plate] according to the pipetting scheme. Next, incubate the plates [Plate] for 90 minutes at 20 to 25°C. The plates can be stacked in the order of pipetting; please cover the top plate.

9.6 Washing

Wash the wells three times in a microplate washer using 350 µl of diluted wash buffer [Wash buffer salt Tween] for each well. Final siphoning should be quantitative, if possible. The plates [Plate] can be tapped out onto absorbable material after washing to remove any residual wash solution.

9.7 Second incubation (conjugate incubation)

Pipette 50 µl of conjugate [Conjugate] into each well. Next, incubate the plates [Plate] for 90 minutes at 20 to 25°C. The plates [Plate] can be stacked in the order of pipetting; please cover the top [Plate].

9.8 Washing

Wash as described in section 9.6.

9.9 Third incubation (substrate incubation)

Pipette 50 µl substrate solution [AllergySub] into each well. Next, incubate the plates [Plate] for 60 minutes at 20 to 25°C. Keep the incubated plates [Plate] out of direct sunlight. The plates [Plate] can be stacked in the order of pipetting; please cover the top plate [Plate].

9.10 Stopping the reaction and performing the measurement

Adding 50 µl of stop reagent [AllergyStop R] into each well will stop the reaction. After the stop reagent [AllergyStop R] is added, gently “tap” the plate [Plate] to achieve thorough mixing within the wells. Measure the extinction at 405 nm against the reference wavelength of 620/630 nm. The plates [Plate] can be remeasured within 24 hours if they are covered and refrigerated.

10. Quality control – Indication of instability or expiration of reagents

10.1 Evaluation of RIDASCREEN® Spec. IgG Foodscreen: A8101, A8121-05/-15, A8221-05/-15, A8321-05/-15, A8621-05/-15, A8823, A8824, A8825, A8826, A8827, A8828, A8821-15, A8822-15, A8829-05/-15

During each test procedure, each plate [Plate] should contain the standards [Standard 1], [Standard 2], [Standard 3], [Standard 4], and the control [Control high]. The test was performed correctly if the following criteria are met (all information is based on a differential measurement of 405 nm to 620 nm (or 630 nm) [OD_{405/620}]):

[Standard 2]	minus	[Standard 1]	≥ 0.08
[Standard 3]	minus	[Standard 2]	≥ 0.30
[Standard 4]	minus	[Standard 3]	≥ 0.70
[Standard 4]			> 1.0

10.2 Evaluation of RIDASCREEN® Spec. IgG Foodscreen: A0301, A8921, A8931

During each test procedure, each plate Plate should contain the cut-off control Cut-Off. The test was done correctly if the following criterion is met (information based on a differential measurement of 405 nm to 620 nm, or 630 nm [OD_{405/620}]):

$$\text{Cut-Off} \geq 0.20$$

10.3 Indication of instability or deterioration of reagents

A deviation from the expected reference value, turbidity of the reagents, or intense yellow coloration of the substrate AllergySub prior to addition to the wells can be an indication of an expired reagent.

If the specified reference values are not met, the following things must be checked before the test is repeated:

- Expiration date of the reagents used
- Functional performance of the equipment used (e.g., calibration)
- Correct test procedure
- Visual inspection of the kit components for contamination or leaks

Do not use substrate solution AllergySub that has an intense yellowish coloration. If the stipulated values are not obtained after repeating the test, please contact the manufacturer.

11. Evaluation and interpretation

11.1 Serum and capillary blood findings

11.1.1 Basis for the calculations

Calculations must be based on a standard curve for the quantitative evaluation and interpretation of the test. For such, the standards must be run on every plate. The extinction values of the standard are plotted as a function of the related concentrations (µg/ml) semi-logarithmically (y-ln-/x-log graph) in a point-to-point graph to generate the standard curve. In duplicate analyses of the standard, the mean extinctions from the duplicate analyses of the standards will be used to generate the standard curve. The concentrations of the specific IgG antibodies in µg/ml can be determined on the basis of the standard curve from the measured OD values and then converted into IgG classes (see Tables 23 and 24). The evaluation can also be done using appropriate software. The qualitative test evaluation and interpretation is explained in the table.

The standard curve of RIDASCREEN® Spec. IgG Foodscreen was calibrated to an international reference preparation: 1st WHO IRP 67/86 for human IgG.

11.1.2 Concentrations, IgG classes, and calculations for RIDASCREEN® Spec. IgG Foodscreen: A0301, A8101, A8121-05/-15, A8221-05/-15, A8321-05/-15, A8621-05/-15, A8821-15, A8822-15, A8823, A8824, A8825, A8826, A8827, A8828, A8829-05/-15, A8921, A8931

Table 23: Quantitative relationship between the determined IgG concentrations [$\mu\text{g/ml}$], IgG classes, and the patient's antigen-specific IgG concentration from a serum sample

IgG concentration [$\mu\text{g/ml}$]	IgG class	Antigen-specific IgG content	Alternative IgG class	Antigen-specific IgG content
< 7.5	0	Negative	0	Negative
7.5 to < 12.5	1	Mild	1	Elevated
12.5 to < 20.0	2	Elevated		
20.0 to < 50.0	3	High	2	Highly elevated
≥ 50.0	4	Very high		

Table 24: Quantitative relationship between the determined concentrations [$\mu\text{g/ml}$], IgG classes, and the patient's antigen-specific IgG concentration from a capillary blood sample

IgG concentration [$\mu\text{g/ml}$]	IgG class	Antigen-specific IgG content	Alternative IgG class	Antigen-specific IgG content
< 8.0	0	Negative	0	Negative
8.0 to < 18.0	1	Mild	1	Elevated
18.0 to < 40.0	2	Elevated	2	Highly elevated
≥ 40.0	3	High		

These values are rough guidelines only. International standards for specific IgGs to foods do not exist.

Table 25: Qualitative evaluation and interpretation of RIDASCREEN® Foodscreen Allergens Customized Plate (A0301) and RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101), Prescreen 1 (A8921), and Prescreen Customized Plate (A8931) from a serum or capillary blood sample

Ratio	Evaluation	Interpretation
≤ 1.0	-	Negative
> 1.0 to < 2.5	+	Mild
≥ 2.5 to < 5.0	++	Elevated
≥ 5.0	+++	High

Ratio = $OD_{405/620} \text{ (sample)} \div OD_{405/620} \text{ (cut-off)}$. The test is valid if $OD_{405/620}$ of cut-off is ≥ 0.20 .

12. Limitations of the method

The IgG concentrations determined using this test system provide information about the level of sensitization of the patient with regard to the investigated food antigens or mixes.

A correlation between the level of a determined IgG concentration and the occurrence or the severity of clinical symptoms cannot be deduced from this. The results which are obtained must always be interpreted in combination with the complete clinical symptoms.

False-positive test results may occur due to cross-reactivity of the tested antigen with epitopes of other antigens.

The manufacture of food extracts and the coating of microtiter plates do not always eliminate antigenic epitopes, which can lead to false-negative results. IgG antibodies against food antigens that do not appear until during industrial preparation, food preparation, or digestion may not be detectable because they are not present in the original food for which the patient is tested.

13. Performance characteristics

13.1. Precision

The reproducibility of the RIDASCREEN® Spec. IgG Foodscreen test was determined using reference samples covering the entire measuring range from negative to high positive. The mean values and the coefficients of variation (CV) were determined for three reagent lots and for two microtiter plates. For the determination inter-assay reproducibility, reference samples from five sequential working days were assayed in

triplicate, with two runs per day. The measurements were performed by two technicians.

Table 26: Intra-assay precision measurement

Mean CV (%)	8.1
-------------	-----

Table 27: Inter-assay precision measurement

Mean CV (%)	11.1
-------------	------

Table 28: Inter-lot precision measurement

Serum samples	1	2	3
Number of reagent kit lots	3	3	3
Number of microtiter plate lots	2	2	2
Number of allergens	44	44	44
Number of measurements	264	264	264
Coefficient of variation	25.0%	16.9%	15.4%
Mean coefficient of variation	19.1%		

The intra-assay, inter-assay, and inter-lot precision measurements suggest that the RIDASCREEN® Spec. IgG Foodscreen test generates accurate and reproducible results in the complete measuring range under constant as well as changing conditions and over several kit lots.

13.2 Stability

13.2.1 Transport stability

The kit components were tested for stability under temperature fluctuations due to transportation. Two simulated transport scenarios (TC 1 and TC 2) and one genuine overseas transport to Australia (TC 3) were investigated. The temperature fluctuations were documented throughout the transport period. Temperatures between 4.5°C and 46°C were recorded. The kit components were compared with the reference kit, which was stored at the prescribed temperature (2 to 8°C).

No impact on the results after transport of the kit components was detected in any scenario, which means that RIDASCREEN® Spec. IgG Foodscreen is stable under transport conditions.

13.2.2 Kit stability after opening

The kit components were tested for their stability after opening of the individual components for a period of 6 weeks. The analysis was done using two positive sera, which, after 2, 4 and 6 weeks, were compared with the test results from week 0, which correspond to 100% (% of REF). The analysis was also done on a microtiter plate left uncovered after 6 weeks combined with a recently opened reagent kit (Ref. Kit) and a recently opened microtiter plate combined with a reagent kit left uncovered (Ref. Plate). The kits were stored at 2 to 8°C.

No losses of stability were detected in the reagents or microtiter plates over a period of at least 6 weeks. The test results obtained with these kit components showed no significant deviations compared to the references either.

13.3 Analytical sensitivity

Analytical sensitivity is a measurement of the accuracy of a test at low concentrations of the analyte.

The mean LOB of RIDASCREEN® Spec. IgG Foodscreen is 2.8 µg/ml and the mean LOD 4.0 µg/ml.

13.4 Interferents

The following potential substances were tested: hemoglobin, bilirubin, triglycerides, HSA.

Relative impacts on the test result of RIDASCREEN® Spec. IgG Foodscreen were not detected for any of the listed substances.

Keep in mind, however, that combinations of substances were not tested. Therefore, an implausible result in the event of hemolytic, hemorrhagic, or lipemic patient samples must always be investigated specifically in regard to the patient's overall clinical picture.

14. Version history

Version number	Section and designation
2016-08-29	Previous version
2020-02-04	General revision: New layout
2020-02-04	1. Intended use 4. Reagents provided 8. Collection and storage of samples 9. Test procedure 10. Quality control – Indication of instability or expiration of reagents 11. Evaluation and interpretation 13. Performance characteristics Test-specific symbols
2020-11-10	13. Performance characteristics

15. Explanation of symbols

General symbols

	For in vitro diagnostic use
	Consult instructions for use
	Batch number
	Use before
	Storage temperature
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

Test-specific symbols

AllergyStop R	Stop reagent
AllergySub	Substrate solution, ready to use
Card	Blood collection card
Plate	Microtiter plate
Conjugate	Conjugate
Control high	Positive control
Sample Buffer IgG	Sample dilution buffer
Standard 1	Standard 1
Standard 2	Standard 2
Standard 3	Standard 3
Standard 4	Standard 4
Wash buffer salt Tween	Wash buffer salt

16. References

Bentz S *et al.* (2010) Clinical relevance of IgG antibodies against food antigens in Crohn's disease: a double-blind cross-over diet intervention study. *Digestion* 81(4):252-64.

Alpay K *et al.* (2010) Diet restriction in migraine, based on IgG against foods: a clinical double-blind, randomised, cross-over trial. *Cephalalgia*. 30(7):829-37.

Aydinlar EI *et al.* (2013) IgG-based elimination diet in migraine plus irritable bowel syndrome. *Headache*. 53(3):514-25.

Uzunismail H *et al.* (2012) The effects of provocation by foods with raised IgG antibodies and additives on the course of Crohn's disease: a pilot study. *Turk J Gastroenterol*. 23:19-27.

G. E. Mullin *et al.* (2010) Testing for Food Reactions: The Good, the Bad, and the Ugly. *Nutr Clin Pract* 25(2):192-198.

A. Aljada *et al.* (2004) Increase in intranuclear nuclear factor B and decrease in inhibitor B in mononuclear cells after a mixed meal: evidence for a proinflammatory effect. *Am J Clin Nutr* 79:682–90.

A. R. Gaby. (1998) The Role of Hidden Food Allergy/Intolerance in Chronic Disease. *Alternative Medicine Review*, 3(2):90-100.

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4. Dodávaná činidla

Tabulka 1: Obsah RIDASCREEN® Foodscreen Allergens Customized Plate (A0301)

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Tabulka 2: Obsah RIDASCREEN® Spec. IgG Foodscreen Reagents Bulk (A8010)

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Tabulka 3: Obsah RIDASCREEN® Spec. IgG Foodscreen Reagents (A8020)

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Tabulka 4: Obsah RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101)

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Tabulka 5: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate 1 (A8121-05/15)

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Tabulka 6: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate 2 (A8221-05/15)

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Tabulka 7: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate 3 (A8321-05/15)

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Tabulka 8: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate 6 (A8621-05/15)

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Tabulka 11: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate ES (A8823)

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Tabulka 12: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate GS (A8824)

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Tabulka 13: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate GB (A8825)

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Tabulka 14: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate GP1 (A8826)

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Tabulka 16: Obsah RIDASCREEN® Spec. IgG Foodscreen Plate GK (A8828)

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Tabulka 18: Obsah RIDASCREEN® Spec. IgG Foodscreen Prescreen 1 (A8921)

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Tabulka 19: Obsah RIDASCREEN® Spec. IgG Foodscreen Prescreen Customized Plate (A8931)

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Tabell 1: Innholdet i RIDASCREEN® Foodscreen Allergens Customized Plate (A0301)

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Tabell 2: Innhold i RIDASCREEN® Spec. IgG Foodscreen Reagents Bulk (A8010)

Wash buffer salt Tween	5 x	Vaskebuffersalt for 5 liter vaskebuffer, 10 mM PBS, 0,05 % TWEEN20
Standard 1	30 ml	Standardserum 1, fortynnet humant serum; konsentrasjon: 2,5 µg IgG/ml; i stabilisert proteinløsning
Standard 2	30 ml	Standardserum 2, fortynnet humant serum; konsentrasjon: 10,0 µg IgG/ml; i stabilisert proteinløsning
Standard 3	30 ml	Standardserum 3, fortynnet humant serum; konsentrasjon: 40,0 µg IgG/ml; i stabilisert proteinløsning
Standard 4	30 ml	Standardserum 4, fortynnet humant serum; konsentrasjon: 200,0 µg IgG/ml; i stabilisert proteinløsning
Control high	30 ml	Kontrollserum, fortynnet humant serum; se sertifikat for konsentrasjon; i stabilisert proteinløsning
Conjugate IgG	1 l	Anti-humant IgG-konjugat (sau); bruksklart, alkalisk fosfatase-konjugert antistoff i stabilisert proteinløsning
AllergySub	1 l	Subtratløsning, bruksklar, inneholder pNPP
AllergyStop R	1 l	Stoppreagens, bruksklar

Tabell 3: Innhold i RIDASCREEN® Spec. IgG Foodscreen Reagents (A8020)

Sample Buffer IgG	120 ml	Prøvefortynningsbuffer, bruksklar
Wash buffer salt Tween	6 x	Vaskebuffersalt for 1 liter vaskebuffer, 10 mM PBS, 0,05 % TWEEN20
Standard 1	2 ml	Standardserum 1, fortynnet humant serum; konsentrasjon: 2,5 µg IgG/ml; i stabilisert proteinløsning
Standard 2	2 ml	Standardserum 2, fortynnet humant serum; konsentrasjon: 10,0 µg IgG/ml; i stabilisert proteinløsning
Standard 3	2 ml	Standardserum 3, fortynnet humant serum; konsentrasjon: 40,0 µg IgG/ml; i stabilisert proteinløsning
Standard 4	2 ml	Standardserum 4, fortynnet humant serum; konsentrasjon: 200,0 µg IgG/ml; i stabilisert proteinløsning
Control high	2 ml	Kontrollserum, fortynnet humant serum; se sertifikat for konsentrasjon; i stabilisert proteinløsning
Conjugate	80 ml	Anti-humant IgG-konjugat (sau); bruksklart, alkalisk fosfatase-konjugert antistoff i stabilisert proteinløsning
AllergySub	80 ml	Substratløsning, bruksklar, inneholder pNPP
AllergyStop R	80 ml	Stoppreagens, bruksklar

Tabell 4: Innhold i RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101)

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Tabell 5: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate 1 (A8121-05/15)

Plate 1	1 x 5, 3 x 5	Mikrotiterplate 1, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
---------	-----------------	---

Tabell 6: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate 2 (A8221-05/15)

Plate 2	1 x 5, 3 x 5	Mikrotiterplate 2, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
---------	-----------------	---

Tabell 7: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate 3 (A8321-05/15)

Plate 3	1 x 5, 3 x 5	Mikrotiterplate 3, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
---------	-----------------	---

Tabell 8: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate 6 (A8621-05/15)

Plate 6	1 x 5, 3 x 5	Mikrotiterplate 6, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
---------	-----------------	---

Tabell 9: Innhold i RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 1 (A8821-15)

Plate	15 x 1	Mikrotiterplate panel 1, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	---

Tabell 10: Innhold i RIDASCREEN® Spec. IgG Foodscreen Microtiter Plate, Panel 2 (A8822-15)

Plate	15 x 1	Mikrotiterplate panel 2, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	---

Tabell 11: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate ES (A8823)

Plate	15 x 1	Mikrotiterplate ES, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	--

Tabell 12: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate GS (A8824)

Plate	15 x 1	Mikrotiterplate GS, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	--

Tabell 13: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate GB (A8825)

Plate	15 x 1	Mikrotiterplate GB, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	--

Tabell 14: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate GP1 (A8826)

Plate	15 x 1	Mikrotiterplate GP1, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	---

Tabell 15: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate GP2 (A8827)

Plate	15 x 1	Mikrotiterplate GP2, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	---

Tabell 16: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate GK (A8828)

Plate	15 x 1	Mikrotiterplate GK, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-------	--------	--

Tabell 17: Innhold i RIDASCREEN® Spec. IgG Foodscreen Plate VEG (A8829-05/15)

Plate VEG	1 x 5, 3 x 5	Mikrotiterplate VEG, belagt med ulike matantigener (se vedlagt allergenmatrise), inkludert standardkurve, kvantitativ
-----------	-----------------	---

Tabell 18: Innhold i RIDASCREEN® Spec. IgG Foodscreen Prescreen 1 (A8921)

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Tabell 19: Innhold i RIDASCREEN® Spec. IgG Foodscreen Prescreen Customized Plate (A8931)

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Tabela 2: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Reagents Bulk (A8010)

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Tabela 3: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Reagents (A8020)

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Tabela 4: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Customized Plate (A8101)

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Tabela 5: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate 1 (A8121-05/15)

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Tabela 6: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate 2 (A8221-05/15)

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Tabela 7: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate 3 (A8321-05/15)

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Tabela 8: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate 6 (A8621-05/15)

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Tabela 12: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate GS (A8824)

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Tabela 16: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Plate GK (A8828)

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Tabela 19: Conteúdo do RIDASCREEN® Spec. IgG Foodscreen Prescreen Customized Plate (A8931)

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