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Testspezifische Symbole

	RIDASCREEN Foodscreen Blood Collection Card
---	---

RIDASCREEN® Foodscreen Blood Collection Card

REF A8025-BCC

1. Intended use

For *in vitro* diagnostic use. The RIDASCREEN® Foodscreen Blood Collection Kit/Card is used for collecting and transporting a human capillary blood sample taken from the fingertip. For professional use only.

2. Summary and explanation of the test

The RIDASCREEN® Foodscreen Blood Collection Card allows for easy collection of capillary blood samples. The samples can then be sent to a laboratory for further analysis and testing for food intolerances.

3. Test principle

The RIDASCREEN® Foodscreen Blood Collection Card is a cellulose card that makes it possible to collect a specifically defined volume of human capillary blood from the fingertip and to transport the sample to a laboratory for further analysis. In combination with the RIDASCREEN® Spec. IgG Foodscreen and RIDA®CHIP FoodGuide tests, the human capillary blood samples can be analyzed for potential food intolerances.

4. Reagents provided

Table 1: Reagents provided

REF	Kit components	Amount	Description
-	Card	100 x	RIDASCREEN® Foodscreen Blood Collection Card for collecting and transporting a human capillary blood sample

5. Storage instructions

Please follow the handling guidelines in Table 2 and store the kit directly after use according to the information specified. The quality guarantee is no longer valid after the expiration date or the recommended storage duration of the opened or diluted reagents. Microbial contamination must be prevented.

Table 2: Storage conditions and information

	Storage temperature	Maximum storage time	Additional notes on storage
unopened	9 °C - 30 °C	Can be used until the printed expiration date	-

6. Reagents required but not provided

- Alcohol swabs
- Safety lancet
- Bandages

7. Warnings and precautions for the users

- For *in vitro* diagnostic use only.
- This test is intended for professional use only, and only qualified laboratory personnel may perform the test. The guidelines for working in medical laboratories must be followed. Always adhere strictly to the operating manual when carrying out this test.
- For safety reasons, the patient samples and all materials that come into contact with the samples should be treated as potentially infectious. Follow the national safety regulations for handling and disposing of potentially infectious material.
- All materials that come into contact with potentially infectious samples must be treated with suitable disinfectants or autoclaved at 121 °C for at least one hour.
- Anyone helping you take the blood sample should wear disposable gloves (not contained in the kit).
- Users are responsible for the proper disposal of all materials after use. For disposal, please adhere to national regulations.
- The RIDASCREEN® Foodscreen Blood Collection Card should be stored dry in the original packaging between 9 °C and 30 °C.
- The RIDASCREEN® Foodscreen Blood Collection Card containing the blood sample must arrive at the laboratory no later than two weeks after blood collection.
- Do not use open or damaged kits to collect blood.
- Do not use the RIDASCREEN® Foodscreen Blood Collection Card after the expiration date listed on the packaging.

- The RIDASCREEN® Foodscreen Blood Collection Card is intended for single use only.
- R-Biopharm AG expressly rejects any liability for the use or sale of the RIDASCREEN® Foodscreen Blood Collection Card or its components if used for anything other than the intended use specified in the instructions.
- The RIDASCREEN® Foodscreen Blood Collection Card is an accessory for the RIDASCREEN® Spec. IgG Foodscreen and RIDA®CHIP FoodGuide test systems and was validated with them. R-Biopharm AG assumes no liability for use with non-validated test systems.
- For users in the European Union: Report all serious adverse events associated with the product to R-Biopharm AG and the appropriate national authorities.

8. Collection and storage of samples

The dried capillary blood sample is stable at room temperature for up to 5 weeks. Do not store capillary blood samples in the refrigerator.

9. Test procedure

9.1

Remove the RIDASCREEN® Foodscreen Blood Collection Card **Card** from the plastic bag containing the desiccant, and record the patient name, date of blood collection, and lab ID, if applicable, on the card. Please do not use a liquid ink roller pen because the ink can “bleed” through the paper. Do not let the perforated rings come in contact with skin.

9.2


Capillary blood should be collected from the finger in accordance with the general hygiene guidelines for blood collection. Use an appropriate safety lancet for blood collection. Apply blood drops to the card **Card** by touching the tip of the blood drop (not the finger) to the first circle. It may take several drops of blood to fill the circles all the way to the edge. Repeat these steps for the other circles.

9.3

The circles should be completely filled up to the edge with blood, and blood should be soaked through to the back side (see illustration).

Table 3: Number of circles needed

RIDACHIP®FoodGuide	RIDASCREEN® Spec. IgG Foodscreen
2 circles	≤ 100 allergens: 4 circles
	≤ 200 allergens: 6 circles
	≤ 300 allergens: 8 circles

9.4

After applying the blood, allow the card **Card** to dry overnight at room temperature.

9.5

Fold the card **Card**. Place the card in the plastic bag with the desiccant for storage or shipping.

9.6

When the sample is processed in the laboratory, the perforated circles filled with blood will be punched out, the blood will be eluted, and the test will be performed.

10. Quality control – Indication of instability or expiration of reagents

Information on quality control cannot be provided because it pertains only to the RIDASCREEN® Spec. IgG Foodscreen or RIDA®CHIP FoodGuide assays. More details on quality control can be found in the respective instructions for use for the RIDASCREEN® Spec. IgG Foodscreen and the RIDA®CHIP FoodGuide.

11. Evaluation and interpretation

Information on evaluation and interpretation cannot be provided because it pertains only to the RIDASCREEN® Spec. IgG Foodscreen or RIDA®CHIP FoodGuide assays. More details on the evaluation and interpretation can be found in the respective instructions for use for the RIDASCREEN® Spec. IgG Foodscreen and the RIDA®CHIP FoodGuide.

12. Limitations of the method

Information on limitations of the method cannot be provided because it pertains only to the RIDASCREEN® Spec. IgG Foodscreen or RIDA®CHIP FoodGuide assays. More details on the limitations of the method can be found in the respective instructions for use for the RIDASCREEN® Spec. IgG Foodscreen and the RIDA®CHIP FoodGuide.

13. Version history

Version number	Section and designation
2022-03-24	IVDR version

14. Explanation of symbols

General symbols

	For in vitro diagnostic use
	Observe operating manual
	Batch number
	Use before
	Storage temperature
	Item number
	Number of tests
	Date of manufacture
	Manufacturer
	Do not reuse
	Keep dry

Test-specific symbols

	RIDASCREEN Foodscreen Blood Collection Card
---	---

RIDASCREEN® Foodscreen Blood Collection Card

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