

Evaluation of the RIDA[®]GENE Bordetella real-time PCR assay for the laboratory diagnosis of Bordetella infections

L. Kastl¹, A. Simons¹, A. Veldenzer², U. Eigner², M. Holfelder²
¹ R-Biopharm AG, Darmstadt, Germany
² Limbach Laboratory, Department of Microbiology, Heidelberg, Germany

Objectives

Bordetella pertussis is the major cause of pertussis (whooping cough). Pertussis can cause a serious illness in people of all age groups, which can be life-threatening particular in infants. Other species, including *Bordetella holmesii* and *Bordetella parapertussis* cause a milder whooping cough-like illness, however it is estimated that 3 - 35 % of the Bordetella infections are caused by *B. parapertussis* and also *B. holmesii* is frequently detected upon Bordetella infections.^{1,2,3} Transmission of pertussis still occurs frequently because protection from vaccination lasts for 5 - 10 years and protection after natural infection wanes after 10 to 15 years. In addition, pertussis vaccination is described to

lack cross-protection against *Bordetella holmesii*.⁴

After the incubation time, the clinical course of a Bordetella infection can be divided into three stages: catarrhal stage (1 - 2 weeks), paroxysmal stage (1 - 2 weeks) and convalescent stage (6 - 10 weeks).⁵ Whereas culture is only appropriate in the first two weeks and serological diagnosis at earliest 2 weeks after infection, real-time PCR allows rapid and sensitive detection in the first four weeks after cough onset. This study aimed to evaluate the performance of the RIDA[®]GENE Bordetella assay for the direct detection of Bordetella infections. The results were compared to the GenoQuick[®] Bordetella assay (Hain Lifescience, Nehren, Germany).

Methods

The RIDA[®]GENE Bordetella multiplex real-time PCR assay is a qualitative assay using fluorogenic target-specific hydrolysis probes for the differential detection of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii* (Table 1). The GenoQuick[®] Bordetella assay allows direct detection and differentiation of *Bordetella pertussis* and *Bordetella parapertussis*.

212 swab specimens from patients with symptoms of a bacterial respiratory infection were isolated with the NucliSENS[®] easyMag[®] automated extraction platform (bioMérieux). DNA was analysed with the RIDA[®]GENE Bordetella assay on the LightCycler[®] 480II (Roche) and by visual inspection of the GenoQuick test strip according to the manufacturer's instructions (Figure 1).


Table 1 Analyte target gene and detection channel overview

Analyte	Target Gene	Detection Channel Light [®] Cycler 480II
<i>B. pertussis</i> / <i>B. holmesii</i>	IS481	465/510
<i>B. parapertussis</i>	IS1001	533/610
<i>B. holmesii</i>	IS1001	618/660
ICD	Synthetic target DNA sequence	533/580

Picture 1 RIDA[®]GENE Bordetella multiplex real-time PCR



Figure 1 Study design



Results

Of the 212 patient samples, 210 were concordant (99 %). 187 samples were negative for a Bordetella infection.

Both tests identified 22 positive samples for *Bordetella pertussis* (10.43 %) and 1 positive sample for *Bordetella parapertussis* (0.47 %). There was a discrepancy in the result for 2 samples between both

assays (Figure 2). Those samples were weak positive for *B. parapertussis* with the GenoQuick[®] Bordetella assay and were not detected by the RIDA[®]GENE Bordetella multiplex real-time PCR assay. The discrepant samples were sent to the national reference center which stated those 2 samples as negative for *Bordetella spp.*, confirming the result obtained by the RIDA[®]GENE Bordetella multiplex real-time PCR assay.

Figure 2 RIDA[®]GENE Bordetella vs. Geno[®]Quick Bordetella

RIDA [®] GENE Bordetella		GenoQuick [®] Bordetella		Total	
		+	-		
+		23	0	23	Pos. agreement 92 %
-		2*	187	189	Neg. agreement 100 %
	Total	25	187	212	Overall agreement 99 %

*two discrepant samples were sent to the national reference center for Bordetella, Germany. Both samples were tested negative for a Bordetella infection

Figure 3 Example of *B. pertussis*/*B. holmesii* run on the LightCycler[®] 480II




Figure 4 Example of *B. parapertussis* run on the LightCycler[®] 480II





Figure 5 Example of *B. holmesii* run on the LightCycler[®] 480II



Acknowledgment

We thank Dr. Marion Riffelmann and Dr. Nicole Kennerknecht from the German consulting laboratory for Bordetella pertussis, Krefeld for analysis of discrepant results.

References:

- Linneman CC et al. Bordetella parapertussis. Recent experience and a review of the literature. Am J Dis Child 1977, 131 (5): 560-563.
- Mertsola, J. Mixed outbreak of Bordetella pertussis and Bordetella parapertussis infection in Finland. Eur J Clin Microbiol 1985, 4(2): 123-128.
- Njamkepo E et al. Significant Finding of Bordetella holmesii DNA in Nasopharyngeal Samples from French Patients with Suspected Pertussis. J Clin Microbiol 2011, 49(12):4347-4348.
- Zhang X et al. Lack of Cross-protection against Bordetella holmesii after Pertussis Vaccination. Emerg Infect Dis 2012, 18(11): 1771-1779.
- Centers for Disease Control and Prevention 2012. Pertussis. Epidemiology and Prevention of Vaccine-Preventable Diseases. The Pink Book: Course Textbook - 12th Edition Second Printing (May 2012). Accessed 21.05.2013.

Conclusion

The RIDA[®]GENE Bordetella assay is a highly sensitive and specific multiplex real-time PCR assay for diagnosis of *Bordetella pertussis*, *Bordetella parapertussis* and *Bordetella holmesii*.

The RIDA[®]GENE Bordetella multiplex real-time PCR assay allows detection and differentiation of *B. pertussis*, *B. parapertussis* and *B. holmesii*.

The real-time PCR format reduces the time to result to 2 hours. The RIDA[®]GENE Bordetella real-time PCR kit contains an internal control DNA that can be used as an amplification control or additionally as an extraction control.