RIDA® LINE Helicobacter IgA, IgG

Art. Nr.: LB2313 (IgA)
LB2323 (IgG)

R-Biopharm AG, An der neuen Bergstraße 17, D-64297 Darmstadt, Germany
Tel.: +49 (0) 61 51 81 02-0 / Telefax: +49 (0) 61 51 81 02-20

1. Allgemeines

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Literatur

1. Marshall, B. F., Royce, H., Annear, D. I., Goodwin, C. S., Taylor, N. S., Edmonds, P., Sly, L. I., Brenner, D. J.: Original Isolation of *Campylobacter pyloridis* from human gastric mucosa. *Microbios Letter* 25, 83 (1982)
2. Warren, J. R.: Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* I, 1273 - 1275 (1983)
3. Parsonnet, J., Friedman, G. D., Vanderstetten, D. P., Chang, Y., Vogelmann, J. H., Orentreich, N., Sibley, R. K.: *Helicobacter pylori* infection and the risk of gastric carcinoma. *N. Engl. J. Med.* 325, 1127 - 1131 (1991)
4. Stolte, M., Eidt, S.: Healing gastric MALT lymphomas by eradicating *H. pylori*. *Lancet* 342, 568 (1993)
5. Stolte, M.: *Helicobacter pylori* and gastric MALT lymphoma. *Lancet* 339, 745 - 746 (1992)
6. Hentschel, E., Brandstatter, G., Dragosics, B., Hirschl, A. M., Nemec, H., Schutze, K., Taufer, M., Wurzer, H.: Effect of Ranitidine and Amoxicillin plus Metronidazole on the eradication of *Helicobacter pylori* and the recurrence of duodenal ulcer [see comments]. *N. Engl. J. Med.* 328, 308 - 312 (1993)
7. Stadelmann, O.: *Helicobacter pylori*: Indikationen und Praxis der Therapie. *Deutsches Ärzteblatt* 92, 39 - 41 (1995)
8. Telford, J. L., Ghiara, P., Dellorco, M., Comanducci, M., Burroni, D., Bugnoli, M., Tecce, M. F., Censini, S., Covacci, A., Xiang, Z. Y., Papini, E., Montecucco, C., Parente, L., Rappuoli, R.: Gene structure of the *Helicobacter pylori* cytotoxin and evidence of its key role in gastric disease. *J. Exp. Med.* 179, 1653 -1658 (1994)
9. Lee, A.: *H. pylori* initiated ulcerogenesis: look to the host. *Lancet* 341, 281 (1993)
10. Schmitt, W., Haas, R.: Genetic analysis of the *Helicobacter pylori* vacuolating cytotoxin - structural similarities with the IgA protease type of exported protein. *Mol. Microbiol.* 12, 307 - 319 (1994)
11. Covacci, A., Censini, S., Bugnoli, M., Petracca, R., Burroni, D., Macchia, G., Massone, A., Papini, E., Xiang, Z., Figura, N., Rappuoli, R.: Molecular characterization of the 128-kDa immunodominant antigen of *Helicobacter pylori* associated with cytotoxicity and duodenal ulcer. *Proc. Natl. Acad. Sci. USA* 90, 5791 - 5795 (1993)
12. Xiang, Z., Censini, S., Bayeli, P. F., Telford, J. L., Figura, N., Rappuoli, R.: Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect. Immun.* 63, 94 (1995)

13. Krakowka, S., Morgan, D. R., Kraft, W. G., Leunk, R. D.: Establishment of gastric *Campylobacter pylori* infection in the neonatal gnotobiotic piglet. *Infect. Immun.* 55, 2789 - 2796 (1987)
14. Karita, M., Li, Q., Cantero, D., Okita, K.: Establishment of a small animal model for human *Helicobacter pylori* infection using germ-free mice. *Am. J. Gastroenterol.* 89, 208 - 213 (1994)
15. Marchetti, M., Arico, B., Burroni, D., Figura, N., Rappuoli, R., Ghiara, P.: Development of a mouse model of *Helicobacter pylori* infection that mimics human disease. *Science* 267, 1655 - 1658 (1995)
16. Jaskowski, T. D., Martins, T. B., Hill, H. R., Litwin, C. M.: Immunoglobulin A Antibodies to *Helicobacter pylori*. *J. Clin. Microbiol.* 35, 2999 - 3000 (1997)
17. Yunoki N, Yokota K, Mizuno M, Kawahara Y, Adachi M, Okada H, Hayashi S, Hirai Y, Oguma K, Tsuji T: Antibody to Heat Shock Protein Can Be Used for Early Serological Monitoring of *Helicobacter pylori* Eradication Treatment. *Clin Diagn Lab Immun*, July 2006, Vol 19, No.3, p449-490.
18. Annibale B, Lahner E, Santucci A, Vaira D, Pasquali A, Severi C, Mini R, Figura N, Fave GD: CagA and VacA are Immunoblot Markers of Past *Helicobacter pylori* in Atrophic Body Gastritis. *J. Compilation* © 2007 Blackwell Publishing Ltd, *Helicobacter* 12:23-30, *Helicobacter* ISSN 1523-5378
19. Gao L, Weck MN, Michel A, Pawlita M, Brenner H: Association between chronic atrophic gastritis and serum antibodies to 15 *Helicobacter pylori* proteins measured by multiplex serology. *Cancer Res* 2009 Apr 1; 69(7): 2973-80

RIDA® LINE Helicobacter IgA, IgG

1. Intended use

For *in vitro* diagnostic use. The RIDA® LINE Helicobacter IgA, IgG is a qualitative line immunoassay for the detection of IgA and IgG antibodies against Helicobacter pylori in human serum or plasma.

In contrast to ELISA, the test principle allows the identification of specific antibodies against the various antigens of Helicobacter pylori through the separate line-up of different single antigens. By using type-specific antigens it is also possible to differentiate between an infection with type I and type II.

The RIDA® LINE Helicobacter IgA, IgG is a confirmatory test and can be used to clarify unclear screening results.

2. General

First descriptions of spirally-formed organisms found in the stomach of humans date back to the turn of the century. At the beginning of the 80's, Marshall and Warren were successful in cultivating the organism from a patient with chronically active gastritis, and therewith were able to give a detailed description of the organism. The spirally-formed bacteria are found especially in the mucosa of the antrum. At first, the organism was named Campylobacter pyloridis or pyloris, but after the creation of a separate genus Helicobacter, it was changed to Helicobacter pylori.

Helicobacter pylori is a Gram-negative, spirally-formed bacterium with up to three spirals. Like all members of the genus Helicobacter, Helicobacter pylori also possesses flagella (unipolar, 2 - 6 flagella). A positive identification is ensured by the very strong urease reaction as well as the positive detection of catalase and oxydase in the corresponding culture morphology.

Different Helicobacter pylori isolates show a high genetic variability, but their phenotypic characteristics are mostly similar. The only phenotypic differences known at the moment restrict themselves to the production or the non-production of a vacuolating cytotoxin, VacA, and a cytotoxin associated protein, CagA. Due to this difference, clinical isolates are divided into 2 types. Helicobacter pylori strains of type I have the ability to produce these proteins. All other strains which do not have this ability belong to type II. It is worth noting that patients with duodenal ulcers are more often infected with the VacA and the CagA expressing type I, which could mean that type I infections lead to more serious illnesses. These clinical diagnoses are reinforced with the observation that lysates of Helicobacter pylori type I or purified VacA can lead to stomach damages in mice.

So far, it has not been determined how *Helicobacter pylori* is transmitted. It is assumed that the pathway of transmission is either oral-oral or fecal-oral. *Helicobacter pylori* is associated with nearly all duodenal and stomach ulcers, as well as with MALT-lymphoma of the stomach and is also associated with a higher risk for gastric adenocarcinoma. *Helicobacter pylori* has been classified as a carcinogen of group 1 by the WHO. Nearly 50% of the world population is infected, whereby the prevalence in the developed countries rises only slowly with advancing age, in contrast to the developing countries. Usually, in the developing countries, infection occurs during childhood and persists over decades. Even though all *Helicobacter pylori* infected people histologically develop a gastritis this leads to more serious illnesses such as peptic ulcers, gastric carcinoma or MALT-lymphoma in only approx. 10%. Next to the different virulent strains of type I and II, host factors are also discussed as possible causes.

Fundamentally, the diagnostic strategies can be differentiated between invasive and non-invasive methods. The invasive method involves the histological detection of the organism from biopsy material which is derived from the upper gastrointestinal tract. This is linked with the detection of urease activity of the organism from the biopsy material. Culturing of the organism can follow. The gastroscopy is usually necessary to exclude the presence of a carcinoma or for therapeutic reasons (e.g. to stop hemorrhages). However, a gastroscopy alone cannot detect or exclude an *Helicobacter pylori* infection. The biopsy can lead to false negative results, since material can be removed from places of the stomach where *Helicobacter pylori* has not colonized.

With the identification of serum antibodies the non-invasive diagnosis by enzyme immunoassay or western or line blot gains importance in view of its advantages for the patients. Herewith, the enzyme immunoassay is distinguished by its easy handling as a screening method. Additionally, the line blot has the advantage that it differentiates between type I and II *H. pylori*. Due to the separation of the antigens, antibodies directed against them can be identified isolated and conclusions about the presence of VacA and CagA can be drawn. Herewith, an analysis of the different *Helicobacter pylori* types with their clinical significance is made possible.

The eradication of *Helicobacter pylori* with a combined antibiotic therapy can heal the gastritis and can also lead to an accelerated healing of an ulcer. Only when *Helicobacter pylori* as the cause of the ulcer has been eradicated will the healed ulcera remain in remission.

3. Test principle

For this line immunoassay, highly purified recombinant antigens of *Helicobacter pylori* (CagA, VacA, GroEL, UreA, HcpC, gGT) are fixed on nitrocellulose membrane strips.

These strips will be incubated with the samples (serum, plasma). Specific antibodies present in the samples will bind to the pathogen antigens on the strip. After removing the unbound antibodies by washing, a second incubation with anti-human IgA or IgG conjugate follows. In the last step (staining reaction), specifically bound antibodies against the different *Helicobacter pylori* antigens become visible as dark bands on the strips.

4. Reagents provided

Table 1: Pack contents

(there are enough reagents in a pack for 20 determinations)

			LB2313 IgA	LB2323 IgG
Strip	20	2 x 10 test strips (numbered); coated with H. pylori antigen	X	X
Wash/Diluent 10x	100 ml	Wash / dilution buffer; 10x concentrated phosphate buffer; contains NaCl, KCl and detergent, preservative : MIT (0.1%) and Oxypyron (0.2%)	X	X
MP	5 g	Milk powder; from skim milk	X	X
Conjugate IgA 100x <i>colourless lid</i>	500 µl	Anti-human IgA conjugate (rabbit); 100x concentrated peroxidase conjugated antibodies; contains < 0.1 % NaN ₃ , MIT (<0.1%), chloracetamide (<0.1%)	X	
Conjugate IgG 100x <i>green lid</i>	500 µl	Anti-human IgG conjugate (rabbit); 100x concentrated peroxidase conjugated antibodies; contains < 0.1 % NaN ₃ , MIT (<0.1%), chloracetamide (<0.1%)		X
Substrate	40 ml	Chromogenic substrate; ready for use; contains tetramethylbenzidine	X	X
Evaluation sheet	1		X	X

5. Reagents required but not provided

5.1. Reagents

- Deionized water

5.2. Accessories

- Incubation trays (please order with Art. No. LB0003 if required)
- Vacuum pump to remove liquid out of the incubation chambers
- Horizontal shaker
- Vortex mixer
- Micropipettes for volumes 20 µl and 1000 µl
- 10 ml pipette
- Measuring cylinder, 50 ml and 1000 ml
- plastic tweezers to handle the teststrips

6. Warnings and precautions for the users

For *in vitro* diagnostic use only.

The test has to be used only by trained and authorized laboratory personal. Please refer to guidelines for safety regulations in medical laboratories. The test protocol must be followed strictly.

Do not mouth pipette samples or reagents. Avoid contact with injured skin or mucous membranes. Wear disposable gloves while handling samples and wash hands after assay is complete. Do not smoke, eat or drink in areas where specimen or kit reagents are handled.

After the addition of patient or control specimens the strips, the material which was in contact with it as well as the patient samples itself should be considered potentially infectious and be treated with respect to the given national safety precautions.

All reagents and materials coming in contact with potential infectious specimens must be treated with disinfectants or autoclaved at 121 °C for at least one hour.

The conjugates contain sodium azides (NaN₃), MIT (methylisothiazolone), oxypyron, chloracetamide and hydrogen peroxide. Contact with skin or mucous membranes must be avoided. Contact with lead or copper pipes sodium azide can form explosive metal azides.

Cross-contamination of patient samples or conjugates can lead to inaccurate test results.

7. Storage instructions

All reagents have to be stored at 2 – 8 °C (**do not freeze**) and can be used up to the expiry date printed on the labels. Substrate, wash/dilution buffer and conjugate from kits with different lot numbers could be used. Please note the shelf-life of these components.

Microbial contamination has to be avoided.

The concentrated conjugates have to be diluted freshly in required amounts. It is not possible to store the ready-for-use conjugate solution.

Ready-for-use wash/dilution buffer is odorless and easily marred. If stored at 2 – 8 °C it can be used up to four weeks.

The tube containing the test strip should be only opened immediately before use to avoid condensation. Leave unused strips in the tube and continue to store at 2 – 8 °C (reseal tube tightly, test strips must not become moist before the test!).

The test kit components must be protected from direct sunlight throughout the entire test procedure. The substrate solution is especially sensitive to light.

8. Indication of instability or deterioration of reagents

The packages have an expiry date after which a quality warranty cannot be given. Substrate solution must not be used if turned blue.

9. Specimen collection and storage

The RIDA[®] LINE Helicobacter can be carried out with serum or plasma.

The use of heat-inactivated, icteric, haemolysed, lipemic or turbid samples is not recommended. Microbial contamination must be avoided. Repeated freezing and thawing of the samples is not recommended. Insoluble substances must be removed from the sample before incubation. The sample material can be stored for up to 2 weeks at 2 – 8 °C if the test cannot be carried out immediately. A prolonged storage of samples is possible at –20 °C or below.

10. Test procedure

10.1. Preliminary comments

Bring all kit components at room temperature (18 – 25 °C) for at least 30 minutes prior to use. The test procedure is carried out at room temperature.

Remove only the volume of reagents that is needed for test procedure. Do not pour reagents back into vials as reagent contamination may occur.

Make sure that incubation solutions do not flow over into other wells. Mix the concentrated reagents and samples thoroughly before use. Avoid a build up of foam. Protect kit components from direct sunlight throughout the entire test procedure. The substrate is especially sensitive to light.

Microwell strips and reagents must not be used if pouch is damaged or vials are leaking.

The strips do not have to be used in their order of numeration. But it is important to keep a record about which strip is incubated with the corresponding sample.

The reproducibility of Blot results depend very strongly on the consistency and regularity of washing the strips. The washing steps described below should always be adhered to.

10.2. Preparation of the ready to use wash / dilution buffer

The diluted, ready to use wash / dilution buffer is used for washing, sample dilution and preparation of conjugate solutions. The required volume for the corresponding number of tests must be determined (Table 2). Ready-for-use wash/dilution buffer is odorless and easily marred. If stored at 2 – 8 °C it can be used up to four weeks.

The skimmed milk powder is first pre-dissolved in wash/dilution buffer concentrate and deionised water is then added to this mixture to reach the final volume (1 part of concentrated buffer is diluted with 9 parts deionized water).

Table 2: Wash / dilution buffer per applied strip (device-specific dead volume is not considered)

Reagent	Formula	Example for 5 strips
Skimmed milk powder (g)	= number of strips x 0,1	0.5 g
Wash-/dilution buffer -concentrate (ml)	= number of strips x 2	10 ml
Deionised water (ml)	= number of strips x 18	90 ml
Ready-for-use wash-/dilution buffer (ml)	= number of strips x 20	100 ml

10.3. First Incubation

For each determination, use one chamber of the incubation plate and fill it with 2 ml diluted wash / dilution buffer. For each determination, one blot strip is taken out of the vial with a pair of tweezers (touch the strips at the end with the numeration only). Dip the strip into the buffer. The strip must be completely moistened with buffer.

The strip numeration must show upwards.

20 µl of an undiluted sample (serum or plasma) are pipetted onto the test strip for each incubation mixture (dilution 1+100).

It must be taken care that the neighboring chambers are not contaminated since this could lead to false positive results. During the incubation, the incubation plates should be covered up with the lids that are enclosed in the test kit.

The strips with the samples are placed on a horizontal shaker and are incubated for 1 hour at room temperature (18 – 25 °C).

10.4. Washing

Remove the solution (via suction), preferably with an attached disinfecting trap. Wash the strips three times for 5 min with 2 ml each of the diluted wash / dilution buffer. Place the incubation plate on the shaker. Between the washing steps the buffer is sucked out of the chambers.

Attention!

The conjugate solution must be prepared during the last washing step.

10.5. Preparation of the conjugate solution

The conjugate enclosed in the kit is provided as concentrate. The concentration is the same for the IgA as well as for the IgG conjugate. For its application to the test, one part of the conjugate concentrate is diluted with 100 parts of the ready-for-use wash/dilution buffer (1+100).

The required volumes for the corresponding number of strips are listed in table 3. Since the diluted conjugate is stable only for a limited period of time, it must be prepared freshly prior to its application.

Table 3: Preparation of the ready-for-use conjugate solutions

Reagent	Formula	Example for 5 strips
Conjugate concentrate (µl)	= number of strips x 20	100 µl
Ready-for-use wash/dilution buffer (ml)	= number of strips x 2	10 ml

The quantities are calculated without dead volume. Depending on handling (manual or on a device) please mix additional conjugate for 1 to 3 strips.

10.6. Second incubation

Add 2 ml of the corresponding freshly prepared conjugate solution to each strip. Incubate the strips for 45 minutes on a horizontal shaker at room temperature.

During the incubation, the incubation plates should be covered up with the lids that are enclosed in the test kit.

10.7. Washing

Washing sequence: see point 10.4.

10.8. Third incubation

Add 1.5 ml of the ready to use substrate to each strip. Incubate the strips for 8 minutes at room temperature with gentle shaking. Thereafter, the substrate is removed via suction. Still in their chambers, the strips should be rinsed three times with deionized water to stop the color reaction.

The strips are carefully taken out of the chambers with tweezers. To dry the strips, they are placed between filter paper (2 h). Stripes should be protected from sun light. After drying they can be evaluated.

Summary of the test procedure

1. Bring all reagents to room temperature
2. Add 2 ml wash / dilution buffer into each chamber
3. Place one strip into each chamber
4. Add 20 µl sample
5. Incubate for 1 hour at room temperature
6. Wash 3 times with 2 ml wash / dilution buffer; 5 min each time
7. Add 2 ml freshly prepared conjugate solution
8. Incubate for 45 minutes at room temperature
9. Wash 3 times with 2 ml wash / dilution buffer; 5 min each time
10. Add 1.5 ml substrate
11. Incubate for 8 minutes at room temperature
12. Rinse strips with deionized water
13. Dry on filter paper (2h)
14. Analysis

11. Analysis

11.1. Preliminary comments

An analysis of the test can be carried out if the following criteria have been fulfilled:

- Reaction control band (upper line) is clearly stained (dark band)
- IgG and IgA conjugate control bands (second and third band) respectively must show a clear staining
- Cut-off control (fourth band) must show a weak, but visible staining

If these conditions are not fulfilled, please check the following before repeating the test:

- expiration date of the reagents
- calibration of the used instruments
- exact test procedure
- visual examination of kit components for signs of contamination, deterioration or leakage; substrate solution must not be used if turned blue

If the quality conditions are not fulfilled after repeating, please contact your local distributor of R-Biopharm.

11.2. Test evaluation and interpretation

For evaluation, stick the corresponding test strip onto the appropriate fields on the evaluation form using a glue stick. Align the test strip with the reaction control bands along the marked lines. Sticking the entire test strip down flat using glue or tape can lead to changes in color. Therefore it is recommended to use a transparent adhesive to attach the test strip to the left of the marked line.

Identify the bands of the developed test strips using the printed control strip on the evaluation form and enter it into the evaluation form. Using table 4, carry out the evaluation of the intensity of the resulting bands separately for each of the corresponding immunoglobulin classes IgA and IgG.

Table 4: Evaluation of band intensity in relation to the cut-off band

Stain intensity of the bands	Evaluation
No reaction	-
Very low intensity (lower than cutoff-band)	+/-
Low intensity (equivalent to cut-off band)	+
Strong intensity (stronger than cut-off band)	++
Very strong intensity	+++

The test result is determined according to the sum of the scores in table 5 of the individual bands reactive over the cut-off. (i.e. with a minimum grade of +).

The positive, borderline or negative assessment of the sample can then be determined directly with table 6.

Table 5: Point evaluation of the antigens at RIDA® LINE Helicobacter IgA, IgG

Antigen	Function	Point evaluation
CagA	Cytotoxin associated protein	2
VacA	vacuole cytotoxin	2
GroEL	Chaperon	2
Urease A	Small subunit of urease	1
HcpC	supposable β -Lactamase	1
gGT	Gamma glutamyl transpeptidase	1

Table 6: Test interpretation

Sum of points	Assessment
0	negative
1	questionable
≥ 2	positive

A serum is **IgA- or IgG- positive** if the total point number is ≥ 2 . It is **IgA- or IgG-negative** if the total point number is 0. Sera with a point number of 1 are considered to be **questionable**. The type differentiation is carried out using CagA and VacA antigens and is only possible in the IgG. If the CagA and/or VacA bands react above the cut-off (i.e. with at least +) there is an infection of Helicobacter pylori type I.

12. Remarks about the test procedure and interpretation

A negative result for Helicobacter pylori cannot exclude an infection with Helicobacter pylori. If a clinical suspicion of an infection with Helicobacter pylori with a negative serological diagnosis subsists, after 4 weeks another patient sample should be taken and tested.

A questionable result could indicate the beginning of an immune response. Antibodies against Helicobacter pylori are present. A progressing sera conversion which points toward an acute illness should be clarified with a renewed test after 2 - 4 weeks. A questionable result can also be observed after therapy and an overcome infection.

A positive result with the RIDA® LINE Helicobacter IgA and/or IgG does not mean that an active illness is at hand for each case. Antibodies can also be detected with a subclinically developing cultivation of Helicobacter pylori. Due to possible pathological changes, an eradication of the organism is advisable.

Serological test results should always be seen in connection with the clinical picture. The direct organism detection is advised if unclear serological diagnoses are present.

Of very high diagnostic value for the interpretation, is the reactivity of antibodies against the vacuolating cytotoxin (VacA) and the cytotoxin associated protein (CagA), whereby the last mentioned is strongly immunogen and in the presence of the corresponding antibodies usually visible as a clearly pronounced band. Due to the low immunogenicity of VacA antibodies against this cytotoxin are rarely detectable or are only visible as weak bands. Normally the VacA band is only detectable in combination with the CagA band. The detection of antibodies against these proteins (CagA and VacA) indicates an infection with highly virulent Helicobacter pylori strains (type I). If the corresponding clinical diagnosis exists, a gastroscopy should be considered, if it has not been already done so, for the clarification of ulcera or carcinoma.

Little research is available which explains the significance of IgA antibodies with an infection of *Helicobacter pylori*. It is advised to carry out this examination if unclear or negative IgG diagnoses are present. A detection of IgA antibodies indicates an acute infection.

Some patient samples can produce a dark, uniform or patterned staining across the entire nitrocellulose strip (e.g. on sera from patients with milk protein allergies). Various factors in each patient are responsible for this. The evaluation of these strips is usually only partly feasible. Thus, “inverse” bands (white bands on dark background), for example, should be evaluated as negative. The respective serum should always be examined using other serological methods.

13. Clinical results

13.1. Diagnostic sensitivity of RIDA® LINE Helicobacter

Correlation between **RIDA® LINE Helicobacter** and histology/culture positive samples

RIDA® LINE Helicobacter	Histology/Culture positive samples (n=120)	
	IgG	IgA
Negative	4	37
Equivocal	3*	8*
Positive	113	75
Sensitivity	96,7%	69,2%

* equivocal results are considered positive

13.2. Diagnostic specificity of RIDA® LINE Helicobacter

Correlation between **RIDA® LINE Helicobacter** and histology/culture negative samples

RIDA® LINE Helicobacter	Histology/Culture negative samples (n=13)	
	IgG	IgA
Negative	13	13
Equivocal	0	0
Positive	0	0
Specificity	100%	100%

13.3. Prevalence blood donors

RIDA® LINE Helicobacter	Blood donors (n=100)	
	IgG	IgA
Negative	48	79
Equivocal	0	0
Positive	52	21
Prevalence	52%	21%

13.4. Analytical specificity

The analytical specificity is defined as the capacity of the test to determine analytical accuracy in the presence of potential interference factors in the sample matrix (e.g. anticoagulants, haemolysis, effects of sample treatment) or cross-reactions with potentially interfering antibodies.

(a) Interferences

Controlled studies of potentially interfering factors have shown that test performance was not influenced by anticoagulants (sodium citrate, EDTA, heparin), haemolysis (up to 1000 mg/dl haemoglobin), lipemia, bilirubinemia (up to 20 mg/dl bilirubin) or the deep-freezing and thawing cycles of the sample.

(b) Cross-reactions

The potential interference of antibodies against other organisms that have similar clinical symptoms to those caused by an infection with helicobacter pylori (e.g. Campylobacter jejuni) was examined in control studies. Additionally, conditions that are attributed to atypical activity of the immune system (antinuclear autoantibodies, rheumatoid factor) were tested. No cross-reactivities were detected.

Literature

1. Marshall, B. F., Royce, H., Annear, D. I., Goodwin, C. S., Taylor, N. S., Edmonds, P., Sly, L. I., Brenner, D. J.: Original Isolation of *Campylobacter pyloridis* from human gastric mucosa. *Microbios Letter* 25, 83 (1982)
2. Warren, J. R.: Unidentified curved bacilli on gastric epithelium in active chronic gastritis. *Lancet* I, 1273 - 1275 (1983)
3. Parsonnet, J., Friedman, G. D., Vanderstetten, D. P., Chang, Y., Vogelman, J. H., Orentreich, N., Sibley, R. K.: *Helicobacter pylori* infection and the risk of gastric carcinoma. *N. Engl. J. Med.* 325, 1127 - 1131 (1991)
4. Stolte, M., Eidt, S.: Healing gastric MALT lymphomas by eradicating *H. pylori*. *Lancet* 342, 568 (1993)
5. Stolte, M.: *Helicobacter pylori* and gastric MALT lymphoma. *Lancet* 339, 745 - 746 (1992)
6. Hentschel, E., Brandstatter, G., Dragosics, B., Hirschl, A. M., Nemec, H., Schutze, K., Taufer, M., Wurzer, H.: Effect of Ranitidine and Amoxicillin plus Metronidazole on the eradication of *Helicobacter pylori* and the recurrence of duodenal ulcer [see comments]. *N. Engl. J. Med.* 328, 308 - 312 (1993)
7. Stadelmann, O.: *Helicobacter pylori*: Indikationen und Praxis der Therapie. *Deutsches Ärzteblatt* 92, 39 - 41 (1995)
8. Telford, J. L., Ghiara, P., Dellorco, M., Comanducci, M., Burroni, D., Bugnoli, M., Tecce, M. F., Censini, S., Covacci, A., Xiang, Z. Y., Papini, E., Montecucco, C., Parente, L., Rappuoli, R.: Gene structure of the *Helicobacter pylori* cytotoxin and evidence of its key role in gastric disease. *J. Exp. Med.* 179, 1653 - 1658 (1994)
9. Lee, A.: *H. pylori* initiated ulcerogenesis: look to the host. *Lancet* 341, 281 (1993)
10. Schmitt, W., Haas, R.: Genetic analysis of the *Helicobacter pylori* vacuolating cytotoxin - structural similarities with the IgA protease type of exported protein. *Mol. Microbiol.* 12, 307 - 319 (1994)
11. Covacci, A., Censini, S., Bugnoli, M., Petracca, R., Burroni, D., Macchia, G., Massone, A., Papini, E., Xiang, Z., Figura, N., Rappuoli, R.: Molecular characterization of the 128-kDa immunodominant antigen of *Helicobacter pylori* associated with cytotoxicity and duodenal ulcer. *Proc. Natl. Acad. Sci. USA* 90, 5791 - 5795 (1993)
12. Xiang, Z., Censini, S., Bayeli, P. F., Telford, J. L., Figura, N., Rappuoli, R.: Analysis of expression of *cagA* and *vacA* virulence factors in 43 strains of *Helicobacter pylori* reveals that clinical isolates can be divided into two major types and that *cagA* is not necessary for expression of the vacuolating cytotoxin. *Infect. Immun.* 63, 94 (1995)

13. Krakowka, S., Morgan, D. R., Kraft, W. G., Leunk, R. D.: Establishment of gastric *Campylobacter pylori* infection in the neonatal gnotobiotic piglet. *Infect. Immun.* 55, 2789 - 2796 (1987)
14. Karita, M., Li, Q., Cantero, D., Okita, K.: Establishment of a small animal model for human *Helicobacter pylori* infection using germ-free mice. *Am. J. Gastroenterol.* 89, 208 - 213 (1994)
15. Marchetti, M., Arico, B., Burrioni, D., Figura, N., Rappuoli, R., Ghiara, P.: Development of a mouse model of *Helicobacter pylori* infection that mimics human disease. *Science* 267, 1655 - 1658 (1995)
16. Jaskowski, T. D., Martins, T. B., Hill, H. R., Litwin, C. M.: Immunoglobulin A Antibodies to *Helicobacter pylori*. *J. Clin. Microbiol.* 35, 2999 - 3000 (1997)
17. Yunoki N, Yokota K, Mizuno M, Kawahara Y, Adachi M, Okada H, Hayashi S, Hirai Y, Oguma K, Tsuji T: Antibody to Heat Shock Protein Can Be Used for Early Serological Monitoring of *Helicobacter pylori* Eradication Treatment. *Clin Diagn Lab Immunol*, July 2006, Vol 19, No.3, p449-490.
18. Annibale B, Lahner E, Santucci A, Vaira D, Pasquali A, Severi C, Mini R, Figura N, Fave GD: CagA and VacA are Immunoblot Markers of Past *Helicobacter pylori* in Atrophic Body Gastritis. *J. Compilation* © 2007 Blackwell Publishing Ltd, *Helicobacter* 12:23-30, *Helicobacter* ISSN 1523-5378
19. Gao L, Weck MN, Michel A, Pawlita M, Brenner H: Association between chronic atrophic gastritis and serum antibodies to 15 *Helicobacter pylori* proteins measured by multiplex serology. *Cancer Res* 2009 Apr 1; 69(7): 2973-80