RIDA[®] CCD-Inhibitor

Art. Nr.: ZA0601

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 2-8 °C

Gebrauchsanweisung
Instructions for use
Instrucciones de uso
Mode d'emploi
Инструкции по использованию
Indicação de uso

R-Biopharm AG, An der neuen Bergstraße 17, D-64297 Darmstadt,
Almanya, Tel.: +49 (0) 61 51 81 02-0 / Faks: +49 (0) 61 51 81 02-20

RIDA® CCD-Inhibitor (ZA0601)

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RIDA® CCD-Inhibitor (ZA0601)

Intended use:

The RIDA® CCD-Inhibitor is an accessory for use in *in vitro* diagnostics for the inhibition of anti-CCD IgE in human serum and plasma.

Package content:

CCD Inhibitor	5x RIDA® CCD-Inhibitor (lyophilised); Instructions for use
---------------	---

Additional necessary reagents and necessary equipment

Test tubes

Vortex mixer

Micropipettes for 10 µl, 55 µl and 400 µl volumes

General information:

A type I allergy is caused by the formation of specific IgE antibodies against allergens. Most allergens are proteins with a molecular weight of 50-500 kDa. However, in combination with a protein, small molecules such as drugs (haptens) can develop into allergens.

In eukaryotic organisms, most proteins are subjected to posttranslational glycosylation and therefore have carbohydrate side chains.

The immune system also develops IgE antibodies against carbohydrate side chains (cross-reactive carbohydrate determinants = CCDs), which are usually of vegetable origin.

The IgE antibodies against these CCDs also lead to cross-reactions with unrelated proteins, which, however, most likely have no clinical relevance and therefore do not trigger allergic symptoms. As these produce positive results in laboratory tests, these positive results must be considered false positive.

In order to correctly distinguish between true positive and false positive results, the anti-CCD IgE antibodies must be inhibited to prevent binding to the CCDs in the laboratory test.

Occurrence of many positive reactions in the single allergen test system is indicative of cross-reactions caused by anti-CCD IgE antibodies. In this case it is recommended to treat the serum with the RIDA® CCD-Inhibitor and to repeat the test.

Test principle:

The RIDA® CCD-Inhibitor bonds to the variable region of IgE antibodies and prevents the antibodies from binding to the carbohydrate side chains of the allergens in the test system.

Precautions:

All reagents and materials coming into contact with potentially infectious specimens must be treated with suitable disinfectants or autoclaved at 121 °C for at least one hour.

The reagent must not be used if the bottle is damaged or has a leak.

Storage:

The RIDA® CCD-Inhibitor must be stored unopened at 2-8 °C and can be used up to the expiration date printed on the label. After the lyophilisate has been dissolved, the ready-to-use solution can be kept for two weeks at 2-8 °C. Microbial contamination must be prevented. After the expiration date, the quality guarantee is no longer valid.

Performing the test:

The RIDA® CCD-Inhibitor must only be used by trained laboratory personnel. The guidelines for work in medical laboratories must be observed. The instructions for performing the test must be strictly complied with.

Dissolve the lyophilisate with 55 µl bidest. H₂O. In order to dissolve the lyophilisate first remove the test tube from the refrigerator and allow it to warm up to room temperature (RT) for approx. 30 minutes. Then add 55 µl H₂O and thoroughly vortex for 30 seconds. Then centrifuge it briefly to ensure that no liquid remains in the cap.

Pipet 10 µl of the dissolved RIDA® CCD-Inhibitor into 400 µl of serum or plasma, shake, and incubate for one hour at RT while shaking.

The treated serum or plasma must be tested immediately after incubation.

Explanation of symbols

	<i>In-vitro</i> diagnostic agent
	Observe the instructions for use
	Batch number
	Use before
	Storage temperature
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

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