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
## **SeraSpot<sup>®</sup> Anti-Treponema-4 IgG** **SeraSpot<sup>®</sup> Anti-Treponema-4 IgM**

Spotimmunoassay for detection of IgG or IgM antibodies to *Treponema pallidum*  
in human serum or plasma

**REF** SP-010-4 G-S6 ▽ 48    **REF** SP-010-4 G-S12 ▽ 96    **REF** SP-010-4 G-S24 ▽ 2x 96

**IVD** *In-vitro*- diagnostic device    **CE**

**REF** SP-010-4 M-S6 ▽ 48    **REF** SP-010-4 M-S12 ▽ 96    **REF** SP-010-4 M-S24 ▽ 2x 96

**IVD** *In-vitro*- diagnostic device    **CE**

**Seramun Diagnostica GmbH** · Sprenhagener Straße 1 · 15754 Heidesee · Germany · [www.seramun.com](http://www.seramun.com)  
phone +49 (0) 33767 79110 · fax +49 (0) 33767 79199 · [info@seramun.com](mailto:info@seramun.com)

## Introduction

*Treponema pallidum ssp. pallidum* is the causative agent of Syphilis, a sexually transmitted globally spread disease in humans. The incidence varies depending on the region and a consistent increase is observed since 2000 in western industrial nations with an actual level of 2-4 cases per 100,000 population per year (1).

The clinical course of Syphilis is characterized by different stages with varying symptoms. Infection occurs by sexual contacts. Skin lesions but also intact mucosa in the genital, oral or anal region are the most likely sites of entry for the pathogen and from there it spreads via blood and lymph vessels. During this primary stage a painless ulcer (ulcus durum) develops at the site of entry and a swelling of the regional lymph nodes (lymphadenopathy). Without treatment the infection will enter into the secondary stage after 9-10 weeks. Clinical signs are patchy papular exanthema or other ulcerative skin irritations, accompanied by several less characteristic symptoms like lymphadenopathy, fatigue, weight loss, subfebrile temperatures etc. The papular skin rashes contain a large number of pathogens. In the course of the secondary stage, which lasts between weeks to months, specific antibodies are produced. Nevertheless reinfections may occur. After decline of the symptoms and a latency period of up to several years the tertiary stage with a wide range of symptoms depending on the respective organ manifestations develops. Involvement of the central nervous system results in the clinical picture of Neurosyphilis or other neurological manifestations (1, 2).

Connatal *T. pallidum* infection by transplacental transmission during pregnancy causes spontaneous abortion, premature birth or stillbirth in 40% of the cases. About 50-70% of infected neonates are clinically unspectacular. Otherwise symptoms like bloody rhinitis, skin and mucosa irritations, fever, oedema, hepatosplenomegaly, lymphadenopathy, enteritis, hydrocephalus and meningitis may occur. In rare cases clinical symptoms do not develop before early childhood (3, 4).

Antibiotic treatment with Penicillin G is the method of choice for Syphilis therapy. Alternative antibiotics are tetracyclines and cephalosporines (3).

*Treponema pallidum* infections are diagnosed either by direct (pathogen detection or PCR) or indirect (serology) methods.

Direct pathogen detection from aspirates or smears of primary lesions, lymph nodes or exanthema by dark-field microscopy or direct immune fluorescence is hampered by the necessity to test the samples immediately after collection. Nucleic acid amplification techniques like PCR from aspirates, smears, blood or liquor are available but the diagnostic relevance is currently questionable (2, 5).

Laboratory diagnosis of Syphilis is preferably performed by detection of *T. pallidum* specific antibodies in serum, plasma or liquor. A two-step algorithm comprising a first screening and a subsequent confirmatory test is recommended (1, 2, 6). *Treponema pallidum* hemagglutination test (TPHA) or *Treponema pallidum* particle agglutination test and a polyvalent, IgG, IgA and IgM antibody detecting enzyme immunoassay (ELISA) are suitable screening methods. A positive or questionable screening result must be checked by a confirmatory test.

The Fluorescence Treponema Antibody Absorption test (FTA-ABS test) is an established confirmatory test. Treponema fixed to the surface of glass slides serve as target antigens for antibody detection from patient samples. In order to reach satisfying test specificity potentially cross reactive antibodies have to be eliminated by absorption with *Treponema phagedenis*. Immunoblots or line immunoassays are also suitable confirmatory tests (2, 4, 7). The major advantage of line immunoassays is the use of defined, highly purified recombinant antigens (Tp47, Tp17, Tp15 and TmpA), that are confirmed as *T. pallidum* specific (8). Therefore cross reactions appear unlikely and an absorption step is not needed (7). For differentiation between acute and recent infections confirmatory testing should always include the separate detection of IgG and IgM antibodies. Either immunoblot or line immunoassay or 19S-IgM-FTA-ABS test can be used to detect specific IgM antibodies for evaluation of the activity of the infection. Since IgM antibodies may persist for months up to years, a positive IgM result must always consider clinical findings for therapeutic decisions (2).

Quantitative detection methods like IgM-ELISA, 19S-IgM-Test or Lipoid antibody detection (e.g. VDRL-test, Cardiolipin-KBR) are recommended for serological follow-up after therapy (2, 4, 6).

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## Intended use

The *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM is an *in vitro* diagnostic device (spotimmunoassay, SIA) for the detection of *Treponema pallidum*-specific IgG or IgM antibodies in human serum or plasma samples. The test can be used as a confirmatory test after a positive or a questionable result of a syphilis screening test (TPHA-TPPA- or TPLA-test).

## Principle of the test

The *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM is a solid phase immunoassay based on the use of recombinant *Treponema* antigens as capture antigens printed in array arrangement (spot array) on the bottom of the wells of 96 well-microtiter plates. The recombinant antigens serve as capture molecules for antibodies against *Treponema* antigens. Bound antibodies are detected by Horseradish Peroxidase-(HRP)-labeled antibodies against human antibodies of IgG- or IgM-type by substrate reaction with hydrogen peroxide and 3,3',5,5'-Tetramethylbenzidine (TMB). At the site of formed immune complexes blue spots are developed by precipitated product from colorless substrate solution. Color intensity is correlated to the antibody concentration. Pale blue to dark blue spots are visible by eye.

**The test procedure for antibody detection is performed in 3 steps:**

### Step 1

Incubation of 1:101 diluted samples for 30 minutes at room temperature in the intended wells. After incubation unbound components are removed by aspirating the samples followed by 3 wash cycles with wash solution.

### Step 2

Incubation of the wells with anti-human IgG- or IgM-HRP-conjugate for 30 minutes at room temperature. After incubation unbound components are removed by aspirating the conjugate followed by 3 wash cycles with wash solution.

### Step 3

Incubation of the wells with the substrate solution *SeramunBlau*<sup>®</sup> spot dark for 30 minutes at room temperature. The reaction is stopped by aspirating the substrate solution, followed by tapping the plate/strip dry onto lint-free absorbent paper.

After finishing step 3, developed assays can be scanned either by using the *Seramun SpotSight*<sup>®</sup> plate or *Seramun SpotSight*<sup>®</sup> strip scanner; the obtained images are interpreted by using the Software *Seramun SpotSight*<sup>®</sup> scan. Taking account of the test-specific evaluation criteria, the test result is determined.

## Used *Treponema* antigens

Nomenclature	Characterization	Specificity
Tpp15	<i>Treponema pallidum</i> protein 15	15 kD
Tpp17	<i>Treponema pallidum</i> protein 17	17 kD
Tpp47	<i>Treponema pallidum</i> protein 47	47 kD
TmpA	<i>Treponema</i> membrane protein A	37 kD

**All listed proteins are highly specific for all stages of infection.**

## Kit components

			For 48 determinations	For 96 determinations	For 2 x 96 determinations
1	<b>WELLS</b>	Wells with arrays	<b>6 single breakable 8-well strips in a frame</b> Vacuum sealed with desiccant	<b>12 single breakable 8-well strips in a frame</b> Vacuum sealed with desiccant	<b>2 x 12 single breakable 8-well strips in a frame</b> vacuum sealed with desiccant
<b>IgG detection: color coding black; IgM detection color coding white</b>					
2	<b>WASHBUF CONC 10X</b>	Wash buffer	<b>100 ml concentrate</b> transparent bottle white cap	<b>100 ml concentrate</b> transparent bottle white cap	<b>2x 100 ml concentrate</b> transparent bottle white cap
3	<b>DIL</b>	Sample diluent	<b>55 ml ready-to-use solution</b> , colored red transparent bottle black cap	<b>2x 55 ml ready-to-use solution</b> , colored red transparent bottle black cap	<b>4x 55 ml ready-to-use solution</b> , colored red transparent bottle black cap
4	<b>CONJ HRP IgG</b>	Anti-human IgG-HRP-conjugate	<b>8 ml ready-to-use solution</b> , colored red transparent bottle red cap	<b>8 ml ready-to-use solution</b> , colored red transparent bottle red cap	<b>2x 8 ml ready-to-use solution</b> , colored red transparent bottle red cap
	<b>CONJ HRP IgM</b>	Anti-human IgM-HRP-conjugate	<b>8 ml ready-to-use solution</b> , colored green transparent bottle green cap	<b>8 ml ready-to-use solution</b> , colored green transparent bottle green cap	<b>2x 8 ml ready-to-use solution</b> , colored green transparent bottle green cap
5	<b>SUBSTR TMB</b>	Substrate 3,3',5,5'-Tetramethylbenzidine and hydrogen peroxide	<b>8 ml ready-to-use solution</b> , black bottle blue cap	<b>8 ml ready-to-use solution</b> , black bottle blue cap	<b>2x 8 ml ready-to-use solution</b> , black bottle blue cap
6	<b>COVER</b>	Adhesive film	<b>2x</b>	<b>2x</b>	<b>4x</b>
7	<b>SWAB</b>	Swab with 70% (v/v) Isopropyl alcohol	<b>3x2</b>	<b>6x2</b>	<b>12x2</b>

## Preparation and storage of samples

Sample collection should be done in a sterile manner. Serum is separated after clotting by centrifugation as soon as possible. Samples can be stored at 2...8 °C for a maximum of 48 hours. For longer storage times samples have to be stored at - 20 °C. Repeated freezing and thawing of samples should be avoided!

Serum or plasma (EDTA, citrate or heparin plasma) can be used. The samples have to be diluted 1 : 101 (10 µl sample and 1000 µl sample diluent) with the ready to use sample diluent.

## Materials required but not provided

- Adjustable single-channel micro-pipettes incl. pipette tips
- Adjustable multi-channel micro-pipettes incl. pipette tips
- Reagent container for multi-channel micro-pipettes
- Distilled or deionized water
- Lint-free absorbent paper
- Test tubes (1 ml) for sample dilution
- Beaker and measuring cylinder
- Washer for 96-well microtiter plates
- Collecting devices for infectious material
- Seramun *SpotSight*<sup>®</sup> plate or Seramun *SpotSight*<sup>®</sup> strip scanner with evaluation software Seramun *SpotSight*<sup>®</sup> scan

## Preparation and storage of reagents

### Kit size and expiry

One kit is designed for 1x 48, 1x 96 or 2x 96 determinations. The expiry date of each component is reported on its respective label; that of the complete kit on the outer box label. Upon receipt, all test components have to be kept at 2...8 °C, preferably in the original kit box. After opening all kit components are stable for at least 3 months, provided proper storage. The ready-to-use wash solution can be used for at least 1 month when stored at 2...8 °C.

## Reagent preparation

Allow all components to reach room temperature prior to use in the assay. The microtitration plate is vacuum-sealed in a foil with desiccant. The plate consists of a frame and strips with breakable wells. Allow the sealed plate to reach room temperature before opening. Unused wells should be stored refrigerated and protected from moisture in the original cover carefully resealed. Prepare a sufficient amount of wash solution by diluting one part of the wash buffer concentrate **WASHBUF CONC 10X** with 9 parts of distilled or deionized water. For Example: 10 ml **WASHBUF 10X** + 90 ml distilled water and vortex thoroughly. The substrate solution **SUBSTR TMB** must be protected from direct light! The test should not be performed with a substrate solution **SUBSTR TMB** that is colored dark or contains colored precipitates!

## Workplace requirements

**Processing the spotimmunoassays requires a clean and lint-free workplace. Adhesion of fibers to the plastic surface of the wells has to be avoided. Fibers may disturb optical image generation of arrays and cause erroneous results. Workplace and kit components should not be exposed to direct sunlight.**

## Assay procedure

Execution at room temperature (18 ... 25 °C) and at the specified incubation times.

### Important notes:

- **Avoid any contact or scratching on the bottom surface of the wells by pipette tips or washer needles. This results in irreparable damages of the spot array.**
- **All liquid reagents (sample, conjugate and substrate) have to be pipetted without causing air bubbles into the wells.**
- **Before starting imaging of the wells it is necessary to remove particles or fibers which may adhere underside the bottom of the wells using the swab **SWAB** provided in the kit.**


### Working steps

1. Warm all reagents and the required number of **WELLS** to room temperature (RT) before use. Mix gently without causing foam.
2. **Add 100 µl** of 1:101 diluted samples to each well.
3. Cover plate with adhesive film **COVER**. Incubate for **30 min** at room temperature.
4. Empty wells (aspirate and discard samples), then wash wells **3 x** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
5. Add **50 µl** of the ready-to-use conjugate solution **CONJ HRP IgG** or **CONJ HRP IgM** to each well.
6. Cover plate with adhesive film **COVER** and incubate for **30 min** at room temperature.
7. Empty wells (aspirate and discard conjugate), then wash wells **3 x** using **400 µl** wash solution per well (made of (2)). Tap plate dry onto lint-free absorbent paper.
8. Add **50 µl** of ready-to-use substrate solution **SUBSTR TMB** to each well.
9. Cover plate and protect from light. Incubate for **30 min** at room temperature.
10. Empty wells (aspirate and discard substrate). Tap plate dry onto lint-free absorbent paper.
11. Clean bottom of wells with Swab **SWAB**.
12. Perform image acquisition using the Seramun *SpotSight*<sup>®</sup> plate or the Seramun *SpotSight*<sup>®</sup> strip scanner and evaluate the results with the Seramun *SpotSight*<sup>®</sup> scan software. If image acquisition is performed with the scanner Seramun *SpotSight*<sup>®</sup> strip, transfer single 8-well microtiter plate strips to the scanners strip holder.

**It is recommended to use a 96-well microtiter plate washer for aspiration and wash steps!**

## Test evaluation

### Array layout

Parameter		Controls	
04	TmpA	01	Positive control (PC)
05	Tpp15	02	Cut-off control (CO)
06	Tpp17	03	Negative control (NC)
07	Tpp47	08	IgG conjugate control (GC)
		09	IgM conjugate control (MC)
		10	Serum control (SC)
		P	Well position marker

### Validity criteria for the test

The *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM tests include the following control spots:

- 1 Positive control (PC). Intensively stained spot, darker than cut-off control. Always stained.
- 2 Cut-off control (CO). Weakly stained spot. Used for result evaluation of parameter specific signals.
- 3 Negative control (NC). Pale spot with intensity lower than cut-off control.
- 4 Function control (IgG-, IgM-conjugate control, GC, MC). Intensively stained spots with different position for IgG and IgM antibody detection. Serving as antibody isotype control.
- 5 Serum control (SC). Intensively stained spot always stained in presence of sample. Absence of spot indicates absence of sample.

The test run is valid if the positive control, the cut-off control, the corresponding function control and the serum control are clearly visible. The color intensity of the negative control should not exceed the intensity of the cut-off control.

### Result interpretation

Evaluation of the test has to be performed by using the Seramun *SpotSight*<sup>®</sup> scan software.

Judgment of spot staining is performed according to the following classification:

Result Interpretation	IgG	IgM
<b>negative</b>	Max. 1 antigen spot > cut-off control	No antigen spot > cut-off control
<b>positive</b>	2 or more antigen spots > cut-off control	1 or more antigen spots > Cut-off control

### Limitations of the procedure

The final result interpretation should always consider clinical findings. Individual cases may require repeated investigations of samples taken at intervals of several weeks, e.g. for samples with borderline result. Impurities and cross contamination of reagents or samples by fungi or bacteria can lead to false positive as well as false negative results.

**Interference**

Hemolytic, lipemic and icteric samples (concentrations up to 500 mg/dl hemoglobin, 1000 mg/dl lipids, 20 mg/dl bilirubin C / bilirubin F) and Rheumatoid factors (500 U/ml) do not interfere with the test.

**Assay procedure**

Incorrect dilution of samples, incorrect washing of the wells or incorrect timing can produce erroneous results.

Air bubbles generated by too forceful pipetting of samples and/or reagents may cause uneven signal formation of arrays. These arrays cannot be used for evaluation.

Damaged arrays e.g. by scratching the well bottom with pipet tips or washer needles are not suitable for evaluation.

Any kind of fibers which may stick underside the bottom of the wells or to the array area cause erroneous results.

## Performance characteristics

### Diagnostic sensitivity

A collection of 231 and 203 samples, characterized by a reference test (line immunoassay) were tested in the *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM (“Initial sensitivity”). Samples with discrepant results were re-analyzed in a second reference assay. The results obtained in the reference test 2 were used for determination of the final status of the samples with discrepant results (“Amended sensitivity”).

IgG	n = 231	Reference test 1		
		positive	borderline	negative
<i>SeraSpot</i> <sup>®</sup> Anti-Treponema-4 IgG	positive	210	2	5
	borderline	0	0	0
	negative	3	0	11

Borderline samples are valued positive.

Initial sensitivity: **98.6 %**

IgM	n = 203	Reference test 1		
		positive	borderline	negative
<i>SeraSpot</i> <sup>®</sup> Anti-Treponema-4 IgM	positive	73	19	7
	borderline	0	0	0
	negative	1	2	101

Borderline samples are valued positive.

Initial sensitivity: **96.8 %**

Amended sensitivity: **97.9 %**

### Diagnostic specificity

A collection of 382 blood donor samples were tested in the *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM in comparison to a reference test (line immunoassay). Samples with discrepant results were re-analyzed in a second reference test and assessed as described in “Diagnostic sensitivity”.

IgG	n = 382	Reference test 1		
		positive	borderline	negative
<i>SeraSpot</i> <sup>®</sup> Anti-Treponema-4 IgG	positive	0	0	0
	borderline	0	0	0
	negative	0	0	382

Specificity: **100.0 %**

IgM	n = 382	Reference test		
		positive	borderline	negative
<i>SeraSpot</i> <sup>®</sup> Anti-Treponema-4 IgM	positive	0	0	8*
	borderline	0	0	0
	negative	0	0	374

Borderline samples are valued positive.

\* 2 samples were confirmed borderline and 1 sample was confirmed positive by a second reference test.

Initial specificity: **97.9 %**

Amended specificity: **98.7 %**

## Cross-reactivity

Serum samples from pregnant women (n=51) were analyzed by *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM. All samples were tested negative for IgG antibodies (specificity 100 %) and 2 samples positive for IgM antibodies (initial specificity 96.1 %). 1 sample was confirmed positive by a second reference test (amended specificity 98.0 %).

## Precision

Samples of known antibody titer were analyzed by *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgG / *SeraSpot*<sup>®</sup> Anti-Treponema-4 IgM. Color intensity (OD) of spots was recorded and the coefficients of variation (CV) were calculated.

IgG detection								
	Intra-assay CV		Inter-assay CV		Interbatch CV		Total-assay CV	
	OD [x̄]	CV[%]	OD [x̄]	CV[%]	OD [x̄]	CV[%]	OD [x̄]	CV [%]
Tpp15	54.6	6.7	111.7	9.7	63.9	8.8	126.2	15.2
Tpp17	120.2	3.1	110.3	76	73.5	9.6	67.3	14.7
Tpp47	107.0	6.4	103.3	7.7	92.0	8.7	118.7	10.4
TmpA	112.6	5.3	98.1	6.1	113.7	11.0	118.3	14.3
CO	37.2	9.4	37.8	6.1	40.5	9.8	44.4	12.5
Expected value		< 10		< 15		< 20		< 25
Procedure	1 operator, 24x determination, 1 batch		1 operator, 10x 2 determination on 3 days, 1 batch		1 operator, 8x 3 determination, 3 batches		2 operators 2x8x3 determination, 3 batches	

IgM detection								
	Intra-assay CV		Inter-assay CV		Interbatch CV		Total-assay CV	
	OD [x̄]	CV[%]	OD [x̄]	CV[%]	OD [x̄]	CV[%]	OD [x̄]	CV [%]
Tpp15	72.0	9.8	48.2	10.7	79.8	13.3	80.7	17.7
Tpp17	111.7	8.2	76.6	8.5	116.8	10.3	110.8	12.2
Tpp47	50.9	9.1	74.4	12.5	56.0	16.6	71.4	18.2
TmpA	69.5	10.9	59.7	12.1	66.1	11.5	74.0	17.5
CO	37.2	9.4	25.5	14.0	30.5	9.8	32.3	11.1
Expected value		< 10		< 15		< 20		< 25
Procedure	1 operator, 24x determination, 1 batch		1 operator, 10x 2 determination on 3 days, 1 batch		1 operator, 8x 3 determination, 3 batches		2 operators 2x8x3 determination, 3 batches	

## Automation

Samples of known antibody titer (n=96, IgG and n=64, IgM) were assayed manually or by use of the microplate processor Dynex DS2<sup>®</sup> (manual sample dilution). Color intensity (OD) of spots was recorded and used for calculation of the coefficient of determination  $r^2$  for selected antibody reactivities and the different methods of operation.

Antigen	Automated (DS2 <sup>®</sup> ) vs. manual procedure	
	$r^2$ IgG	$r^2$ IgM
Tpp15	0.979	0.968
Tpp17	0.979	0.990
Tpp47	0.975	0.980
TmpA	0.979	0.970

## Common advices and precautions

This kit is for *in-vitro* use only. Follow the working instruction carefully. The test should be performed by trained technical staff only. The test instructions have to be followed strictly. The expiry dates stated on the respective labels are to be observed. Do not use reagents from damaged packages or bottles.

**Do not use or mix reagents from different lots except for sample diluent, washing buffer and TMB substrate solution. These can be used across batches and parameters.**

Do not use reagents from other manufacturers to complete the kit. Some of the reagents (2, 3, 4) contain small amounts of Thimerosal (0.01 % v/v (2)), Kathon (0.2 % v/v (3)), or Bronidox (0.1 % v/v (4)) as preservative. They must not be swallowed or allowed to come into contact with skin or mucous membranes. In case of contact immediately remove by rinsing with water. All reagents should be stored at 2...8°C. Handle all components and all patient samples as if potentially hazardous.

In particular the following precautions should be observed:

- **Do not smoke, eat or drink while handling kit material!**
- **Never pipette by mouth!**
- **Always use protective gloves!**
- **Note safety precautions of the single test components!**

## Incubation scheme SeraSpot® Anti-Treponema-4 IgG / IgM

-  100 µl  
30 min  
 3 x wash  
diluted sample (1 : 101)  
incubation (room temperature)  
with washing solution, each 400 µl per well
-  50 µl  
30 min  
 3 x wash  
conjugate **CONJ HRP IgG** or **CONJ HRP IgM** (4)  
incubation (room temperature)  
with washing solution, each 400 µl per well
-  50 µl  
30 min  
 aspiration  
substrate **SUBSTR TMB** (5)  
incubation (room temperature, protected from light)
- imaging of wells and image analysis  
scanner Seramun SpotSight® plate / Seramun SpotSight® strip and software Seramun SpotSight® scan

SeramunBlau®, SeraSpot® und Seramun SpotSight® are registered trademarks of Seramun Diagnostica GmbH.