## **RIDA®CHIP FoodGuide**

**REF** A8501 RIDA®CHIP FoodGuide 25

**REF** A8502 RIDA®CHIP FoodGuide 50

**REF** A8503 RIDA®CHIP FoodGuide 100

Deutsch .....	3
English.....	19
Español.....	36
Français.....	52
Italiano .....	69

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







**REF** A8501, A8502, A8503

### 1. Intended use

For *in vitro* diagnostic use. The RIDA<sup>®</sup>CHIP FoodGuide test is a quantitative microarray enzyme immunoassay (M-EIA) for the determination of specific IgG antibodies against food allergens. The test is designed for manual processing. The test is suited for the use of human serum and human capillary blood as sample matrices. The test is a diagnostic aid in cases of suspected IgG-mediated food hypersensitivities.

### 2. Summary and explanation of the test

The human immune system provides important protection against harmful exogenous substances. It mediates the identification and inactivation of, for example, bacteria, viruses, and parasites as part of a complex immune response. If the immune response erroneously targets non-pathogenic molecules, such as food, this is referred to as a hypersensitivity reaction or an allergy. The delayed type of such allergic reactions to food, a type III allergy, is mediated by IgG antibodies. Due to the high permeability of the intestinal wall, food proteins may cross the intestine barrier unnaturally, coming into contact with the immune system and ultimately triggering inflammatory reactions. The associated symptoms occur in a delayed manner and range from gastrointestinal symptoms and joint pain to migraines. Because the symptoms are delayed, it is often difficult to establish a causal relationship between the food and the onset of symptoms. Detecting the food-specific IgG antibodies in human serum, however, allows the triggering allergens to be easily identified, and the symptoms can be alleviated by avoiding the corresponding food.

### 3. Test principle

The test is a microarray-based enzyme immunoassay (M-EIA) to detect IgG antibodies against food allergens. Food extracts are elaborately extracted and purified and deposited on a microarray. Each well of the RIDA<sup>®</sup>CHIP FoodGuide contains several food allergens. Each well also has two standard series for quantification as well as positive and negative controls. The patient samples (serum or capillary blood eluate) are pipetted into the wells and incubated at room temperature. During incubation, specific IgG antibodies bind to the corresponding

adsorbed food antigens. Unbound material is removed by washing. An anti-human IgG antibody conjugated with horseradish peroxidase is then added. During a second incubation, this antibody enzyme conjugate binds to the human IgG antibodies of the patient sample. Any unbound conjugate is removed by washing. When the substrate is added, it is oxidized by the horseradish peroxidase and turns into a blue, insoluble product. The amount of blue precipitate which is formed is proportional to the quantity of antigen-specific antibodies in the serum and can be detected photographically. The RIDASOFT® FoodGuide software can be used to quantify these data and create comprehensive reports.

#### 4. Reagents provided

**Table 1:**

Slide	3 x	Slide with 1 x 8 wells Enables the detection of 101 antigens for a maximum of 12 patients (A8503) or the detection of 25 or 50 antigens for a maximum of 24 patients (A8501 and A8502, respectively).
Sample Buffer	2 x 15 ml	Ready-to-use sample dilution buffer to dilute the patient sample and elute capillary blood samples
Activation Buffer	5 ml	Ready-to-use activation buffer for CHIP preparation
Conjugate IgG	5 ml	Ready-to-use, diluted anti-human IgG conjugate Horseradish peroxidase-conjugated antibodies in stabilized protein solution
Substrate	5 ml	Ready-to-use substrate solution. Contains TMB.
Wash buffer salt Tween	2 x	Wash buffer salt for producing 2 x 1 liter of wash buffer. 10 mM PBS, 0.05% TWEEN20.

#### 5. Storage instructions

The test kit must be stored at 2 to 8°C and can be used until the expiration date printed on the label. The diluted wash buffer has a maximum shelf life of 4 weeks when stored at 2 to 8°C. Microbial contamination must be prevented. After the expiration date, the quality guarantee is no longer valid.

Contamination of the substrate solution **Substrate** with conjugate **Conjugate IgG** must be avoided as this will cause the substrate to deteriorate. Direct exposure of the substrate to light must also be avoided to prevent denaturing or discoloration. See 10.2, Indication of instability or expiration of reagents.

## 6. Reagents required but not provided

### 6.1 Reagents provided

- Distilled or deionized water

### 6.2 Equipment

- Slide carrier **Slide Carrier** for RIDA®CHIP FoodGuide slide **Slide**
- Reader for RIDA®CHIP FoodGuide slide **Slide**  
(R-Biopharm product number: ZCHIP READER SSCL)
- Microtiter plate washer **Washer** for 8 wells
- Microliter single-channel fine pipette (e.g., 1-25 µl)
- Microliter single-channel dispensing pipette or microliter multi-dispensing aid
- Suitable sample diluent vials

## 7. Warnings and precautions for the users

- For *in vitro* diagnostic use only.
- This test must be carried out only by trained laboratory personnel. The guidelines for working in medical laboratories must be followed. Always adhere strictly to the user instructions for carrying out this test.
- The product is intended for use by professionals working in hospital laboratories, reference laboratories, private laboratories, or public laboratories.
- Samples or reagents must not be pipetted by mouth, avoid contact with injured skin or mucous membranes. Wear disposable gloves when handling samples and wash the hands after the test.
- Do not smoke, eat, or drink in areas where samples or test reagents are being processed.
- Activation buffers, sample dilution buffers, conjugates, and wash buffer concentrates contain sodium azide as a preservative. Even though the concentration of sodium azide is below the classification limit for identifying hazardous materials, contact with skin and mucous membranes should be avoided. Contact with lead or copper pipes may cause the formation of explosive metal azides.
- Do not allow reagents to touch the skin, eyes, or clothing! In case of skin contact, wash area thoroughly with soap and water immediately. In case of eye contact, immediately flush under running water for 15 minutes with the eyelids open. Consult an eye specialist. If swallowed, drink plenty of water, avoid vomiting, and seek immediate medical attention.
- All components of the kit must be appropriately disposed of by the user after use (see safety datasheets for the components in the kit).
- All reagents and materials that come into contact with potentially infectious samples must be treated with suitable disinfectants and/or autoclaved at 121°C for at least one hour.

## 8. Collection and storage of samples

The RIDA<sup>®</sup>CHIP FoodGuide test was developed for testing human serum or capillary blood.

After venous blood collection, the serum should be separated from the blood clots (after complete clotting) as quickly as possible to prevent hemolysis. The samples must be stored cold (2 to 8°C) or frozen (-20°C) until they are tested. Avoid repetitive freezing and thawing of the serum as well as microbial contamination. The use of heat-inactivated, lipemic, hemolytic, icteric, or turbid sera may lead to false results.


When capillary blood samples are used in the RIDA<sup>®</sup>CHIP FoodGuide test, appropriate equipment must be used (A8025/A8025-IMU/A8025-PIM RIDASCREEN<sup>®</sup> Foodscreen Blood Collection Kit, A8025-BCC RIDASCREEN<sup>®</sup> Foodscreen Blood Collection Card). Only this equipment was validated for this test. The dried capillary blood samples on the Blood Collection Card (component in A8025 and A8025-BCC) are stable for a maximum of 3 weeks when stored in a dry location at room temperature (20°C to 25°C). Do not store capillary blood samples in the refrigerator.

**Table 2:** Sample stability

Undiluted serum		Diluted serum	Dried capillary blood sample	Eluate of capillary blood sample
2 to 8°C	-20°C	20 to 25°C	20 to 25°C	20 to 25°C
1 month	> 3 months	Max. 6 hours	3 weeks	Max. 6 hours

### Important note for processing capillary blood samples:

The circles on the blood card must be completely filled with blood. Also check the back of the card to ensure the circles are soaked through.

## 9. Test procedure

### 9.1 General information

**Important:** All reagents, patient sera/capillary blood samples, and the RIDA®CHIP FoodGuide slide [Slide] must be brought to room temperature before use. Mix the reagents well immediately before use.

The RIDA®CHIP FoodGuide slide [Slide] cannot be re-used. The components of the kit must not be used if the packaging is damaged or the vials are leaking.

**Any deviation from the specified incubation times and temperatures will cause the standard series values to shift in comparison with the certificate.**

**Significant differences in the values of the standard series can lead to invalid test results.**

**R-Biopharm assumes no warranty for any modification of incubation time and/or temperature that deviates from this description, for example, for automated processing of the test.**

Do not conduct the test in direct sunlight. It is recommended to cover the slide [Slide]. The manual test procedure takes a total of approximately 2 hours and 50 minutes.

### 9.2 Preparing the wash buffer

Dissolve one package of wash buffer salt in 1 liter of distilled or deionized water. Prepare only the amount of wash buffer needed for the particular assay. The diluted wash buffer has only a limited shelf life of max. 4 weeks at a storage temperature of 2 to 8°C.

### 9.3 Preparation of samples

#### 9.3.1 Sample dilution; serum samples

Dilute all serum samples with the sample dilution buffer [Sample Buffer] in a 1:101 ratio. To avoid systematic errors, we recommend making 1010 µl of diluted sample consisting of 10 µl serum and 1000 µl sample dilution buffer [Sample Buffer]. Make sure that the sample is mixed thoroughly. The sample must be pipetted several times and inverted and vortexed several times to ensure that it is thoroughly mixed. The diluted serum sample is stable for a maximum of 6 hours and must be processed within this time.

#### 9.3.2 Sample dilution; capillary blood samples

The dried capillary blood sample will be shaken with sample dilution buffer [Sample Buffer] for 1 (to 14) hour(s) at 20 to 25°C on a horizontal shaker in a suitable vial. For this step, mix 1 circle filled with blood [Card] with 2 ml sample dilution buffer [Sample Buffer]. Wear disposable gloves to pop out the dried circles filled with blood [Card] using your finger. A suitable sterile object, such as the end of

a pipette tip or a laboratory spatula, may also be used to pop out the circles **Card**. The eluate of the capillary blood sample already corresponds to the diluted patient sample and is pipetted into the wells **Wells** of the slide **Slide** according to the instructions.

#### 9.4 Preparation of the CHIP

The RIDA®CHIP FoodGuide slide **Slide** must be activated prior to use. For this purpose, place the desired number of slides **Slide** in the slide carrier **Slide Carrier**. 12 slides **Slide** can be used per slide carrier **Slide Carrier**. Next, pipette 100 µl of activation buffer **Activation Buffer** into each well and incubate for 5 minutes at 20 to 25°C.

#### 9.5 Washing

Use of a microplate washer is expressly recommended. R-Biopharm assumes no warranty for accurate processing if manual washing systems are used. The wells **Wells** are washed in overflow mode with 500 µl wash buffer **Wash buffer salt Tween** each. Final siphoning should be as thorough as possible; programming of a longer washing time is recommended.

**Important: Do not tap out the slide carrier **Slide Carrier** onto cellulose or similar materials as may normally be performed with other ELISA tests.**

#### 9.6 First incubation (sample incubation)

Pipette 100 µl each of the patient sample (diluted patient serum or capillary blood eluate) into the wells **Wells** of the slide **Slide** according to the pipetting scheme. Important: Depending on the test variant, either every well **Wells** (A8501, A8502) or 2 wells **Wells** each (A8503) must be used for one patient sample (see Figure 1). The slides **Slide** are then incubated for 45 minutes at 20 to 25°C. If multiple slide carriers **Slide Carrier** are used, they must be marked accordingly. The slide carrier **Slide Carrier** should be covered during incubation if possible.

**Fig. 1:** Pipetting scheme of individual test variants of RIDA®CHIP FoodGuide, S = patient sample

### 9.7 Washing

Wash as described in section 9.5.

### 9.8 Second incubation (conjugate incubation)

Pipette 100 µl (ready-to-use) conjugate **Conjugate IgG** into each well **Wells** of the slide **Slide**. Incubate the slides **Slide** filled with conjugate **Conjugate IgG** for 30 minutes at 20 to 25°C. The slide carrier **Slide Carrier** should be covered during incubation if possible.

### 9.9 Washing

Wash as described in section 9.5.

### 9.10 Third incubation (substrate incubation)

Pipette 100 µl (ready-to-use) substrate **Substrate** into each well **Wells** of the slide **Slide**. Incubate the slides **Slide** filled with substrate **Substrate** for 15 minutes at 20 to 25°C. The slide carrier **Slide Carrier** should be covered during incubation if possible.

### 9.11 Stopping the color reaction and performing the measurement

Stop the reaction by washing the wells **Wells** of the slide **Slide** twice each with 1000 µl deionized H<sub>2</sub>O in overflow mode. Omit the final siphoning in this step. **Leave the remaining deionized H<sub>2</sub>O in the wells **Wells** of the slide **Slide**. The remaining water will not affect the evaluation of the test. The evaluation must be performed within 2 hours, otherwise the results will be compromised.** The slides **Slide** can be retained for a second evaluation; for this purpose, the fluid must be siphoned off and the slide **Slide** must be completely dried in the drying cabinet at 37°C to 45°C (duration: about one hour). The slides **Slide** are placed in

the reader **Reader** and documented as specified in the instructions for the reader **Reader** and the software manual.

## 10. Quality control – Indication of instability or expiration of reagents

### 10.1 Quality control of the results

Each well **Wells** of the slide **Slide** contains two full standard series, one negative control (in duplicate), and two different positive controls (each in duplicate). The evaluation software (RIDASOFT® FoodGuide) evaluates the standard series and the individual controls. All specifications are stored in the software and checked automatically. The stored specifications are listed in Table 3.

The test has run correctly if the following specifications are met:

**Table 3:** Specifications for standards and controls

Standards	Min. intensity [%]	Max. intensity [%]
Standard 1		< Standard 2
Standard 2	-	< Standard 3
Standard 3	-	< Standard 4
Standard 4	-	< Standard 5
Standard 5	65.0	100

Controls	Min. concentration	Max. concentration
	µg/ml	µg/ml
Negative control	-	< 3.5
Positive control 1	12.9	23.0
Positive control 2	28.9	45.4

## 10.2 Indication of instability or expiration of reagents

Deviation from the stipulated values, turbidity of the reagents, or a blue precipitate in the reagent vial of the substrate may indicate that the reagent has expired.

If the stipulated values are not met, the following things must be checked before the test is repeated:

- Kit shelf life
- Functional performance of the equipment used (e.g., calibration)
- Correct test procedure according to the instructions for use
- Visual inspection of the kit components for contamination or leaks

**If the stipulated values are not obtained after repeating the test, please contact the manufacturer.**

## 11. Evaluation and interpretation

### 11.1 Basis for the calculation

To evaluate the test, a calculation based on the standard series with the aid of the software (RIDASOFT® FoodGuide) is required. The concentrations of the specific IgG antibodies in µg/ml are determined on the basis of the standard series from the measured intensities and then converted into IgG classes (see Table 4).

The standard series for RIDA®CHIP FoodGuide is calibrated based on an international reference preparation: 1st WHO IRP 67/86 for human IgG.

### 11.2 Concentrations, IgG classes, and calculations for RIDA®CHIP FoodGuide

**Table 4:** Relationship between the determined µg/ml, IgG classes, and antigen-specific IgG concentration of the patient from a serum or capillary blood sample

µg/ml	IgG class	Antigen-specific IgG content
< 10.0	0	Negative
≥ 10.0 < 20.0	1	Elevated
≥ 20.0	2	Highly elevated

These values are rough guidelines only. International standards for specific IgGs to foods do not exist.

## 12. Limitations of the method

The IgG concentrations determined with this test system provide information about the level of sensitization of the patient with regard to the investigated food antigens or antigen mixtures.

A correlation between the level of a determined IgG concentration and the occurrence or the severity of clinical symptoms cannot be deduced from this. The results obtained must always be interpreted in combination with the complete clinical symptoms and are only a diagnostic aid.

Negative results do not rule out IgG-mediated food intolerance and should not be used as the sole basis for diagnosis.

False positive test results may occur due to cross-reactivity of the tested allergen with epitopes of other allergens.

It cannot be ruled out that antigenic epitopes are missing due to the production processes, e.g., extraction of the food or coating of the microtiter plates. Potentially missing epitopes may lead to false negative results. IgG antibodies against food antigens that appear only during industrial preparation, food preparation, or during digestion may not be detectable because they are not present in the original food for which the patient is tested.

### 13. Performance characteristics

All performance characteristics of RIDA®CHIP FoodGuide were determined for all 101 antigens according to the current CLSI guidelines. The evaluations were conducted for each antigen, but here, the data are presented in a summarized form for greater clarity.

#### 13.1 Precision

The precision or reproducibility of the RIDA®CHIP FoodGuide test was determined using eight reference samples that cover the entire measurement range from negative to highly positive. Intra-assay reproducibility was tested 12 times for each reference sample. For inter-assay reproducibility, references from six sequential working days were assayed in triplicate, with two runs per day. Each measurement was performed by two technicians. Inter-lot reproducibility was calculated using the average values and the coefficients of variation (CV) for three reagent lots; the data were calculated using single-factor analysis of variance (ANOVA).

**Table 5:** Intra-assay precision analysis

Allergen	CV (%)	Allergen	CV (%)	Allergen	CV (%)
Almond	6.5	Gluten	5.0	Pineapple	6.8
Apple	7.0	Goat milk and cheese	15.2	Pistachio	8.2
Apricot	12.2	Grape	8.8	Pollock	7.7
Aspergillus niger	11.4	Green bean	9.5	Poppy	5.8
Banana	9.4	Guar flour (E412)	10.5	Pork	8.2
Barley	8.5	Hazelnut	8.5	Potato	9.2

Basil	9.2	Honey	9.2	Pumpkin seed	9.4
Beef	6.1	Horseradish	8.2	Quinoa	8.6
Beetroot	7.4	Kiwi	10.2	Raspberry	10.4
Black pepper	15.7	Kohlrabi (turnip cabbage)	6.8	Red cabbage	8.9
Black tea	18.5	Lamb	8.6	Rennet cheese	19.6
Broccoli	8.7	Lamb's lettuce	11.3	Rice	7.8
Buckwheat	6.7	Leek	7.8	Rosemary	12.8
Butterhead lettuce	6.9	Lemon	10.5	Rye	8.2
Candida albicans	13.1	Linseed	8.9	Salmon	6.6
Cane sugar	6.4	Lobster	9.6	Sesame	9.5
Carrot	6.5	Maize, sweet corn	7.2	Sheep milk and cheese	14.2
Cashew kernel	9.7	Meadow mushroom	10.1	Sour-milk products (cow)	19.6
Celeriac, knob celery	7.6	Millet	17.8	Soybean	6.8
Cherry	8.5	Mustard seed	27.2	Spelt	9.8
Chicken	10.8	Nectarine	9.2	Strawberry	8.1
Chicken egg-white	5.4	Nutmeg	8.5	Sunflower seed	8.4
Chicken egg-yolk	12.3	Oats	12.7	Sweet pepper	5.6
Chili cayenne	9.9	Ocean perch	9.3	Thyme	6.3
Cinnamon	15.6	Olive	11.1	Tomato	8.5
Cocoa bean	8.7	Onion	8.4	Tunafish	12.1
Cod, codling	9.9	Orange	7.3	Turkey	10.9
Coffee	13.1	Oregano	10.2	Vanilla	8.5
Cow's milk	3.7	Oyster mushrooms	21.6	Walnut	9.7
Crayfish	22.2	Paprika, spice	16.1	Watermelon	7.9
Cucumber	12.6	Parsley	15.9	Wheat	7.2
Curry	11.3	Pea	13.0	Yeast	11.1
Eggplant	10.1	Peanut	9.4	Zucchini (Courgette)	6.6
Garlic	9.3	Peppermint tea	11.3		

**Table 6:** Summary of intra-assay precision analysis

Number of sera	8
Replicates (per serum)	12
<b>Average CV</b>	<b>10.1%</b>

**Table 7:** Inter-assay precision analysis

Allergen	CV (%)	Allergen	CV (%)	Allergen	CV (%)
Almond	8.6	Goat milk and cheese	15.4	Pollock	17.9
Apple	11.2	Grape	17.9	Poppy	8.5
Apricot	18.7	Green bean	6.6	Pork	7.2
Aspergillus niger	11.8	Guar gum (E412)	9.2	Potato	11.4
Eggplant	19.2	Hazelnut	11.0	Pumpkin seeds	9.2

Banana	15.3	Honey	8.8	Quinoa	10.6
Barley	15.0	Horseradish	10.1	Raspberry	6.5
Basil	11.3	Kiwi	9.3	Red cabbage	8.3
Beef	6.8	Kohlrabi (turnip cabbage)	13.2	Rennet cheese (cow)	17.2
Beetroot	14.4	Lamb	14.2	Rice	8.6
Black pepper	17.9	Lamb's lettuce	6.7	Rosemary	21.6
Broccoli	8.0	Leek	10.4	Rye	8.3
Buckwheat	8.1	Lemon	16.6	Salmon	11.7
Butterhead lettuce	8.6	Linseed	7.0	Sesame	15.3
Candida albicans	11.5	Lobster	12.2	Sheep milk and cheese	15.6
Cane sugar	12.7	Maize, sweet corn	13.6	Sour-milk products (cow)	13.7
Carrot	7.6	Meadow mushroom	11.4	Soybean	8.3
Cashew nut	12.1	Millet	12.2	Spelt	12.9
Cayenne pepper	13.1	Mustard seed	22.6	Strawberry	15.6
Celeriac, knob celery	7.8	Nectarine	12.2	Sunflower seeds	10.0
Cherries	9.7	Nutmeg	10.4	Sweet pepper	8.3
Chicken	13.7	Oats	8.0	Tea, black	14.6
Chicken egg white	6.4	Ocean perch	13.3	Tea, peppermint	9.1
Chicken egg yolk	11.7	Olive	11.0	Thyme	12.3
Cinnamon	11.8	Onion	9.9	Tomato	9.6
Cocoa beans	15.8	Orange	10.5	Tunafish	16.4
Cod, codling	8.6	Oregano	18.0	Turkey	11.2
Coffee	14.0	Oyster mushrooms	16.5	Vanilla	11.0
Cow's milk	5.9	Paprika, spice	25.1	Walnut	10.4
Crayfish	19.8	Parsley	13.4	Watermelon	7.9
Cucumber	8.4	Pea	11.4	Wheat	10.7
Curry	12.1	Peanut	7.2	Yeast	9.0
Garlic	7.1	Pineapple	12.4	Zucchini (Courgette)	10.8
Gluten	9.9	Pistachio	8.7		

**Table 8:** Summary of inter-assay precision analysis

Number of sera	8
Replicates (per serum)	36 (6 days, 2 runs, each in triplicate)
<b>Average CV</b>	<b>11.8%</b>

**Table 9:** Summary of inter-lot precision analysis

Number of lots	3
Number of sera	8
Replicates (per serum)	3
<b>Average CV</b>	<b>18.0%</b>

## **13.2 Stability**

### **13.2.1 Transport stability**

The kit components were tested for stability under temperature fluctuations due to transportation. Two simulated transport scenarios and one genuine overseas transport were investigated. Temperature fluctuations between 3°C and 43°C were induced during the transport simulations. In the genuine transport investigation, the temperature fluctuations during the entire transportation period were documented. Temperatures between 1°C and 24°C were recorded. The stressed kits were compared with the reference kit from the same lot. The reference kits were stored at 2 to 8°C. A total of 4 sera were used for the analysis.

**No significant impact on results was detected after the kit components were transported, which means RIDA®CHIP FoodGuide is stable under the analyzed transport conditions.**

### **13.2.2 Kit stability after accidental freezing**

The kit components were tested for stability after freezing for 24 hours. The kit components were stored at -16°C for 2 days. The stressed kits were then compared with the reference kit from the same lot. The reference kits were stored at 2 to 8°C. A total of 4 sera were used for the analysis.

**Note: Freezing of kit components is not a permissible method of storage.**


**No significant impact on results was observed after 2-day freezing of the kit components.**

### **13.2.3 Kit stability after opening**

**Note:** The test is intended for single use only so it does not contain any components that may be reused. For this reason, “kit stability after opening” was not determined for RIDA®CHIP FoodGuide.

## **13.3 Analytical sensitivity**

Analytical sensitivity is a measurement of the accuracy of a test at low concentrations of the analyte.

**Fig. 2:** Box-whisker plot shows the distribution of all limit of blank (LOB) and limit of detection (LOD) values determined; green = LOB, red = LOD. The boxes represent 50% of the totality of data. The whiskers show overall distribution of the data. The horizontal line in the box is the median.

**For RIDA®CHIP FoodGuide, the median LOB is 2.4 µg/ml and the median LOD is 3.2 µg/ml.**

### 13.4 Interferents

The following potential interfering substances (interferents) were tested in the sample: hemoglobin, bilirubin, triglycerides.

**None of the listed substances were observed to impair the test results of the RIDA®CHIP FoodGuide.**

### 14. Version history

Version number	Section and designation
2017-10-12	Previous version
2020-04-30	General revision: new layout








2020-04-30	<ol style="list-style-type: none"> <li>1. Intended use</li> <li>4. Reagents provided</li> <li>7. Warnings and precautions for the users</li> <li>8. Collection and storage of samples</li> <li>9. Test procedure</li> <li>10. Quality control – Indication of instability or expiration of reagents</li> <li>11. Evaluation and interpretation</li> <li>12. Limitations of the method</li> <li>13. Performance characteristics</li> <li>15. Test-specific symbols</li> </ol>
2022-06-13	Revision
2022-06-13	<ol style="list-style-type: none"> <li>4. Reagents provided</li> <li>9. Test procedure</li> </ol>

## 15. Explanation of symbols

### General symbols

	In vitro diagnostic use
	Consult instructions for use
	Batch number
	Use before
	Store at
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

### Test-specific symbols

## 16. References

A. Aljada et al. (2004) Increase in intranuclear nuclear factor B and decrease in inhibitor B in mononuclear cells after a mixed meal: evidence for a proinflammatory effect. *Am J Clin Nutr* 79:682–90.

A. R. Gaby. (1998) The Role of Hidden Food Allergy/Intolerance in Chronic Disease. *Alternative Medicine Review*, 3(2):90-100.

El Aydinlar et. al (2013) IgG-based elimination diet in migraine plus irritable bowel syndrome. *Headache*. 53(3):514-25.

G. E. Mullin et al. (2010) Testing for Food Reactions: The Good, the Bad, and the Ugly. *Nutr Clin Pract* 25(2):192-198.

H. Uzunısmail et. al (2012) The effects of provocation by foods with raised IgG antibodies and additives on the course of Crohn's disease: a pilot study. *Turk J Gastroenterol*. 23:19-27.

K. Alpay et. al (2010) Diet restriction in migraine, based on IgG against foods: a clinical double-blind, randomised, cross-over trial. *Cephalalgia*. 30(7):829-37.

S. Bentz et al.(2010) Clinical relevance of IgG antibodies against food antigens in Crohn's disease: a double-blind cross-over diet intervention study. *Digestion* 81(4):252-64.

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