RIDA[®] PRECISION ABCB1
real-time PCR

Art. Nr.: Y5500
50 Reaktionen

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 -20 °C

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real-time PCR

1. Intended use

For *in-vitro* diagnostic use. RIDA[®]PRECISION ABCB1 is a multiplex real-time PCR for the direct qualitative detection of ABCB1 gene variants (SNPs: rs2032583, rs2235015) in human EDTA whole blood samples.

2. Summary and explanation of the test

Depression is one of the most serious mental disorders being estimated to affect about 350 million people worldwide.¹ Current treatment of moderate and severe depression includes antidepressants in combination with psychotherapy.² Less than one third of patients suffering from depression experience complete remission of their illness when initially treated with antidepressant agents.³ Often it may take several months to achieve a successful treatment which can cause unnecessary side effects. Approximately only two third of patients suffering from severe depression experience remission after several treatment trials.⁴

A key factor for the effectiveness of antidepressants is the ability to pass the blood-brain barrier, which is regulated by transporter molecules such as P-glycoproteins.⁵ The P-glycoprotein, encoded by the ABCB1 gene (also known as multidrug-resistance gene 1 (MDR1)), is a transmembrane protein located within the luminal membrane of endothelial cells which form the blood-brain barrier.⁶ The ABCB1 genotyping assay, developed at Max Planck Institute of Psychiatry in Munich, allows detection of specific sequence variations of the ABCB1 gene encoding the P-glycoprotein.⁷ Due to the gene polymorphisms patients can be identified which have a blood-brain barrier that is less permeable for many of the common antidepressants and which respond insufficient to therapy therefore. Previous studies from Breitenstein et al., 2014 and Sarginson et al., 2010 showed that use of the ABCB1 genotyping assay led to a positive remission outcome and enables the effective treatment for patients suffering from depression.^{6,8}

3. Test principle

RIDA®PRECISION ABCB1 is a multiplex real-time PCR for the direct, qualitative detection of ABCB1 gene variants (SNPs: rs2032583, rs2235015) in human EDTA whole blood samples.

After DNA isolation, amplification of the DNA fragments, which include the target SNPs, occur. Genotyping of the different alleles occur with specific fluorescence-marked probes (see Tab. 7).

4. Reagents provided

Tab. 1: Reagents provided (Reagents provided in the kit are sufficient for 50 determinations)

Kit Code	Reagent	Amount	Lid Color
1	Reaction Mix	1x 1100 µl	yellow
2	Taq-Polymerase	1x 40 µl	dark red
A	Control A	1x 100 µl	dark blue
B	Control B	1x 100 µl	dark green

5. Storage instructions

- Protect all reagents from light and store at -20 °C. All reagents can be used until the expiration date. After expiry the quality guarantee is no longer valid.
- Carefully thaw and fully defrost reagents before using (e.g. in a refrigerator at 2 – 8 °C).
- Reagents can sustain up to 20 freeze/thaw cycles without influencing the assay performance (e.g. after the first thawing separate it in aliquots and freeze immediately).
- During PCR preparation all the reagents should be stored cold in an appropriate way (2 – 8 °C).

6. Additional necessary reagents and necessary equipment

- The RIDA®PRECISION ABCB1 real-time PCR Assay is suitable for use with following extraction platforms and real-time PCR instruments:
- Extraction platform:
 - Maxwell® RSC Instrument (Promega)
- Real-time PCR instruments:

Roche:	LightCycler® 480II
Roche:	cobas z 480 Analyzer
Agilent Technologies:	Mx3005P
Bio-Rad:	CFX96™

If you want to use other extraction platforms or real-time PCR instruments please contact R-Biopharm at ridaprecision@r-biopharm.de.

- RIDA®GENE Color Compensation Kit IV (PG0004) for use with the LightCycler® 480II and cobas z 480 Analyzer
- Real-time PCR consumables (plates, tubes, foil)
- Centrifuge with a rotor for the reaction vials
- Vortexer
- Pipettes (0.5 – 20 µl, 20 – 200 µl, 100 – 1000 µl)
- Filter tips
- Powder-free disposal gloves

7. Precautions for users

For *in-vitro* diagnostic use.

This test must be carried out only by trained laboratory personnel. The guidelines for working in medical laboratories have to be followed. The instruction manual for the test procedure has to be followed. Do not pipet samples or reagents by mouth. Avoid contact with bruised skin or mucosal membranes. During handling reagents or samples wear appropriate safety clothing (appropriate gloves, lab coat, safety goggles) and wash your hands after finishing the test procedure. Do not smoke, eat or drink in areas where samples or reagents are being used.

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.
- Samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.
- Do not use the kit after expiration date.

All reagents and materials used have to be disposed properly after use. Please refer to the relevant national regulations for disposal.

For more details see Material Safety Data Sheets (MSDS) at www.r-biopharm.com

8. Sample collection and storage

8.1 Sample storage

The test can be carried out with human EDTA whole blood samples. Prior to DNA extraction blood samples can be stored at room temperature for up to 24 hours or at 2 – 8 °C for up to 72 hours.⁹ Microbial contamination of the samples must be avoided. Using heat-inactivated, lipaemic, haemolytic, icteric or turbid samples can lead to false results.

8.2 Sample preparation

8.2.1 DNA Isolation from EDTA whole blood

To isolate DNA from human EDTA whole blood, we recommend using a commercially available DNA extraction system (e.g. Maxwell[®] RSC Instrument (Promega)). Extract DNA according to manufacturer's instructions.

9. Test procedure

9.1 Master-Mix preparation

Calculate the total number of PCR reactions (sample and control reactions) needed. One control A and control B must be included in each assay run.

We recommend to calculate an additional volume of 10 % to compensate imprecise pipetting (see Tab. 2). Thaw, mix gently and centrifuge briefly the **Reaction Mix**, the **Taq-Polymerase**, the **Control A** and **Control B** before using. Keep reagents appropriately cold during working step (2 – 8 °C).

Tab. 2: Calculation and pipetting example for 10 reactions of the Master-Mix

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
	Total	20.0 µl	220 µl

Mix the components of the Master-Mix gently and briefly spin down.

9.2 Preparation of the PCR-Mix

Pipette 20 µl of the Master-Mix in each reaction vial (tube or plate).

Sample: Add 5 µl DNA-Extract to the pre-pipetted Master-Mix.

One control A and B must be included in each assay run.

Control A: Add 5 µl Control A as control A to the pre-pipetted Master-Mix.

Control B: Add 5 µl Control B as control B to the pre-pipetted Master-Mix.

Cover tubes or plate. Spin down and place in the real-time PCR instrument.

9.3 PCR Instrument Set-up

The PCR reaction should be started according to the PCR instrument Set-up (see Tab. 3 for LightCycler® 480II and cobas z 480 Analyzer; see Tab. 5 for Mx3005P and CFX96™). The RIDA®PRECISION ABCB1 and other R-Biopharm assays can be combined on one PCR instrument by using the Set-up universal profile (see Tab. 4 or see Tab. 6).

Tab. 3: Real-time PCR DNA profile for LightCycler® 480II and cobas z 480 Analyzer

Initial Denaturation	1 min, 95 °C
<u>Cycles</u>	45 Cycles
PCR Denaturation	10 sec, 95 °C
Annealing / Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

Tab. 4: Real-time PCR universal profile for LightCycler® 480II and cobas z 480 Analyzer

Reverse Transcription	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
<u>Cycles</u>	45 Cycles
PCR Denaturation	10 sec, 95 °C
Annealing / Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

Tab. 5: Real-time PCR DNA profile for Mx3005P and CFX96™

Initial Denaturation	1 min, 95 °C
<u>Cycles</u>	45 Cycles
PCR Denaturation	15 sec, 95 °C
Annealing / Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

Tab. 6: Real-time PCR universal profile for Mx3005P and CFX96™

Reverse Transcription	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
<u>Cycles</u>	45 Cycles
PCR Denaturation	15 sec, 95 °C
Annealing / Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

Note: Annealing and Extension occur in the same step.

9.4 Detection Channel Set-up

Tab. 7: Selection of detection channels

Real-time PCR Instrument	Detection	Detection Channel	Designation Detection Channel	Note
Roche LightCycler® 480II	rs2032583 T-allele	465/510	83T	RIDA®GENE Color Compensation Kit IV (PG0004) is required
	rs2032583 C-allele	533/580	83C	
	rs2235015 G-allele	533/610	15G	
	rs2235015 T-allele	618/660	15T	
Roche cobas z 480 Analyzer	rs2032583 T-allele	465/510	83T	RIDA®GENE Color Compensation Kit IV (PG0004) is required
	rs2032583 C-allele	540/580	83C	
	rs2235015 G-allele	540/610	15G	
	rs2235015 T-allele	610/670	15T	
Agilent Techn. Mx3005P	rs2032583 T-allele	FAM	83T	-
	rs2032583 C-allele	HEX	83C	
	rs2235015 G-allele	ROX	15G	
	rs2235015 T-allele	Cy5	15T	
Bio-Rad CFX96™	rs2032583 T-allele	FAM	83T	-
	rs2032583 C-allele	VIC	83C	
	rs2235015 G-allele	ROX	15G	
	rs2235015 T-allele	Cy5	15T	

10. Quality control

The analysis of the samples is done by the software of the used real-time PCR instrument according to manufacturer`s instructions. Positive and negative controls have to show correct results (see Tab. 8) in order to determine a VALID run.

Tab. 8: For a VALID run, the following conditions must be met:

	Gene polymorphisms			
	rs2032583 T-allele	rs2032583 C-allele	rs2235015 G-allele	rs2235015 T-allele
Designation Detection Channel	83T	83C	15G	15T
Control A	positive	negative	positive	negative
Control B	negative	positive	negative	positive

See Quality Assurance Certificate (QAC) for target Ct values of control A and B.

The PCR run must be rerun, if one of the controls, control A or B, show no amplification signal in the detection system according to specifications.

If the required criteria are not met, following items have to be checked before repeating the test:

- Expiry of the used reagents
- Functionality of the used instrumentation
- Correct performance of the test procedure

11. Result interpretation

The result interpretation is done according to table 9.

Tab. 9: Sample interpretation

	Gene polymorphisms				Result
	rs2032583 T-allele	rs2032583 C-allele	rs2235015 G-allele	rs2235015 T-allele	
Designation Detection Channel	83T	83C	15G	15T	
Control A	positive	negative	positive	negative	
Control B	negative	positive	negative	positive	
E.g. sample 1	positive	negative	positive	negative	rs2032583 TT rs2235015 GG
E.g. sample 2	positive	negative	negative	positive	rs2032583 TT rs2235015 TT
E.g. sample 3	positive	negative	positive	positive	rs2032583 TT rs2235015 GT
E.g. sample 4	negative	positive	positive	negative	rs2032583 CC rs2235015 GG
E.g. sample 5	positive	positive	positive	negative	rs2032583 TC rs2235015 GG
E.g. sample 6	negative	positive	negative	positive	rs2032583 CC rs2235015 TT
E.g. sample 7	negative	positive	positive	positive	rs2032583 CC rs2235015 GT
E.g. sample 8	positive	positive	negative	positive	rs2032583 TC rs2235015 TT
E.g. sample 9	positive	positive	positive	positive	rs2032583 TC rs2235015 GT
E.g. sample 10	negative	negative	negative	negative	Invalid (Target genes are not detectable)

A signal is evaluated negative, if the controls A and B show an amplification signal in the detection system according to specifications and the fluorescence of the sample is $\leq 10\%$ of the fluorescence of the controls A and B.

The PCR run is not evaluable, if either controls A and B or control A or B show no amplification signal in the detection system. The rerun of the PCR run is recommended.

If the sample and controls A and B show no amplification signal in the detection system (see Tab. 9, E.g. sample 10) either a component was not added to the PCR-Mix or a component was non-functional.

If controls A and B show an amplification signal in the detection system according to specifications but the sample show no amplification signal in the detection system according to specifications (see Tab. 9, E.g. sample 10) either the sample was not added to the PCR-Mix or unsuitable template DNA (Quality, PCR inhibitors) was used.

12. Limitations of the method

1. This assay is validated for human EDTA whole blood samples only.
2. The presence of PCR inhibitors may cause invalid results.
3. The result of this molecular genetic analysis gives information about the permeability of antidepressants through the blood-brain barrier only and not about its clinical effect.
4. The exact causal relation between genetic variants and the response to pharmacotherapy is not fully clarified.
6. According to Gendiagnostikgesetz (Genetic Diagnostic Act) patients get a comprehensive genetic counselling before the genetic analysis. Also genetic testing will not be conducted without the patient's active written agreement to the test.

13. Performance characteristics

13.1 Interfering substances

The RIDA®PRECISION ABCB1 real-time PCR is validated for human EDTA whole blood samples only. No interference could be detected with the following substances (see Tab. 10):

Tab. 10: List of substances and the concentration of the substances used for testing

Substances	Concentration
Heparin	15 U/ml
Cholesterol	3.0 mg/ml
Bilirubin	0.1 mg/ml
Haemoglobin	0.2 mg/ml
K ₂ EDTA	1.8 mg/ml

Explanation of Symbols

	For <i>in vitro</i> diagnostic use
	Consult instructions for use
	Lot number
	Expiry
	Store at
	Article number
	Number of test
	Date of manufacture
	Manufacturer

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