



# PROCEED<sup>TM</sup>

## PRODUCT INFORMATION SHEET

**Product Name:** PROCEEDx<sup>TM</sup> SARS-CoV-2 Negative Sample

**Catalogue Number:** VP-99-M3

**Storage:** Store this preparation at 2 °C to 8 °C.

**Description:** This product contains human fibroblast cells in a preservative solution.

**Intended Use:** The product is to be used as a negative test sample in polymerase chain reaction (PCR) based tests for SARS-CoV-2 detection. The product is to be used for nucleic acid test instrument (device) verification/validation, manufacturing QMS, educational and R&D purposes.

**NOT FOR USE IN DIAGNOSTIC PROCEDURES. FOR RESEARCH USE ONLY.**

**Handling:** Mix by micropipette aspiration of 200 µL from the sample volume and vortex the PROCEEDx<sup>TM</sup> SARS-CoV-2 Negative Sample immediately prior use. It is recommended that the product is tested according to kit Manufacturer's instructions for testing of clinical samples in a liquid sample transport medium.

**Expected Results:** Testing of the product in PCR-based assays for detection of SARS-CoV-2 should result in negative finding, including such assays where the presence of human DNA is required for the determination of a positive result.

**Hazards:** The reagent has been inactivated and should contain no infectious material: however, it should be handled as if potentially infectious. Generally accepted good laboratory practices appropriate to biological reagents should be employed. The buffer solution contains methanol (CAS no. 67-56-1) as a preservative. Please refer to the Safety Data Sheet (SDS).

**Cell line:** MRC-5

**Assistance:** If you have any questions regarding the production, testing or use of this material, please send them together with any relevant data to Microbix Technical Services, at [customer.service@microbix.com](mailto:customer.service@microbix.com) or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

## Quality Control Information

**Product Name:** PROCEEDx™ SARS-CoV-2 Negative Sample

**Lot Number:** 99M30003A

**Product Code:** VP-99-M3

qPCR Test Kit	Target	Result
CDC 2019-Novel Coronavirus (2019-nCoV) based RT-PCR test for detection of SARS-Cov-2	SARS-CoV-2	Negative
	Human DNA	Positive
Xpert Flu/RSV	Influenza A	Negative
	Influenza B	Negative
	RSV	Negative

**Inactivation:** Product was inactivated.

**Quality Control Summary:** This lot meets specifications and is released for shipment.

**Expiry date:** 29 April, 2021

Quality Assurance Signature:



Date: 17 June, 2020

