



RIDA®QUICK Rotavirus

Fast enteric human Rotavirus detection



Product specifications



Intended use:

For *in vitro* diagnostic use. The RIDA®QUICK Rotavirus Test is a quick immunochromatographic test for the qualitative determination of rotaviruses in stool samples.



Kit storage:

2 - 30 °C



Kit shelf life:

Max. 24 months



Sample type:

Untreated human stool



Analytes:

Rotavirus (red T-line)



Integrated control:

Control (blue C-line)



Relative sensitivity:

Rotavirus: 100 % (compared to ELISA)
100 % (compared to PCR)



Relative specificity:

Rotavirus: 99 % (compared to ELISA)
94.4 % (compared to PCR)



Accessory:

RIDA®QUICK Rotavirus/Adenovirus Sample diluent (Art. No. ZN1004): prefilled sample buffer tubes; Simplifies sample dilution and application on the test cassette

Why is rotavirus rapid testing relevant?

- Rotavirus infections remain among the main causes of viral gastroenteritis
- Although vaccination against rotavirus gastroenteritis is available, the disease is still a relevant cause of acute gastroenteritis
- Rapid and reliable testing improves patient care, prevents unnecessary antibiotic treatments, and controls the spread of rotavirus infections

Diagnostic solution

The RIDA®QUICK Rotavirus for professional laboratory use

- detects Rotavirus Group A in clinical stool samples
- is rapid and efficient
- provides results after 5 minutes test run
- offers reliant quality

Your benefits

Quality

- Integrated control lines to evaluate test run validity
- Easy interpretation with color-coded result lines, and interpretation marks on the test cassette
- Available compatible control to check workflow performance and result quality
- Accurate Rotavirus Group A detection¹ due to an elaborated test design and monoclonal antibodies

Efficiency

- One step immunochromatographic assay with minimized hands on time saves time and resources
- Fast results assist health care professionals in immediate decisions on patient care
- Simple workflow helps to reduce errors

Flexible

- Individual determinations for low and irregular sample throughput
- Flexible storage conditions relax the storage capacities

1: LABQUALITY data 2021/2022