


RIDA[®] GENE Color Compensation Kit II

Art. No.: PG0002

3 reactions

 -20 °C



R-Biopharm AG, An der neuen Bergstraße 17, D-64297 Darmstadt, Germany

Tel.: +49 (0) 61 51 81 02-0 / Telefax: +49 (0) 61 51 81 02-20

1. Intended use

RIDA[®]GENE Color Compensation Kit II is intended for generating a Color Compensation File for duplex and triplex real-time PCR experiments on the LightCycler[®] 1.5 and 2.0. The generated Color Compensation File can be applied to analyze multiplex real-time PCR experiments of RIDA[®]Gene real-time PCR kits on the LightCycler[®] 1.5 and 2.0.

2. Explanation of the test

In a multiplex real-time PCR, the wavelengths of light emitted by the reporter dyes may overlap, causing one channel to pick up signals (crosstalk) from a dye measured by another channel. This crosstalk of fluorescence signal can result in incorrect data unless a correction is made by using a Color Compensation File. Color compensation is used to subtract fluorescence crosstalk from a reporter dye into inappropriate channels outside of its dominant emission channel.

3. Kit components

Tab.1: Reagents provided (Reagents provided in the kit are sufficient for 3 Color Compensation experiments)

Kit Code	Reagent	Volume	Lid color
1	Blank	1x 80 µl	white
2	Dye 1	1x 80 µl	green
3	Dye 2	1x 80 µl	yellow
4	Dye 3	1x 80 µl	red

4. Storage instructions

- Protect RIDA[®]GENE Color Compensation Kit II from light and store at -20 °C.
- RIDA[®]GENE Color Compensation Kit II can be used until the expiration date printed on the label
- After expiry the quality guarantee is no longer valid.
- Carefully thaw RIDA[®]GENE Color Compensation Kit II reagents before using (e.g. in a refrigerator at 2 - 8 °C).
- During Color Compensation preparation all the reagents should be stored cold in an appropriate way (2 - 8 °C).

5. Additional equipment and materials required

- LightCycler® 1.5 or 2.0 (Roche)
- Real-time PCR consumables (LightCycler® Capillaries)
- Pipettes (0.5 – 20 µl, 20 – 200 µl, 100 – 1000 µl)
- Filter tips

6. Precautions for users

- This test must only be performed by laboratory personnel trained in molecular biology methods.
- Strictly follow the working instructions.
- When handling samples, wear disposable gloves. After finishing the test, wash your hands.
- Do not smoke, eat or drink in areas where samples or test reagents are being used.
- Do not use the kit after the expiration date.

7. Protocol for creating a Color Compensation File on the LightCycler® 1.5

7.1 Preparation of the Color Compensation

For a color compensation experiment it is necessary to pipette 20 µl of each dye and also of the background (blank) into a LightCycler® capillary. Add 1 µl of Dye 1 to the pre-pipetted 20 µl of Dye 3 in LightCycler® capillary 4 (see Tab.2).

Thaw, mix gently and centrifuge briefly the reagents before use. Keep reagents appropriately cold during working step (2 - 8 °C).

Tab.2: Preparation of the Color Compensation


Kit Code	Reagent	Volume per reaction			
		Capillary Position 1	Capillary Position 2	Capillary Position 3	Capillary Position 4
1	Blank	20 µl	-	-	-
2	Dye 1	-	20 µl	-	1 µl
3	Dye 2	-	-	20 µl	-
4	Dye 3	-	-	-	20 µl

Cap and briefly spin down the LightCycler® capillaries after pipetting the reagents.


Note: Do not change the order of the LightCycler® capillaries on the LightCycler®-Carousel.

7.2 LightCycler® 1.5 Set-up

1. Open the LightCycler® 1.5 software. Click the “**Run**” button to open a new “**LightCycler Experiment**”.



2. The following window opens.




3. Program the LightCycler[®] according to the real-time PCR profile (see Tab.3). Click the “Add” button to program the 4 protocol steps.


Tab.3: LightCycler[®] real-time PCR profile

Program	Cycles / Analysis Mode	Temperature targets			
		Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp rate [°c/s]
Initial Denat.	1 / none	95	none	00:00:05	20
Cycling	5 / Quantification	95	none	00:00:05	20
		60	single	00:00:05	20
TM-Analyse	1 / Color Compensation	95	none	00:00:05	20
		50	none	00:00:05	20
		70	continuous		0.2 (Acquisitions per °C = 1)
Cooling 40	1 / none	40	none	00:00:10	20

Note: Ensure the correct setting of the number of “Cycles” and of the “Analysis Mode”



4. The final screen should look like the screenshot below.



5. Place the LightCycler[®]-Carousel with the capillaries in the LightCycler[®] 1.5. Click the “Run” button and save the experiment in the proper folder.

The following window will open to program the PCR layout. Enter in the “Maximum Position” dialog field the number of capillaries. Click the “Done” button to start the experiment.


File Edit Help

LC Carousel: Please edit the sample data for this run.

#	Sample Name	Type	Replicate of	Concentration	Notes	#	Sample Name	Type	Replicate of	Concentration	Notes
1	Sample 1	Unknown	0	0,00E+0		17	Sample 17	Unknown	0	0,00E+0	
2	Sample 2	Unknown	0	0,00E+0		18	Sample 18	Unknown	0	0,00E+0	
3	Sample 3	Unknown	0	0,00E+0		19	Sample 19	Unknown	0	0,00E+0	
4	Sample 4	Unknown	0	0,00E+0		20	Sample 20	Unknown	0	0,00E+0	
5	Sample 5	Unknown	0	0,00E+0		21	Sample 21	Unknown	0	0,00E+0	
6	Sample 6	Unknown	0	0,00E+0		22	Sample 22	Unknown	0	0,00E+0	
7	Sample 7	Unknown	0	0,00E+0		23	Sample 23	Unknown	0	0,00E+0	
8	Sample 8	Unknown	0	0,00E+0		24	Sample 24	Unknown	0	0,00E+0	
9	Sample 9	Unknown	0	0,00E+0		25	Sample 25	Unknown	0	0,00E+0	
10	Sample 10	Unknown	0	0,00E+0		26	Sample 26	Unknown	0	0,00E+0	
11	Sample 11	Unknown	0	0,00E+0		27	Sample 27	Unknown	0	0,00E+0	
12	Sample 12	Unknown	0	0,00E+0		28	Sample 28	Unknown	0	0,00E+0	
13	Sample 13	Unknown	0	0,00E+0		29	Sample 29	Unknown	0	0,00E+0	
14	Sample 14	Unknown	0	0,00E+0		30	Sample 30	Unknown	0	0,00E+0	
15	Sample 15	Unknown	0	0,00E+0		31	Sample 31	Unknown	0	0,00E+0	
16	Sample 16	Unknown	0	0,00E+0		32	Sample 32	Unknown	0	0,00E+0	


Seek Temperature: Maximum Position: Done Enter Samples later Clear Sample List Default Sample List Concentration Units:

1. After the experiment ends the following window opens.



2. Click **“Select a Program”** button to select segment 3 (TM-Analyse) for Color Compensation.

Ensure that the green cursors are placed around the start and end of the melting curve segment. Then select **“Calibration”** from the Color Compensation pull-down menu.



Save the Color Compensation File. The created Color Compensation File can now be applied for LightCycler® experiments.

To apply Color Compensation for a multiplex PCR experiment click the **“Select CC Data”** button to select and import the appropriate Color Compensation File. Select the appropriate channel (e.g. F1 or F2) and click the **“Quantification”** button to analyze the experiment.

Note: The Color Compensation File is specific for every LightCycler® instrument. A new Color Compensation File has to be created if the LightCycler® instrument is changed or after the optical system has been repaired.

8. Protocol for creating a Color Compensation File on the LightCycler® 2.0

8.1. Preparation of the Color Compensation

For a color compensation experiment it is necessary to pipette 20 µl of each dye and also of the background (blank) into a LightCycler® capillary. (see Tab.4).

Thaw, mix gently and centrifuge briefly the reagent before using. Keep reagents appropriately cold during working step (2 - 8 °C).

Tab.4: Preparation of the Color Compensation


Kit Code	Reagent	Volume per reaction	Capillary position
1	Blank	20 µl	1
2	Dye 1	20 µl	2
3	Dye 2	20 µl	3
4	Dye 3	20 µl	4

Cap and briefly spin down the LightCycler® capillaries after pipetting the reagents.


Note: Do not change the order of the LightCycler® capillaries on the LightCycler® Carousel.

8.2 LightCycler® 2.0 Set-up

1. Open the LightCycler® 2.0 software. Click the “**File**” menu, and select “**New**” to open the “**Create New Object**” dialog box. Double-click the “**LightCycler Experiment**” icon to open a new LightCycler experiment file.



2. The following window opens.



3. Program the the LightCycler[®] according to the Set-up (see Tab. 5) and the protocol steps of the real-time PCR profile (see Tab. 6)


Tab.5: LightCycler[®] Set-up

Parameter	Setting
Default Channel	530
Seek Temperature	70 °C
Max. Seek Pos.	4
Instrument Type	6 Ch.
Capillary Size	20µl


Tab.6: LightCycler[®] real-time PCR profile

Program	Cycles / Analysis Mode	Temperature targets			
		Target [°C]	Acquisition Mode	Hold [hh:mm:ss]	Ramp rate [°c/s]
Initial Denat.	1 / none	95	none	00:00:05	20
Cycling	5 / Quantification	95	none	00:00:05	20
		60	single	00:00:05	20
TM-Analyse	1 / Color Compensation	95	none	00:00:05	20
		50	none	00:00:05	20
		70	continuous		0.2 (Acquisitions per °C = 1)
Cooling 40	1 / none	40	none	00:00:10	20


Note: Ensure the correct setting of the number of “Cycles” and of the “Analysis Mode”



4. The final screen should look like the screenshot below.




5. Click the “**Samples**” icon on the left-hand side of the window. Then click the “**Analysis Type**” menu, and select “**Color Compensation**”.



6. Define in the “**Dominant Channel**” dialog fields the dominant channel for the selected dyes for each position (see Tab.7). Select for the reagent “**Blank**” the dominant channel “**Water**”. It is not necessary to specify a “**Sample Name**”.

Tab.7: Detection Channel Set-up

Reagent	Detection Channel
Blank	Water
Dye1	530
Dye2	560
Dye3	705



Place the LightCycler[®]-Carousel with the capillaries in the LightCycler[®] 2.0. Click the “Run” icon and save the experiment in the proper folder.

8.3 Evaluation and creation of a Color Compensation File

After the experiment ends, click the “Analysis” button to open the “Create New Analysis” dialog box. Under “Other Methods”, select “Color Compensation”. Then click “OK”

In the “Color Compensation” window which appears, click the “Save CC Object” button to save the Color Compensation File in the “CCC” folder. The created Color Compensation File can now be applied for LightCycler[®] experiments.

To apply Color Compensation for a multiplex PCR experiment open the experiment and click the “Analysis” button. Click the “Color Compensation (Off)” button, and select the appropriate Color Compensation File. The “Color Compensation (Off)” button switches to “Color Compensation (On)” to confirm that color compensation is active. The multiplex real-time PCR experiment can now be analyzed.

Note: The Color Compensation File is specific for every LightCycler[®] instrument. A new Color Compensation File has to be created if the LightCycler[®] instrument is changed or after the optical system has been repaired.