## RIDA® GENE Gut Balance

**REF** PG0105

R-Biopharm AG, An der neuen Bergstraße 17, 64297 Darmstadt, Germany  
Phone: +49 (0) 61 51 81 02-0 / Fax: +49 (0) 61 51 81 02-20

Deutsch .....	3
English.....	27
Español .....	49
Français.....	71
Italiano.....	93

# RIDA<sup>®</sup>GENE Gut Balance

**REF** PG0105

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### Testspezifische Symbole

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## 16. Literatur

1. Matsuki. Real-time PCR Analysis of human intestinal microflora with 16S rRNA-gene-targeted Genus and species-specific primers. XI Congreso Nacional de Biotecnología y Bioingeniería.
2. Mai V *et al.* Associations between dietary habits and body mass index with gut microbiota composition and fecal water genotoxicity: an observational study in African American and Caucasian American volunteers. *Nutr. Journal* 2009, 8: 49 – 59.
3. Turnbaugh P *et al.* A core gut microbiome in obese and lean twins. *Nature* 2009, 457(7228): 480 – 484.
4. Vaarala O. Gut Microbiota and Type 1 Diabetes. *Rev. Diab. Stud.* 2013, 9(4): 251 - 259.

**English**

## **RIDA<sup>®</sup>GENE Gut Balance**

**REF** PG0105

### **1. Intended use**

For *in vitro* diagnostic use. RIDA<sup>®</sup>GENE Gut Balance is a multiplex real-time PCR for the direct qualitative and quantitative detection of *Bacteroides*- and Cluster XIVa DNA from human stool samples.<sup>1</sup>

### **2. Summary and explanation of the test**

90 % of the normal human gut flora is populated by two phylogenetic groups which exist in a symbiotic balance. *Bacteroides* are anaerobic, gram-negative bacteria which are part of the normal gut flora of the intestinal tract. In the large intestine, approximately 10<sup>11</sup> *Bacteroides*/g stool exist and are therefore the dominant bacteria in terms of numbers. The second phylogenetic group is the *firmicutes*. *Clostridium* Cluster XIVa are a class of *firmicutes* to which, besides others, *Eubacterium* spp. and *Roseburia* spp. belong.

Different sources associate a disbalance of the composition of the gut flora (dysbiosis) with obesity.<sup>2,3</sup> Here, a decreased number of *Bacteroides* corresponds to presence of obesity whereas at the same time an increasing number of *Eubacterium rectale* was detected in patients with obesity.<sup>3,4</sup>

### **3. Test principle**

RIDA<sup>®</sup>GENE Gut Balance is a multiplex real-time PCR for the direct, qualitative and quantitative detection of *Bacteroides*- and Cluster XIVa-DNA in human stool samples. After DNA isolation, amplification of the gene fragment (if present) specific for *Bacteroides* and Cluster XIVa (16S-rRNA) occurs. The amplified targets are detected with hydrolysis probes, which are labeled at one end with a quencher and at the other end with a fluorescent reporter dye (fluorophore). In the presence of a target the probes hybridize to the amplicons. During the extension step, the **Taq-Polymerase** breaks the reporter-quencher proximity. The reporter emits a fluorescent signal, which is detected by the optical unit of a real-time PCR instrument. The fluorescence signal increases with the amount of formed amplicons. With the standards, **Standard A**, **Standard B** and **Standard C**, included in the kit, it is possible to quantify the results. The identified DNA amount in the sample (copies/reaction) is

converted into the concentration unit cells / g stool with a correction factor (K, see also table 12). The RIDA®GENE Gut Balance multiplex real-time PCR kit contains an **Internal Control DNA** (ICD) that detects PCR inhibition, monitors reagent integrity and confirms that nucleic acid extraction was sufficient.

#### 4. Reagents provided

**Tab. 1: Reagents provided (Reagents provided in the kit are sufficient for 100 determinations)**

Kit Code	Reagent	Amount		Lid Color
1	Reaction Mix	2x	1050 µl	yellow
2	Taq-Polymerase	1x	80 µl	red
D	Internal Control DNA	2x	1700 µl	orange
N	No Template Control	1x	450 µl	white
P	Positive Control	1x	200 µl	blue
10 <sup>2</sup>	Standard A	1x	100 µl	dark blue
10 <sup>4</sup>	Standard B	1x	100 µl	dark blue
10 <sup>6</sup>	Standard C	1x	100 µl	dark blue

#### 5. Storage instructions

- Protect all reagents from light and store at -20 °C. All reagents can be used unopened until the expiration date. After expiry the quality guarantee is no longer valid.
- Carefully thaw reagents before using (e.g. in a refrigerator at 2 - 8 °C).
- Reagents can sustain up to **20 freeze/thaw cycles** without influencing the assay performance (e.g. after the first thawing separate it in aliquots and freeze immediately).
- During PCR preparation all the reagents should be stored cold in an appropriate way (2 - 8 °C).

## 6. Additional necessary reagents and necessary equipment

The RIDA®GENE Gut Balance multiplex real-time PCR assay is suitable for use with following extraction platforms and real-time PCR instruments:

**Tab. 2:** Necessary equipment

Extraction platform	
R-Biopharm	RIDA® Xtract
Promega	Maxwell® RSC
Real-time PCR instrument:	
Roche	LightCycler® 480II
Agilent Technologies	Mx3005P
Applied Biosystems	ABI 7500
Bio-Rad	CFX96™
QIAGEN	Rotor-Gene Q

**Note: Only use 0.1 ml tubes on the Rotor-Gene Q (QIAGEN).**

If you want to use other extraction platforms or real-time PCR instruments please contact R-Biopharm at [mdx@r-biopharm.de](mailto:mdx@r-biopharm.de).

- RIDA®GENE Color Compensation Kit IV (PG0004) for use with the LightCycler® 480II

- Real-time PCR consumables (plates, tubes, foil)
- Centrifuge with a rotor for the reaction vials
- Vortexer
- Pipettes (0.5 – 20 µl, 20 – 200 µl, 100 – 1000 µl)
- Filter tips
- Powder-free disposal gloves
- PCR water (nuclease-free)

## 7. Precautions for users

For *in-vitro* diagnostic use.

This test must only be carried out by trained laboratory personnel. The guidelines for working in medical laboratories have to be followed. The instruction manual for the test procedure has to be followed. Do not pipet samples or reagents by mouth.

Avoid contact with bruised skin or mucosal membranes. During handling reagents or samples, wear appropriate safety clothing (appropriate gloves, lab coat, safety goggles) and wash your hands after finishing the test procedure. Do not smoke, eat or drink in areas where samples or reagents are being used.

- Extraction, PCR preparation and the PCR run should be separated in different rooms to avoid cross-contaminations.

- Samples must be treated as potentially infectious as well as all reagents and materials being exposed to the samples and have to be handled according to the national safety regulations.

- Do not use the kit after the expiration date.

All reagents and materials used have to be disposed properly after use.

Please refer to the relevant national regulations for disposal.

For more details see Safety Data Sheets (SDS) at [www.r-biopharm.com](http://www.r-biopharm.com)

## 8. Collection and storage of samples

### 8.1 Sample preparation from stool samples

For DNA isolation of human stool samples, use a commercially available DNA isolation kit (e.g. RIDA® Xtract (R-Biopharm)) or DNA extraction system (e.g. Maxwell® RSC (Promega)). Extract DNA according to the manufacturer's instructions.

We recommend to dilute the stool samples before extraction 1:3 with water. Vortex the diluted stool sample intensely and centrifuge at 1000 x g for 30 sec. Use from the supernatant the appropriate volume according to the manufacturer's instruction.

The RIDA®GENE Gut Balance assay contains an Internal Control DNA that detects PCR inhibition, monitors reagent integrity and confirms that nucleic acid extraction was sufficient. The Internal Control DNA can either be used as PCR inhibition control or as extraction control for the sample preparation procedure and as PCR inhibition control.

If the Internal Control DNA is used only as a PCR inhibition control, 1 µl of the Internal Control DNA should be added to the Master-Mix (see Tab.4).

If the Internal Control DNA is used as an extraction control for the sample preparation procedure **and** as PCR inhibition control, 20 µl of the Internal Control DNA has to be added during extraction procedure. The Internal Control DNA should always be added to the specimen-lysis buffer mixture and must **not** be added directly to the specimen. We also recommend to add 1 µl of the Internal Control DNA to the negative control and positive control PCR Mix.

## 9. Test procedure

### 9.1 Master-Mix preparation

Calculate the total number of PCR reactions (sample and control reactions) needed. One positive control and one negative control must be included in each assay run.

We recommend to calculate an additional volume of 10 % to compensate imprecise pipetting (see Tab. 3, Tab. 4). Thaw, mix gently and centrifuge briefly the **Reaction Mix**, the **Taq-Polymerase**, the **Positive Control**, the **No Template Control**, the **Internal Control DNA** and **Standard A**, **Standard B** and **Standard C** before using. Keep reagents appropriately cold during working step (2 - 8 °C).

**Tab. 3:** Calculation and pipetting example for 10 reactions of the Master-Mix (ICD as extraction and PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
	<b>Total</b>	<b>20 µl</b>	<b>220 µl</b>

Mix the components of the Master-Mix gently and briefly spin down.

**Tab. 4:** Calculation and pipetting example for 10 reactions of the Master-Mix (ICD only as PCR inhibition control)

Kit code	Master-Mix components	Volume per reaction	10 reactions (10 % extra)
1	Reaction Mix	19.3 µl	212.3 µl
2	Taq-Polymerase	0.7 µl	7.7 µl
D	Internal Control DNA	1.0 µl	11 µl
	<b>Total</b>	<b>21.0 µl</b>	<b>231.0 µl</b>

Mix the components of the Master-Mix gently and briefly spin down.

## 9.2 Preparation of the PCR-Mix

Pipette 20 µl of the Master-Mix in each reaction vial (tube or plate).

**Negative control:** Add 5 µl **No Template Control** to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the negative control.

**Sample:** Add 5 µl eluate to the pre-pipetted Master-Mix.

**Positive control:** Add 5 µl **Positive Control** to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the positive control.

**Standard (A, B, C):** Add 5 µl **Standard** (A, B, C) to the pre-pipetted Master-Mix.

**Note:** If the **Internal Control DNA** is used as extraction control for the sample preparation procedure and as PCR inhibition control, we recommend to add 1 µl of the **Internal Control DNA** to the PCR-Mix of the standards.

**Note:** Using the following cyclers requires to include a standard curve in each run: ABI 7500 (Applied Biosystems) and CFX96™ (Bio-Rad).

**For all other cyclers, only one sample of the standard curve (**Standard B**) has to be included in the experimental set-up as calibrator for each new real-time PCR run. Here, the application of a standard curve is only required to run once per lot number.**

Cover tubes or plate. Spin down and place in the real-time PCR instrument. The PCR reaction should be started according to the PCR instrument set-up (see Tab. 5, Tab. 6, Tab. 7, Tab. 8).

## 9.3 PCR instrument set-up

### 9.3.1 DNA real-time PCR profile

**Tab. 5:** DNA real-time PCR profile for LightCycler® 480II and Rotor-Gene Q

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and extension occur in the same step.

**Tab. 6:** DNA real-time PCR profile for Mx3005P, ABI7500, CFX96™

Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and extension occur in the same step.

**Note:** The total copy number per reaction of **Standard A**, **Standard B** and **Standard C** has to be typed in into the Setup File of the software program of the respective real-time PCR cycler. A total volume of 5 µl DNA is used resulting in following concentrations:

**Standard A:**  $5 \times 10^2$  copies/reaction

**Standard B:**  $5 \times 10^4$  copies/reaction

**Standard C:**  $5 \times 10^6$  copies/reaction

**Note:** The standard curve can be saved on the real-time PCR cycler for each parameter. Apart from the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler, the standard curve is only required to run once per lot number. Using the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler requires to include a standard curve in each run. For all other cyclers, only one sample of the standard curve (**Standard B**) has to be included in the experimental set-up as calibrator for each new real-time PCR run.

### 9.3.2 Universal real-time PCR profile

**Note:** The universal real-time PCR profile should only be used for DNA assays when combining RIDA®GENE DNA and RNA real-time PCR assays in one run.

**Tab. 7:** Universal real-time PCR profile for LightCycler® 480II

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	10 sec, 95 °C
Annealing/Extension	15 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Tab. 8:** Universal real-time PCR profile for Mx3005P, ABI7500, CFX96™ and Rotor-Gene Q

<u>Reverse Transcription</u>	10 min, 58 °C
Initial Denaturation	1 min, 95 °C
Cycles	45 Cycles
<u>PCR</u> Denaturation	15 sec, 95 °C
Annealing/Extension	30 sec, 60 °C
Temperature Transition Rate / Ramp Rate	Maximum

**Note:** Annealing and Extension occur in the same step.

**Note:** The total copy number per reaction of **Standard A**, **Standard B** and **Standard C** has to be typed in into the Setup File of the software program of the respective real-time PCR cycler. A total volume of 5 µl DNA is used resulting in following concentrations:

**Standard A:** 5 x 10<sup>2</sup> copies/reaction

**Standard B:** 5 x 10<sup>4</sup> copies/reaction

**Standard C:** 5 x 10<sup>6</sup> copies/reaction

**Note:** The standard curve can be saved on the real-time PCR cycler for each parameter. Apart from the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler, the standard curve is only required to run once per lot number. Using the ABI 7500 (Applied Biosystems) and the CFX96™ (Bio-Rad) cycler requires to include a standard curve in each run.

For all other cyclers, only one sample of the standard curve (**Standard B**) has to be included in the experimental set-up as calibrator for each new real-time PCR run.

#### 9.4 Detection channel set-up

**Tab. 9:** Selection of appropriate detection channels

Real-time PCR instrument	Detection	Detection channel	Note
Roche LightCycler® 480II	<i>Bacteroides</i>	465/510	<b>RIDA® GENE Color Compensation Kit IV (PG0004) is required</b>
	ICD	533/580	
	Cluster XIVa	618/660	
Agilent Technologies Mx3005P	<i>Bacteroides</i>	FAM	<b>Check that reference dye is none</b>
	ICD	HEX	
	Cluster XIVa	Cy5	
ABI 7500	<i>Bacteroides</i>	FAM	<b>Check that passive reference option ROX is none</b>
	ICD	VIC	
	Cluster XIVa	Cy5	
Bio-Rad CFX96™	<i>Bacteroides</i>	FAM	-
	ICD	VIC	
	Cluster XIVa	Cy5	
Qiagen Rotor-Gene Q	<i>Bacteroides</i>	Green	<b>The gain settings have to be set to 5, according to the default settings</b>
	ICD	Yellow	
	Cluster XIVa	Red	

## 10. Quality control

The analysis of the samples is done by the software of the used real-time PCR instrument according to the manufacturer's instructions. Positive control and negative control have to show correct results (see Table 10, Fig. 1, Fig. 2) in order to determine a valid run.

The **Positive Control** has a concentration of  $10^3$  copies/ $\mu$ l. In each PCR run it is used in a total amount of  $5 \times 10^3$  copies.

**Tab. 10:** For a valid run, the following conditions must be met:

Sample	Assay result	ICD Ct	Target Ct
Positive control	Positive	NA *1	See Quality Assurance Certificate
Negative control	Negative	Ct > 20	Not detectable


*\*1 No Ct value is required for the ICD to make a positive call for the positive control.*

If the positive control is not positive within the specified Ct range but the negative control is valid, prepare all new reactions including the controls.


If the negative control is not negative but the positive control is valid, prepare all new reactions including the controls.

If the required criteria are not met, following items have to be checked before repeating the test:

- Expiry of the used reagents
- Functionality of the used instrumentation
- Correct performance of the test procedure

**Fig. 1:** Correct run of the positive control (red) and negative control (green) (*Bacteroides*) on the LightCycler® 480II

**Fig. 2:** Correct run of the positive control (red) and negative control (green) (Cluster XIVa) on the LightCycler® 480II

## 10.1 Validity of quantitative detection

For the validity of a quantitative diagnostic test run, all control conditions of a valid qualitative diagnostic test run must be met. Furthermore, for accurate quantification results a valid standard curve has to be generated. For a valid quantitative diagnostic test run, the following control parameter values of the standard curve should be achieved.

	Control parameter	Valid value
<b>Roche LightCycler® 2.0</b>	Efficiency	1,9 – 2,1
<b>Roche LightCycler® 480II</b>	Efficiency	1,9 – 2,1
	Slope	-3,1 – -3,6
<b>Agilent Techn. Mx3005P</b>	Rsq	> 0,98
	Y	-3,1 – -3,6
<b>ABI 7500</b>	R <sup>2</sup>	> 0,98
	Slope	-3,1 – -3,6
<b>Bio-Rad CFX96™</b>	R <sup>2</sup>	> 0,98
	Slope	-3,1 - -3,6
<b>Qiagen Rotor-Gene Q</b>	R <sup>2</sup>	> 0,98
	M	-3,1 – -3,6

## 11. Result interpretation

The result interpretation is done according to Table 11.

**Tab. 11:** Sample interpretation

Target genes			
<i>Bacteroides</i>	Cluster XIV	ICD	Result
positive	negative	positive / negative	<b><i>Bacteroides</i> detected*</b>
negative	positive	positive / negative	<b>Cluster XIVa detected*</b>
positive	positive	positive / negative	<b><i>Bacteroides</i> and Cluster XIVa detected</b>
negative	negative	positive	<b>Target genes not detected*</b>
negative	negative	negative	<b>Invalid</b>

A sample is evaluated positive, if the sample DNA and the Internal Control DNA show an amplification signal in the detection system.

A sample is also evaluated positive, if the sample DNA shows an amplification signal but none for the Internal Control DNA in the detection system. The detection of the internal amplification control is not necessary because high concentrations of the amplicon can cause a weak or absent signal of the Internal Control DNA.

A sample is evaluated negative, if the sample DNA shows no amplification signal, but an amplification signal for Internal Control DNA in the detection system. An inhibition of the PCR reaction can be excluded by the detection of the Internal Control DNA.

A sample is invalid, if the sample DNA and Internal Control DNA show no amplification signal in the detection system. The sample contains a PCR inhibitor. The extracted sample needs to be further diluted with PCR water (1:10) and re-amplified, or the isolation and purification of the sample has to be improved.

**\*Note: A double-negative result for *Bacteroides* and Cluster XIVa DNA is unlikely since both bacterial groups are human commensal bacteria. Accordingly, this is also valid for a negative result for only one of the both bacteria groups.** If a double-negative result occurs for *Bacteroides* and Cluster XIVa DNA, it is likely that, upon use of the ICD as inhibition control, the sample extraction was not successful. If a double-negative result occurs for *Bacteroides* and Cluster XIVa DNA, it is recommended to improve isolation and purification of the sample and repeat amplification of the sample.

## 11.1 Quantification of samples

To quantify *Bacteroides* and Cluster XIVa positive samples, a standard curve with the **Standard A**, **Standard B** and **Standard C** has to be performed separately. The standard curve measurement has to be saved separately. However, the same standard curve measurement can be used in all runs with products from the same lot number by importing the saved experiment.

**Note: This is not valid for the following cyclers: ABI 7500 (Applied Biosystems), and CFX96™ (Bio-Rad). Here, a standard curve has to be measured with each run. For all other cyclers, one sample of the standard curve (Standard B) has to be included in the experimental set-up as calibrator for each new real-time PCR run.**

To quantify *Bacteroides* and Cluster XIVa positive samples, all standard samples (A, B and C), the positive and the negative control as well as the unknown samples to be quantified, have to be selected and analyzed according to the instructions of the cycler manufacturer.

With the quantitative RIDA®GENE Gut Balance multiplex real-time PCR the amount of DNA in copies/reaction of the parameter is calculated. The conversion into the concentration unit cells/g stool sample is done with a correction factor K and takes into account the dilutions of the extraction procedure (dependent on the extraction kit used) and the PCR set-up as well as the number of target sequences in the whole genome.

The conversion of the result of the quantitative RIDA®GENE Gut Balance multiplex real-time PCR in cells/g stool is calculated with following formula:

$$C \text{ [cells/g stool]} = c \text{ [copies/reaction]} \times K$$

C [cells/g stool] - bacterial concentration of sample in cells/g stool

c [copies/reaction] - DNA concentration in PCR reaction  
(result of quantitative PCR)

K - correction factor

For the calculation of the correction factor, following information has to be considered:

- Sample dilution
- Starting volume of sample for DNA extraction
- DNA extract from total eluate used for PCR reaction
- Number of target sequence in the whole genome

**Tab. 12:** Example calculation of correction factor using Maxwell® RSC (Promega) for sample preparation of a 1:3 diluted sample

Description	Factor
Sample dilution 1:3 before extraction	x 3
300 µl sample for extraction*	x 3.33
5 µl DNA extract into PCR reaction**	x 20
a. Target sequence contained 6x in total <i>Bacteroides</i> genome or	a. x 0,167 ( <i>Bacteroides</i> )
b. Target sequence contained 5x in total Cluster XIVa genome	b. x 0,2 (Cluster XIVa)
Correction factor K for <i>Bacteroides</i>	<b>0,33 x 10<sup>2</sup></b>
Correction factor K for Cluster XIVa	<b>0,40 x 10<sup>2</sup></b>

\* Result corresponds to 1 g stool

\*\* Corresponding to a total eluate of 100 µl (= 1/20)

**Note:** For further information on quantification of please contact [mdx@r-biopharm.de](mailto:mdx@r-biopharm.de).

## 12. Limitations of the method


1. The result of molecular analysis should not lead to the diagnosis, but always be considered in the context of medical history and symptoms of the patient.
2. This assay is only validated for human stool samples.
3. Inappropriate specimen collection, transport, storage and processing or a pathogen load in the specimen below the analytical sensitivity can result in false negative results.
4. The presence of PCR inhibitors may cause invalid results.
5. Mutations or polymorphisms in primer or probe binding regions may affect detection of new variants resulting in a false negative result with the RIDA®GENE Gut Balance assay.
6. As with all PCR based *in vitro* diagnostic tests, extremely low levels of target below the limit of detection (LoD) may be detected, but results may not be reproducible.
7. A positive test result does not necessarily indicate the presence of viable organisms. However, a positive result is indicative for the presence of the target genes (16S-rRNA).

## 13. Performance characteristics


### 13.1 Analytical sensitivity

The RIDA®GENE Gut Balance multiplex real-time PCR has a detection limit of  $\geq 10$  DNA copies per reaction for *Bacteroides* and Cluster XIVa.

The following figures 3 and 4 show dilution series of *Bacteroides* and Cluster XIVa ( $10^6 - 10^2$  DNA copies per  $\mu\text{l}$ ) on the LightCycler® 480II.

**Fig. 3:** Dilution series *Bacteroides* ( $10^6 - 10^2$  DNA copies per  $\mu\text{l}$ ) on the LightCycler® 480II

**Fig. 4:** Dilution series Cluster XIVa ( $10^6 - 10^2$  DNA copies per  $\mu\text{l}$ ) on the LightCycler<sup>®</sup> 480II

The detection limit of the whole procedure depends on the sample matrix, DNA extraction and DNA concentration.

### 13.2 Analytical specificity

The RIDA<sup>®</sup>GENE Gut Balance multiplex real-time PCR is specific for *Bacteroides* and Cluster XIVa. No cross-reaction could be detected for the following species (see Tab. 13). All strains tested were detected by RIDAGENE Gut Balance multiplex real-time PCR and sequence matching (\*).

**Tab. 13:** Cross-reactivity testing

Adenovirus 1, human, strain Adenoid 71	-	<i>Campylobacter lari</i> subsp. <i>lari</i>	-	<i>E. coli</i> (O157:H7)	-	Norovirus GG* I	-
Adenovirus 7, human, strain Gomen	-	<i>Campylobacter upsaliensis</i>	-	<i>E. coli</i> (O26:H-)	-	Norovirus GG II*	-
Adenovirus 40, human, strain Dugan	-	<i>Candida albicans</i>	-	<i>E. coli</i> (O6)	-	<i>Proteus vulgaris</i>	-
Adenovirus 41, human, strain Tak	-	<i>Citrobacter freundii</i>	-	<i>Entamoeba histolytica</i>	-	<i>Pseudomonas aeruginosa</i>	-
<i>Aeromonas hydrophila</i>	-	<i>Clostridium bifermentans</i>	-	<i>Enterobacter cloacae</i>	-	Rotavirus	-
<i>Arcobacter butzleri</i>	-	<i>Clostridium difficile</i>	-	<i>Enterococcus faecalis</i>	-	<i>Salmonella enteritidis</i>	-
Astrovirus Type 2	-	<i>Clostridium novyi</i>	-	<i>Enterococcus faecium</i>	-	<i>Salmonella typhimurium</i>	-
Astrovirus Type 8	-	<i>Clostridium perfringens</i>	-	<i>Giardia intestinalis</i> Portland 1	-	<i>Serratia liquefaciens</i>	-
<i>Bacillus cereus</i>	-	<i>Clostridium septicum</i>	-	<i>Giardia intestinalis</i> WB Clone C6	-	<i>Shigella flexneri</i>	-
<i>Bifidobacterium bifidum</i> *	-	<i>Clostridium sordellii</i>	-	<i>Giardia lamblia</i>	-	<i>Staphylococcus aureus</i>	-
<i>Campylobacter coli</i>	-	<i>Clostridium sporogenes</i>	-	<i>Klebsiella oxytoca</i>	-	<i>Staphylococcus epidermidis</i>	-
<i>Campylobacter fetus</i> subsp. <i>fetus</i>	-	<i>Cryptosporidium muris</i>	-	<i>Lactobacillus ruminis</i>	-	<i>Vibrio parahaemolyticus</i>	-
<i>Campylobacter jejuni</i>	-	<i>Cryptosporidium parvum</i> *	-	<i>Lactobacillus salivaris</i>	-	<i>Yersinia enterocolitica</i>	-

### 13.3 Analytical reactivity

The reactivity of the RIDA®GENE Gut Balance multiplex real-time PCR was evaluated exemplarily with different *Bacteroides* and Cluster XIVa strains (see Tab. 14). All tested strains were detected by the RIDA®GENE Gut Balance multiplex real-time PCR assay or by sequence alignment (\*).

**Tab. 14:** Analytical reactivity testing

<b>Bacteroides</b>					
<i>Bacteroides acidifaciens</i> *	+	<i>Bacteroides fragilis</i>	+	<i>Bacteroides rodentium</i> *	+
<i>Bacteroides caccae</i> *	+	<i>Bacteroides gallinarum</i> *	+	<i>Bacteroides stercorisoris</i> *	+
<i>Bacteroides caecimuris</i> *	+	<i>Bacteroides helcogenes</i> *	+	<i>Bacteroides stercoris</i> *	+
<i>Bacteroides cellulosilyticus</i> *	+	<i>Bacteroides intestinalis</i> *	+	<i>Bacteroides thetaiotaomicron</i> *	+
<i>Bacteroides clarus</i> *	+	<i>Bacteroides mediterraneensis</i> *	+	<i>Bacteroides timonensis</i> *	+
<i>Bacteroides dorei</i> *	+	<i>Bacteroides nordii</i> *	+	<i>Bacteroides uniformis</i> *	+
<i>Bacteroides eggerthii</i> *	+	<i>Bacteroides oleiciplenus</i> *	+	<i>Bacteroides xylanisolvens</i> *	+
<i>Bacteroides fingoldii</i> *	+	<i>Bacteroides ovatus</i> *	+		
<b>Cluster XIVa</b>					
<i>Acetivomaculum ruminis</i>	+	<i>Clostridium aminophilum</i>	+	<i>Eubacterium rectale</i>	+
<i>Clostridium aerotolerans</i>	+	<i>Clostridium clostridioforme</i>	+		

## 14. Version history

Version number	Chapter and designation
2019-07-15	General revision 3. Test principle 4. Reagents provided 5. Storage instructions 6. Additional necessary reagents and necessary equipment 7. Precautions for users 8. Collection and storage of samples 9. Test procedure 10. Quality control 11. Result interpretation 13. Performance characteristics 14. Version history 15. Explanation of symbols

## 15. Explanation of symbols

### General symbols

	For in vitro diagnostic use
	Consult instructions for use
	Lot number
	Expiry
	Store at
	Article number
	Number of tests
	Date of manufacture
	Manufacturer

### Testspecific symbols

Not applicable

## 16. Literature

1. Matsuki. Real-time PCR Analysis of human intestinal microflora with 16S rRNA-gene-targeted Genus and species-specific primers. XI Congreso Nacional de Biotecnología y Bioingeniería.
2. Mai V *et al.* Associations between dietary habits and body mass index with gut microbiota composition and fecal water genotoxicity: an observational study in African American and Caucasian American volunteers. *Nutr. Journal* 2009, 8: 49 – 59.
3. Turnbaugh P *et al.* A core gut microbiome in obese and lean twins. *Nature* 2009, 457(7228): 480 – 484.
4. Vaarala O. Gut Microbiota and Type 1 Diabetes. *Rev. Diab. Stud.* 2013, 9(4): 251 - 259.

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No aplicable

## 16. Bibliografía

1. Matsuki. Real-time PCR Analysis of human intestinal microflora with 16S rRNA-gene-targeted Genus and species-specific primers. XI Congreso Nacional de Biotecnología y Bioingeniería.
2. Mai V *et al.* Associations between dietary habits and body mass index with gut microbiota composition and fecal water genotoxicity: an observational study in African American and Caucasian American volunteers. *Nutr. Journal* 2009, 8: 49 – 59.
3. Turnbaugh P *et al.* A core gut microbiome in obese and lean twins. *Nature* 2009, 457(7228): 480 – 484.
4. Vaarala O. Gut Microbiota and Type 1 Diabetes. *Rev. Diab. Stud.* 2013, 9(4): 251 -259.

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