RIDASOFT®WIN.NET

REF Z9996

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PC-Software for the evaluation of data

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Please contact for questions and further information:

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RIDA®, RIDASCREEN® and RIDASOFT®
are registered trademarks of R-Biopharm AG.
Manufacturer: R-Biopharm AG, Darmstadt, Germany

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Content

	page
1. Intended use	22
2. Minimum System requirements	22
3. Installation	22
4. Use of the software	23
4.1. Starting the program	23
4.1.1. Allergy	23
4.1.2. Clinical	23
4.1.3. Most-recently accessed method file	23
4.2. Changing a method	23
4.3. Preparation of the reader for your assay measurements	25
4.4. Plate layout	26
4.4.1. Specific IgE	26
4.4.2. Specific IgG	26
4.4.3. All other methods	26
4.5. Measuring	27
4.6. Evaluation	27
4.7. Printing and exporting data	27
5. Explanation of the command line	27
5.1. File menu	27
5.2. Edit menu	29
5.3. Results menu	30
5.4. View menu	32
5.5. Windows menu	32
5.6. Help menu	32
6. Troubleshooting	33

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RIDASOFT®WIN.NET – User Manual

1. Intended use

For *in vitro* diagnostic use. The RIDASOFT®WIN.NET is a software for calculation and interpretation of concentrations based on optical density from photometers. These are carried out based on assay specific parameters of the validated qualitative and quantitative clinical RIDASCREEN® enzyme linked immunosorbent assay (ELISA) as well as ELISA distribution products.

For professional use only.

2. Minimum System requirements

- OS Win XP SP3 or higher
- .NET Framework from Version 4.0
- RIDASOFT®WIN.NET Version 1.107.1.0247
- For the installation of the software: administration rights

3. Installation

- Insert the installation disk into the disk drive
- Open *setupDotNet.exe*
- The setup program will install RIDASOFT®WIN.NET on your computer
- Start RIDASOFT®WIN.NET
- Select your preferred language
- You will be asked if you would like *step-by-step instructions*
- If you confirm you will be introduced with the program functionality step by step
- If you do not wish to be asked again each time the program starts please click the check button
- All help topics including the test implementation can be accessed by using the *help index* command in the help menu

4. Use of the software

4.1. Starting the program

4.1.1. Allergy

Click on the Windows start button and select the menu item *Program/ RIDASOFT WIN.NET/RIDASOFT WIN.NET*. In the **Select item** window that opens, select the *Allergy* item and confirm your choice by clicking on the *OK* button.

All the methods stored for the selected item are listed in the window that follows. The following standard methods are already stored in the software for the *Allergy* item:

spezige.met: for specific IgE
spezigg.met: for specific IgG
totalige.met : for total IgE

These methods cannot be changed.

Select the method you require and confirm your choice with *OK*. See point 4.2. for changing the settings of the methods.

Please note that after the settings have been changed, the methods must be stored under a different name.

4.1.2. Clinical

Click on the Windows start button and select the menu item *Program/ RIDASOFT WIN.NET/RIDASOFT WIN.NET*. In the **Select item** window that opens, select the *Clinical* item and confirm your choice by clicking on the *OK* button.

All the subfolders stored for the selected item are listed in the window that follows.

Select the subfolder and then the method you require and confirm your choice with *OK*. See point 4.2. for changing the settings of the methods.

4.1.3. Most-recently accessed method file

Click on the Windows start button and select the menu item *Program/ RIDASOFT WIN.NET/RIDASOFT WIN.NET*. In the **Select item** window that opens, select the *most-recently accessed method file* item and confirm your choice by clicking on the *OK* button. See point 4.2. for changing the settings of the methods.

4.2. Changing a method

Change a method by clicking on this item or on the menu item *File/Open/Method*. Select the method you want to change from the list and confirm your choice with *OK*.

In order to change this method, click on this icon or on the menu item *Edit/Method*. You can choose the specifications for the method in the window that then opens:

Standards

- Amounts:** Enter the desired number of standards. You can now switch between the standards one by one with the cursor buttons and define their properties.
- Id (optional):** Use a maximum of 20 characters to name the standard. Initially standards are named Std1..Stdn. The standards' names are visible in the method view and in the print output.
- Replicates:** Enter the desired number of replicates for each standard. If you would like to use the same number of replicates for all standards, please press the auto button after putting in the desired number.
- Concentration:** Enter the desired concentration for each standard. Zero concentrations are only allowed for the first standard, which is the zero standard. Not available in *Clinical*.
- Serial dilution:** If you put in a dilution factor greater than zero, the concentration of each standard beginning with the highest is calculated. Enter the concentration for the highest standard in the concentration field. The concentration of the remaining standards is calculated using following formula:
Conc. Stdi = Conc. Stdi+1/dilution factor.

Properties

- Name:** The name states the name of the method to be changed. When saving edited methods you have to use new filenames. The supplied method files are write protected and not changeable.
- Dimension:** Input the unit in which you want the results displayed, usually IU/ml or kU/L for *Allergy*. Not available for *Clinical*.
- Evaluation:** The evaluation indicates the calculation for the method. Depending on the selected item there are different calculation methods available. In *Allergy* the standard curve is compiled by joining the individual standards in linear fashion (point-to-point). Some of the settings may not be changed.

Samples

- Amounts:** Reflects the remaining maximum amount of available samples. Switch the samples one by one with the cursor buttons and define their properties.
- Id (optional):** Use a maximum of 20 characters to name the samples. Initially samples are named P1 ... Pn. The sample names are visible in the method view and in the print output.
- Replicates:** Enter the desired number of replicates for each sample. If you would like to use the same number of replicates for all samples, please press the auto button after inputting the desired number.
- Factor: (optional)** Enter the desired dilution factor for each sample. If you would like to use the same dilution factor for each sample, please press the auto button after inputting the desired factor. If you would like to change only the samples' dilution factor for the current method it is more suitable to use the Edit Samples command.

Controls

- Amount:** Enter the desired number of controls. You can now switch between the samples one by one with the cursor buttons and define their properties.
- Id (optional):** Use a maximum of 20 characters to name the controls. 2 controls are normally used and are named QC1 and QC2. The control names are visible in the method view and in the print output.
- Replicates:** Enter the desired number of replicates for each control. If you would like to use the same number of replicates for all controls, please press the auto button after inputting the desired number.

4.3. Preparation of the reader for your assay measurements

Before performing your first assay measurement, you have to adjust RIDASOFT® WIN.NET to your photometer.

To determine the photometer and the serial port use the *Reader Setup* command. Select your reader in the dialog box. The following dialog box determines the serial port settings:

The dialog controls are initialized with the selected device specifications. These specifications refer to the reader manual. Please compare the settings. It is possible that your reader uses different specifications. In this case you have to change the settings:

Port:

Select the serial port connected to your reader.

Baud Rate:

Select the speed the reader transfers its data.

Data Bits:

Select the number of data bits.

Parity:

Select the kind of parity bit.

Stop Bits:

Select the number of stop bits.

Flow (optional):

Select the kind of data flow control. Normally readers use a hardware handshake for communication with the computer. Please change this parameter only if you know that your reader uses a different kind of flow control than the one already checked.

If the device is remotely controlled by RIDASOFT®WIN.NET, please select the filter wavelengths in the following step.

These settings are saved on disk and are used every time you start a new measurement.

Important:

Be sure your reader works in computer mode. Readers obtained from the manufacturer are in some cases set to printer mode, thus the data is transferred to a printer connected to the reader. In this case you have to switch to computer mode. Please refer to the reader manual. Be sure the serial cable is RS-232C standard. In most cases you must obtain the serial cable from the reader manufacturer.

4.4. Plate layout

4.4.1. Specific IgE

You have selected the **spezige.met** method at program start and are now in the main menu. If

another method is loaded, you can select the method for specific IgE by clicking on this icon or by clicking on the menu item *File/Open/Method* and confirming with *OK*.

Occupy the plate by clicking on the menu item *Edit/Layout*. The Layout window opens.

Input the first patient by clicking on the *New Patient* button. You can now input the number and the name (optional).

Allergens are selected either by editing the allergen code directly in the *Allergen* field or by selecting an allergen from the code list or the name list. Each time an allergen is edited, this needs to be confirmed with the *Enter* key. Every allergen that is confirmed is immediately inserted into the occupancy matrix.

When all the allergens of the first patient have been entered and you wish to input further patients, click on the *New Patient* button and proceed as with the first patient. After all patients have been entered, terminate the input process by clicking on the *Complete* button.

The counter shows you how many allergens can still be laid on this plate whilst you are entering the allergens.

4.4.2. Specific IgG

In order to measure specific IgG, the method **spezigg.met** must have been selected when the software was started or must be selected as detailed under 4.4.1. Further occupancy is made in the same way as for specific IgE as described under 4.4.1.

4.4.3. All other methods

The layout of the plate with standards and controls is already stipulated. The layout shifts if you change the characteristics of the method (number of standards, number of repeats, etc., see point 4.2.). If you want to alter the layout made automatically, mark the *BLK* entry in the small right-hand window. Then position the mouse over all entries (standards and controls) one after the other and press the left-hand mouse button each time.

In order to now occupy the plate with the standards and controls according to the stipulations of the selected method, mark the desired entry in the list and then position the mouse over the desired position in the layout window and then press the left-hand mouse button.

Select *P* in the list to occupy your samples and click on as many fields in the layout window as the number of samples you want to specify. In order to allocate an ID to the samples, select the menu item *Edit/Samples*. A window then opens in which you can allocate an ID to each sample, which then appears in the printout of results. You can delete an existing entry with *BLK*. You can occupy the plate in one step with *Auto*. Mark any assignment of fields and the standards/sample are occupied according to their number and replicates in the desired allocation in chronological sequence on the plate.

4.5. Measuring

Before you can measure the occupied and processed plate, you must ensure that the photometer is connected and switched on and that the correct reader is set in the software (see point 4.3.).

Place the plate in your photometer and start the measurement either by clicking on this icon or by clicking on the menu item *File/New/Assay*.

4.6. Evaluation

Click on this icon or click on the menu item *Results/Default* to start the evaluation. You can view the measured values, the plate occupancy, the method used, the standard curve (quantitative evaluations) and the sample results on the lower edge of the window by clicking on the respective tab.

Save the test by clicking on this icon or click on the menu item *File/Save/Assay* and enter a name under which you want to save the test.

Stored tests can be reopened by clicking on this icon or by clicking on the menu item *File/Open/Assay* and then on the name under which the test was stored. As only the raw data were stored, the evaluation must be made again to enable the results to be visualised.

4.7. Printing and exporting data

Click on this icon or click on the menu item *File/Print* to print the results. Select *Results* if you want to print out all results, including the standard curve of the plate layout etc., in list form. Click on the desired printouts if you do not want to print out all data. If you want to print out the results oriented by patient as individual findings, click on *report*. This print function is only available when measuring specific IgE and specific IgG. You can select all measured patients or individual patients.

In order to make the results available to other database programs, click on the menu item *File/Save Sample Results*. The data is then stored in the ALLERGDB.TXT file for *Allergy* or RIDADB.TXT file for all other items and can be read in by other database programs.

5. Explanation of the command line

5.1. File menu

New Method command (same as *Edit Method* command)

Use this command to create a new method based on the current method. You can open an existing method with the *Open Method* command.

Shortcut
Keys: CTRL+N

New Assay command

Use this command to start a new measurement.

You need to start the reader after activating this command for manual reader operation.

The data can be seen in the plate values view. You can open an existing assay with the Open Assay command.

Shortcut
Keys: CTRL+M

Open Method command

Use this command to open an existing method. You can create a new method with the New Method command.

Shortcut
Keys: CTRL+O

Open Assay command

Use this command to open an existing assay. You can start a new assay measurement with the New Assay command.

Shortcut
Keys: CTRL+A.

Save Method command

Use this command to save and name the active method along with the current plate layout. RIDASOFT® WIN.NET displays the Save As dialog box so you can name your method. The supplied method files are write protected.

Shortcut
Keys: CTRL+S

Save Assay command

Use this command to save and name the active assay along with the current plate layout.

RIDASOFT® WIN.NET displays the Save As dialog box so you can name your assay. Existing assay files are write protected.

Shortcut
Keys: CTRL+R

Save Sample Results command

Use this command to save the current samples' data. The file RIDADB.txt is automatically updated. You can import this file to a database application to create your own database.

Classify Samples command

Use this command to classify selected samples (for example suppliers, customers.).

Add (a category):

Press this button after filling the edit box and the list will be updated.

Remove (a category):

Press this button after selecting a list item and the item will be removed from the list.

Classify:

Press this button after selecting the desired samples in the layout view and after selecting the desired list item. The samples will be classified with the selected list item. Successful classification will be confirmed with an 'ok' message.

This command is unavailable if the layout view is not active.

The classification of samples is displayed in the method view of the method layout. Use the Save Sample Result command to add the data to the samples list.

Print command

Use this command to print a selection of method specific data. This command presents the Selection dialog box, where you may select the desired information to be printed.


Shortcut
Keys: CTRL+P

5.2. Edit menu

Cut command

Use this command to remove the currently selected data from the document and copy it to the clipboard. This command is unavailable if there is no data currently selected.

Cutting data to the clipboard replaces the contents previously stored there.

Shortcut
Keys: CTRL+X.

Copy command

Use this command to copy selected plate values onto the clipboard. This command is unavailable if there is no plate value currently selected.


Copying data to the clipboard replaces the contents previously stored there.

Shortcut
Keys: CTRL+C

Paste command


Use this command to insert a copy of the clipboard contents at the insertion point. This command is unavailable if the clipboard is empty.

Shortcut
Keys: CTRL+V

Method command (same as New Method command)

Use this command to edit the current method. You can open an existing method with the Open Method command.

Shortcut
Keys: CTRL+N

Samples command

Use this command to change the sample dilution factor and/or the samples id for the current method.


Factor:

Put in the desired dilution factor for each sample. If you like to use the same dilution factor for each sample, please press the auto button after putting in

the desired factor. Use the cursor buttons to switch between samples.

Id (optional):

Use a maximum of twenty characters to name the samples. Initially samples are named P1 ... Pn. The sample names are visible in the method view and in the print output.

5.3. Results menu


Default command

Use this command to start the method specific calculation. The view will be switched to method layout, and you will have access to the results in the graph view (quantitative evaluations) and in the samples view.


4P-Method command

Use this command to calculate the 4P-curve from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.


Linear Regression command

Use this command to calculate the regression line from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

Logit/Log command


Use this command to calculate the Logit/log from the standard OD values. The view will be switched to method layout, you have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.

Shortcut

Cubic Spline command

Use this command to calculate the cubic spline from the standard OD values. The view will be switched to method layout, and you will have access to the results in the graph view and in the samples view. This command is unavailable if the current method uses a qualitative or point-to-point calculation.


Shortcut

5.4. View menu

Toolbar command

Use this command to display and hide the Toolbar, which includes buttons for some of the most common commands in RIDASOFT®WIN.NET, such as File Open. A checkmark appears next to the menu item when the Toolbar is displayed.

Status Bar command


Use this command to display and hide the Status Bar, which describes the action to be executed by the selected menu item or depressed toolbar button, and keyboard latch state. A checkmark appears next to the menu item when the Status Bar is displayed.

5.5. Windows menu

Plate Layout command


Use this command to switch to the Plate Layout mode. This command is unavailable if the mode is already active.

Shortcut
Keys: CTRL+L

Method Layout command

Use this command to switch to the Method Layout mode. In this mode you can check the template of the current method in tabbed views. When performing a calculation you can check the results in the updated view. This command is unavailable if the mode is already active.

Shortcut
Keys: CTRL+H

5.6. Help menu

Index command

Use this command to display the opening screen of *Help*. From the opening screen, you can jump to step-by-step instructions for using RIDASOFT® WIN.NET and various types of reference information.


Once you open Help, you can click the Contents button whenever you want to return to the opening screen.

Using Help command

Use this command for instructions about using Help.

Information

Use this command to display the copyright notice, version number and the serial number of your copy of RIDASOFT®WIN.NET.

6. Troubleshooting

Error Description	Possible Source	Solution
Error message during measurement “Reader not found or reader not started!”	Reader not correctly connected or not turned on	Connect and/or turn on the reader correctly (see 4.3.)
Error message during measurement „Data transfer failed.“	Reader is not working properly	Check functionality of the reader
Error message in reader properties: „Connection failed.“	Reader not correctly connected or not turned on	Connect and/or turn on the reader correctly (see 4.3.)
Error during reading/initializing/writing of files	No reading/writing permission to the program directory	Grant reading/writing permissions
Connection issues with Biotek ELX800	Wrong baudrate set and/or data flow control DTR/DSR not active	Set baudrate DTR/DSR 2400 and/or activate data flow control DTR/DSR (see 4.3)
Connection issues with Multiskan EX	Data flow control active	Disable all data flow control options (see 4.3)
Connection issues with Metertech Accureader	Data is not provided to the PC	Reader: SETUP->COMPUTER ->RAW-DATA Measurement: MODE->END POINT(S) ->CONTINUOUS

Further product information and applications, please contact your local distributor or R-Biopharm at this address: info@r-biopharm.de.

Version overview

Version number	Chapter and title
2015-02-20	Release version
2021-02-03	Adaptation to new master template and new intended use. All references to food & feed analytics have been removed. Added Software Version requirement.

Explanation of symbols

- General symbols:

Follow the instructions for use

Batch number

Expiry date (YYYY-MM)

Storage temperature

Article number

Number of test determinations

Manufacturing date (YYYY-MM)

Manufacturer + address