Instructions for use/Testanleitung

Pankrin™ ELISA

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English:

Intended Use

With the Pankrin™ ELISA from BIOSERV Diagnostics an acute pancreatitis or an acute attack of chronic pancreatitis is diagnosed in patients with abdominal pain.

Clinical Relevance

An acute pancreatitis, the acute inflammation of the pancreas, is mostly caused by the dystopic protease activation in the azinus cells, the obstruction of the luminal secretion and the formation of intracellular vacuoles (Schneider 1999). The yearly incidence rate is 50 to 100 cases per 100.000 persons. The clinical picture is determined by pronounced abdominal pain and general symptoms like nausea, vomiting, and circulatory collapse (Selberg *et al.* 1995). The mild form of an acute pancreatitis, the endematous pancreatitis (ca. 80% of all cases of an acute pancreatitis) progresses without complications and under adequate therapy patients normally recover within 72 hours. The severe form of an acute pancreatitis, the hemorrhagic narcotising pancreatitis is characterized by necrosis and by a protracted course with frequent failure of organs or organ systems (Asanuma *et al.* 1999, Singer *et al.* 1988). Morphologically a thickening and irregular delimitation of the pancreas can be observed as well as an accumulation of liquid in the peritoneal space. The main cause of an acute pancreatitis are an increased alcohol consumption (Blank *et al.* 1999, Lankisch 2000) and trapped gallstones which possibly get into the pancreatic duct via a reflux (Niederau *et al.* 1997, UHL *et al.* 1999, Sharma *et al.* 1999). Other pathogenetic factors are more rare. In a considerable percentage of the patients suffering from acute pancreatitis no risk factor can be detected. Cholelithiasis (gallstone disease) peaks between 40 and 60 years, the age maximum for alcohol abusus lies between 20 and 40 years.

Fields of Application

In the functional diagnostics of the exocrine pancreas the BIOSERV Pankrin™-ELISA is a sensitive and clinically easily practicable tool ideally complementing the sonographically established initial diagnosis.

Principles of the Assay Method

The proteins detected by the BIOSERV Pankrin™-ELISA are absolutely pancreas-specific. In addition to pancreatic elastase also other pancreatic enzymes are detected to increase the diagnostic sensitivity and specificity of the test. In marked contrast to clinical chemistry tests measuring only the enzymatic activity with the BIOSERV Pankrin™-ELISA the pancreatic proteins are determined immunologically, their immunological properties depending on the molecular structure. Therefore not only an absolute pancreas-specificity can be guaranteed but it is also possible to measure already inactivated proteins in the serum. The pancreatic proteins getting into the serum during an acute attack can thus be detected for a much longer period of time than with conventional clinical chemistry test systems.

The BIOSERV Pankrin™-ELISA (**E**nzyme **L**inked **I**mmuno**S**orbent **A**ssay) is a solid phase enzyme immunoassay based on the double sandwich technique. The wells are coated with polyclonal antibodies specifically recognizing pancreatic proteins. Pancreatic proteins from the serum samples or standards are immobilized by these antibodies. Then a secondary polyclonal antibody solution also specifically reacting with pancreatic proteins is added and binds to the pancreatic proteins immobilized on the first antibody layer in the well. In the following incubation step the biotin-residue of the bound biotin-antibody-complex reacts with an added peroxidase-streptavidin-complex which then enzymatically oxidizes the added substrate TMB (3,3',5,5'-tetramethylbenzidine). This enzymatic activity is stopped after a defined time period by adding 0.25 mol/l H₂SO₄ and the amount of oxidized TMB is determined photometrically at 450 nm. The use of a reference measurement with a wavelength ≥550 nm is recommended.

Reagents

(sufficient for 96 determinations)

- | | |
|---------------------------------------------------------------------------------------------------------|-----------|
| 1. Pankrin™ reference standard set (lyophilised) 2 x - per vial | 0.5 ml |
| – Standard 1 (100 Pankrin™ units/ml serum – colourless screw cap) | |
| – Standard 2 (200 Pankrin™ units/ml serum – white screw cap) | |
| – Standard 3 (400 Pankrin™ units/ml serum – yellow screw cap) | |
| – Standard 4 (800 Pankrin™ units/ml serum – blue screw cap) | |
| 2. Positive control, lyophilised (corresponding to 200 Pankrin™ units/ml serum ±20%, – green screw cap) | 0.5 ml |
| 3. Biotinylated secondary polyclonal anti-pancreas enzyme antibodies (red screw cap) | 0.12 ml |
| 4. Biotin conjugate dilution buffer | 8 ml |
| 5. Streptavidin horseradish peroxidase conjugate (ready for use) | 8 ml |
| 6. Washing solution (10x concentrated) | 2 x 50 ml |
| 7. Substrate solution (solution of TMB, ready for use) | 13 ml |
| 8. Stop solution (0.25 mol/l H ₂ SO ₄) | 13 ml |
| 9. Microtiter strips coated with polyclonal anti-pancreas enzyme antibodies | 96 wells |
| 10. Holder for single strips | 1 x |

Materials Required but not Included

1. Microplate reader with 450 nm filter, optionally with a reference filter ≥550 nm.
2. Microliter pipettes with disposable tips: 5 µl, 50 µl, 100 µl and 1000 µl.
3. Tubes for the dilution of the samples.
4. Distilled or deionised water.
5. Absorbent paper.
6. Please use only calibrated pipettes and instruments.

Warnings and Precautions

1. This kit is intended for *in vitro* use only.
2. Avoid contact with the stop solution, it may cause skin irritations and burns.
3. Do not pipette reagents by mouth.
4. Please regard all samples as potentially infectious and handle them with utmost care.
5. Handling and disposal should be in accordance with the procedures defined by an appropriate national biohazard safety guideline or regulation where this exists.

Instructions for Reagent Preparation

1. The components of this kit are intended for use as an integral unit and should not be interchanged with the components of other kits.
2. All reagents and specimens must be brought to room temperature before use.
3. All reagents have to be mixed without foaming.
4. Once the test procedure has been started, all steps should be continued without interruption.
5. Pipette all reagents and samples onto the bottom of the wells. Mixing or shaking after pipetting is not required.
6. Use new disposable tips for each specimen.
7. Before starting the assay, all reagents to be used should be prepared and ready for immediate use, all needed strips should be secured in the holder etc. This will ensure equal time periods for each pipetting step without interruption.
8. For optimal results it is important to wash the wells thoroughly after incubation and to remove even the last water drops by hitting the plate on absorbent paper or cloth.
9. Since the kinetics of the enzymatic reaction depends on the surrounding temperature different extinctions correlating with the respective room temperature may be observed. The optimum laboratory room temperature is 21,5 °C (71 °F).
10. It is recommended to effect all tests in double determination in order to minimize the consequences of pipetting or handling errors.

Storage Instructions and Shelf Life Information

1. Store the reagents at 2 – 8 °C (36 °F – 46 °F).
2. The reagents remain stable until the expiration date of the kit.
3. The diluted washing solution is stable for 4 weeks at refrigerator temperatures (2 °C – 8 °C / 36 °F – 46 °F).
4. Put caps back on the vials immediately after use.
5. Store the microtiter strips in a dry bag with desiccants. The remaining strips must be stored in the tightly resealed bag together with the desiccants. Under these storage conditions, they are stable at least for 4 weeks after opening of the sealed bag.

Sample Material

Serum

Specimen Collection and Preparation

Collect blood by venipuncture, allow to clot, and separate serum by centrifugation at room temperature; avoid haemolysis. Avoid repeated freezing and thawing. Store tubes closed as they may be a danger of contamination or alteration of concentration.

1. Handle all samples with utmost care since they may be infectious.
2. There are no known interferences with extrinsic factors or other substances.
3. Samples may be stored at different temperatures for the following time-spans:
 - Environmental temperature up to 30 °C (86 °F): up to three days
 - Refrigerator temperature (2 – 8 °C / 36 °F – 46 °F): up to one week
 - Household freezer temperature (-10 °C – -20 °C / 14 °F – -4 °F): up to one year

ATTENTION! There are no test methods available which may guarantee that Hepatitis B virus, Human Immunodeficiency Virus (HIV/HTLV-III/LAV), or other infectious agents are absent from the reagents in this kit. Therefore, all human blood products, including patient samples, should be considered potentially infectious.

Assay Procedure

1. Preparation of the washing solution(10x): Dilute the concentrated washing solution (50ml) with 450 ml distilled or deionised water. **Attention:** Do not use unpurified tap water!
2. Preparation of the lyophilised standards: Reconstitute each lyophilised standard with 500 µl aqua bidest. Do not use buffer! The reconstituted standards can be stored at 2 °C – 8 °C (36 °F – 46 °F) for four weeks.
3. Warm all reagents to room temperature and mix thoroughly before use.
4. Fix the required number of coated wells or strips in the strip holder. Fill in a paper form accordingly. Use up the microtiter strips within four weeks.
5. Dilute the serum samples 1:101 (1 + 100) with washing solution (e.g. 25 µl serum + 2.5 ml washing solution). The reconstituted standards are ready for use.
6. Pipette 50 µl each of the zero standard (blank, washing solution), the standards, positive control and the diluted serum samples with new disposable tips into the respective wells.
7. Incubate for 60 minutes at room temperature (18 – 25 °C).
8. Briskly shake out the contents of the wells and then rinse the wells 3 times with 200 µl diluted (see step 1) washing solution. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
9. Dilute the biotinylated anti-pancreatic protein antibody solution 1:201 with biotin conjugate buffer: (1 part biotinylated antibody + 200 parts conjugate buffer). Examples:
 - 5 µl anti-pancreatic protein antibody solution + 1 ml biotin conjugate buffer
 - 10 µl anti-pancreatic protein antibody solution + 2 ml biotin conjugate buffer
 - 15 µl anti-pancreatic protein antibody solution + 3 ml biotin conjugate buffer
 - 25 µl anti-pancreatic protein antibody solution + 5 ml biotin conjugate bufferDispense 50 µl into each well.
10. Incubate for 30 minutes at room temperature (18 – 25 °C).
11. Briskly shake out the contents of the wells and then rinse the wells 3 times with 200 µl diluted (see step 1) washing solution.
12. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
13. Dispense 50 µl of streptavidin peroxidase conjugate (ready for use) into each well.
14. Incubate for 30 minutes at room temperature (18 – 25 °C).
15. Briskly shake out the contents of the wells and then rinse the wells 3 times with 200 µl diluted (see step 1) washing solution.
16. Knock the residual water out of the wells by hitting them (in the holder) on absorbent paper or cloth.
17. Dispense 100 µl of substrate solution into each well.
18. Incubate for 20 min at room temperature (18 – 25 °C) Attention: please take the time from start of pipetting the substrate.
19. Stop the enzymatic reaction by adding 100 µl stop solution to each well, in the same sequence and time interval as dispensing the substrate.
20. Measure the extinction of the samples at 450 nm. It is recommended to carry out the measurement of the extinction within 10 minutes after stopping the reaction.
21. The value of the positive control should be between 160 to 240 Pankrin units/ml serum. If the positive control should be out of this range, please repeat the test. In this case please check once again all incubation steps and times.
 - 21 a) If the positive control included to this kit repeats giving results higher than 240 Pankrin units/ml serum, the kit must be discarded.
 - 21b) If the positive control included to this kit repeats giving results lower than 160 Pankrin units /ml serum, the kit must be discarded.

As a general rule the enzymatic reaction is linearly proportional to time and temperature. This makes interpolation possible for fixed physico-chemical conditions.

If in a test run the absorbance of the 800 U/ml-standard is lower than 1.0 the incubation time of the final enzymatic reaction may be extended.

If, on the other hand, the absorbance of the 800 U/ml-standard is above the upper performance limit of the microplate spectrophotometer used the enzymatic reaction time may be reduced.

Since calibrators are assayed in each run, absorbance fluctuations do not affect the absolute results.

Pipetting Scheme for the Pankrin™ Serum ELISA from BIOSERV Diagnostics

	1	2	3	4	5	6	7	8	9	10	11	12
A	BLANK	BLANK	P	3	P	11	P	19	P	27	P	35
B	S	1	P	4	P	12	P	20	P	28	P	36
C	S	2	P	5	P	13	P	21	P	29	P	37
D	S	3	P	6	P	14	P	22	P	30	P	38
E	S	4	P	7	P	15	P	23	P	31	P	39
F	P	C	P	8	P	16	P	24	P	32	P	40
G	P	1	P	9	P	17	P	25	P	33	P	41
H	P	2	P	10	P	18	P	26	P	34	P	42

In this pipetting scheme the recommended positions for the blank, standards (S1 – S4), positive control (PC) and for the patient samples (P1 – P42) are shown as double determinations.

Calculation of the Results

1. Calculate the average absorbance values for each set of reference standards, controls and patient samples
2. The optical density of each standard value is plotted as y value (y-axis), the corresponding Pankrin™ value is drawn in as the x-value (x-axis). The resulting calibration curve is used to determine the values of the patient samples. The OD values of the serum samples are correlated with the corresponding Pankrin™ concentration values by interpolation.
3. Using the mean absorbance value for each sample determine the corresponding concentration of Pankrin™ in U/ml from the standard curve.

Limitations of the Assay

- At temperatures higher than 30 °C (86 °F) the samples should be transported cooled or refrigerated. The time to stop the (enzymatic colour) reaction may have to be adjusted (shortened).
- Severely haemolytic or lipaemic sera or sera from patients with liver diseases should not be used. Results may be adversely affected by certain pathologic conditions, such as poly- and monoclonal gammopathies, autoimmune diseases or by an altered immune status.

Expected Values

- Normal: < 160 Pankrin™–Units/ml serum
- Near Cut-Off: 160 – 190 Pankrin™–Units/ml serum
- Pathological: > 190 Pankrin™–Units /ml serum

These ranges were determined in a normal collective of 558 blood donors and 177 clinically clarified patients with a specificity of 96% at a sensitivity of 98%.
The cut-off of the Pankrin ELISA is at 175 Pankrin™–Units/ml serum. In case of a value in the range near the cut-off (160 to 190 Pankrin™–Units/ml serum) we recommend a follow-up determination using a new serum sample taken within 24 hours. Other additional diagnostic examinations should follow, then.

Assay Performance Characteristics

1. Diagnostic Specificity: 96%

Samples from 558 healthy test persons were investigated. The test persons were either healthy blood donors or patients for whom pancreatic ailments had been excluded by other diagnostic methods.

2. Diagnostic Sensitivity for Chronic Pancreatitis: 98%

Samples from 177 patients suffering from acute pancreatitis were investigated.

3. Intraassay variation coefficient: 7.48% (5.07 – 9.10%)

For the determination of the intraassay variation coefficient 6 kits from 6 different batches (produced on different days) were used.

One patient sample (optical density about 1.0) was applied 96 times per testing procedure.

4. Interassay variation coefficient: 6.51% (5.40 – 8.70%)

In order to determine the coefficient of interassay variation one strip each of 12 kits stemming from 6 different batches (produced on different days) were used.

One patient sample (optical density about 1.0) was applied 72 times per testing procedure.

Bibliography

Especially for Pankrin®-Serum-ELISA from BIOSERV Diagnostics

Lankisch PG, Weber-Dany B, Doobe C, Finger T, Maisonneuve P, Lowenfels AB, Keim V (2006): Pankrin: A new parameter for the diagnosis of acute pancreatitis in cases of late clinical presentation. *Pancreas* Volume 32/3, pages 330-331.

Keim, V, Teich, N, Bodeker H, Mossner, J (2003): Evaluation of Pankrin®, a new serum test for diagnosis of acute pancreatitis. *Clinica Chimica Acta* 332 , pages 45 – 50.

Generally for acute Pancreatitis

Asanuma Y, Fukuya T, Tanaka J-I, Sato T, Shibata S, Koyama K (1999): The Application of Immobilized Polymyxin B Fiber in the Treatment of Septic Shock Associated with Severe Acute Pancreatitis: Report of Two Cases. *Jpn. J. Surg* 29, pages 1177-1182.

Bank S, Indaram A (1999): Causes of acute and recurrent pancreatitis. *Gastroenterology Clinics of North America* 28, pages 571-589.

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Sharma VK, Howden CW (1999): Metaanalysis of randomised controlled trials of endoscopic retrograde cholangiography and endoscopic sphincterotomy for the treatment of acute biliary pancreatitis. *Am. J. Gastroenterol.* 94, pages 3211-3214.

Singer MV, Gyr KE, Sarle H (1988): Revised classification of pancreatitis. Report of the Second International Symposium on the Classification of Pancreatitis in Marseille, France, March 1984. *Gastroenterology* 89 (1988); pages 683-690.

Uhl W, Müller, CA, Krähenbühl L, Schmid SW (1999) Schölzel, St., Büchler, W.: Acute gallstone pancreatitis. *Surgical Endoscopy* 13, pages 1070-1076.

Also available from BIOSERV Diagnostics:

Please contact us for more information: info@bioserv-diagnostics.com

Gastroenterology	Order Code	Determinations per Kit
Human Pancreatic Elastase ELISA (pancreatic insufficiency)	BS - 86 - 01	96
Accessories Kit: Stool Sample Preparation Kit with 45 Tubes	BS - 00 - 03	
Pankrin™ ELISA (acute pancreatitis)	BS - 86 - 02	96
Infertility		
Anti-Spermatozoa Antibody Latex Agglutination Test	BS - 10 - 10	50
Anti-Spermatozoa Antibody ELISA for Serum	BS - 10 - 20	96
Anti-Spermatozoa Antibody ELISA for Seminal Plasma	BS - 10 - 21	96
Anti-Spermatozoa Antibody Haemagglutination Test	BS - 10 - 30	40
Anti-Spermatozoa Antibody ELISA - Ig Classifying for Serum	BS - 10 - 50	96
Anti-Zona Pellucida Antibody Latex Agglutination Test	BS - 20 - 10	50
Anti-Zona Pellucida Antibody ELISA	BS - 20 - 20	96
Anti-Zona Pellucida Antibody Haemagglutination Test	BS - 20 - 30	40
Anti-Zona Pellucida Antibody ELISA - Ig Classifying	BS - 20 - 50	96
Anti-Ovary Antibody Latex Agglutination Test	BS - 40 - 10	50
Anti-Ovary Antibody ELISA	BS - 40 - 20	96
Anti-Ovary Antibody Haemagglutination Test	BS - 40 - 30	40
Anti-Ovary Antibody ELISA - Ig Classifying	BS - 40 - 50	96
Monitoring of Risk Pregnancies		
IGF-BP1 ELISA (PP12)	BS - 30 - 10	96
Glycodelin ELISA (PP14)	BS - 30 - 20	96

BIOSERV Diagnostics GmbH
Doberaner Str. 151
18057 Rostock
Tel.: 0381 3758 2090
Fax: 0381 3758 2099
E-mail: info@bioserv-diagnostics.com
Web: www.bioserv-diagnostics.com