

Need for IFX rapid testing? Therapeutic drug monitoring (TDM)

The lateral flow assay RIDA®QUICK IFX Monitoring is a rapid point-of-care assay, which allows the quantitative determination of infliximab trough levels (TLs) within 20 minutes. It facilitates the use of therapeutic drug monitoring in hospitals and infusion centers and the decision making in regards to therapy adaptation.

RIDA®QUICK IFX Monitoring, Art. No. GN3041



- Fast and reliable drug monitoring of infliximab (IFX, Remicade®) and its biosimilars Remsima® & Inflectra®
- Quantitative lateral flow test
- 20 minutes time to result
- High correlation to RIDASCREEN® IFX Monitoring due to identical monoclonal antibody (MA-IFX6B7)
- 25 test cassettes per package
- 2 controls available as accessory


RIDA®QUICK IFX Monitoring is read on the RIDA®QUICK SCAN II for precise and accurate measurement.

Scan to view video clip about
RIDA®QUICK IFX Monitoring




Products for therapeutic drug monitoring (TDM)

Infliximab lateral flow assay

Product	Description	Tests	Matrix	Art. No.
Lateral flow test and device				
 RIDA®QUICK IFX Monitoring	Immunochromatographic lateral flow assay for the quantitative determination of infliximab (IFX, Remicade®) and its biosimilars Remsima® and Inflectra®	25	Serum/plasma	GN3041
RIDA®QUICK IFX Monitoring Control Set	Positiv controls, accessory for GN3041			GP3041
RIDA®QUICK SCAN II	Lateral flow reader			ZRQS2-KD

ELISAs based on assays of KU Leuven, Belgium

Product	Description	Tests	Matrix	Art. No.
Enzyme immunoassay				
RIDASCREEN® IFX Monitoring	Enzyme immunoassay for the quantification of infliximab (IFX, Remicade®) and its biosimilars Remsima® und Inflectra®	96	Serum/plasma	G09041
RIDASCREEN® Anti-IFX Antibodies	Enzyme immunoassay for the quantification of anti-infliximab antibodies	96	Serum/plasma	G09042
RIDASCREEN® ADM Monitoring	Enzyme immunoassay for the quantification of adalimumab (ADM, Humira®)	96	Serum/plasma	G09043
 RIDASCREEN® Anti-ADM Antibodies	Enzyme immunoassay for the quantification of anti-adalimumab antibodies	96	Serum/plasma	G09044

For further details and information visit our website www.r-biopharm.com, contact your local distributor or Clinical Sales International.

R-Biopharm contacts:

Clinical Sales International:
 Phone: +49 (0) 61 51 - 81 02-0
 Fax: +49 (0) 61 51 - 81 02-688
 E-mail: clinical.sales@r-biopharm.de

Order Department:
 Phone: +49 (0) 61 51 - 81 02-0
 Fax: +49 (0) 61 51 - 81 02-20
 E-mail: orders@r-biopharm.de