



Your competent partner in gastroenterology

R-Biopharm AG

ELISA tests with diagnostic focus on:

Colon cancer prevention

RIDASCREEN® Haemoglobin
RIDASCREEN® Haemo-/Haptoglobin Complex

Inflammatory bowel disease

RIDASCREEN® Calprotectin
RIDASCREEN® α_1 -Antitrypsin
RIDASCREEN® sIgA

Pancreatic diagnostics

Pancreatic Elastase ELISA
Pankrin® ELISA

Therapeutic drug monitoring

RIDASCREEN® IFX Monitoring
RIDASCREEN® Anti-IFX Antibodies
RIDASCREEN® ADM Monitoring
RIDASCREEN® Anti-ADM Antibodies
RIDASCREEN® VDZ Monitoring
RIDASCREEN® GLM Monitoring
RIDA®QUICK IFX Monitoring
RIDA®QUICK ADM Monitoring

Haemoglobin and Haemo-/Haptoglobin Complex
Markers for colon cancer screening. A combined testing of both parameters increases the probability to detect adenomas.

Calprotectin

Indicator of gastrointestinal inflammation, used for differentiation of IBD and IBS and for treatment monitoring.

α_1 -Antitrypsin and sIgA

Indicators for the integrity of the gut wall, relevant for monitoring and diagnosing a Leaky Gut Syndrome.

Therapeutic Drug Monitoring (TDM)

TDM of Infliximab (IFX) and its biosimilars as well as Adalimumab (ADM), Vedolizumab (VDZ) and Golimumab (GLM) supports the therapy optimization and allows a treatment cost reduction at the same time.

Pancreatic elastase and pankrin

Markers for the detection of pancreas insufficiency and pancreatitis.

Accessories:

Stool collection tubes

Allow an easy, rapid and clean collection of stool samples.

R-Biopharm's Gastroenterology product portfolio is dedicated to the advanced diagnosis of gut-related diseases. The ELISAs are suited for single plate or medium-high-throughput analysis on automated ELISA systems.

ELISA

| Product | Description | Matrix | Tests | Art. No. |
|---|---|--------|-------|----------|
| Enzyme immunoassays | | | | |
| RIDASCREEN® Haemoglobin | Enzyme immunoassay for the quantitative detection of haemoglobin | Stool | 96 | G09030 |
| RIDASCREEN® Haemo-/Haptoglobin Complex | Enzyme immunoassay for the quantitative detection of haemoglobin/haptoglobin complex | Stool | 96 | G09031 |
| RIDASCREEN® α ₁ -Antitrypsin | Enzyme immunoassay for the quantitative detection of α ₁ -Antitrypsin | Stool | 96 | G09034 |
| RIDASCREEN® sIgA | Enzyme immunoassay for the quantitative detection of secretory IgA | Stool | 96 | G09035 |
| RIDASCREEN® Calprotectin | Enzyme immunoassay for the quantitative detection of calprotectin | Stool | 96 | G09036 |
| Pancreatic Elastase ELISA | Enzyme immunoassay for the quantitative detection of pancreatic elastase | Stool | 96 | G09038 |
| Pancreatic Elastase ELISA (SK15) | Enzyme immunoassay for the quantitative detection of pancreatic elastase; additional standard (SK15) | Stool | 96 | G09040 |
| Pankrin® ELISA | Enzyme immunoassay for the quantitative detection of pancreatic elastase and other pancreatic enzymes | Serum | 96 | G09039 |

Therapeutic Drug Monitoring (TDM)

| | | | | |
|---------------------------------------|---|--------------|----|--------|
| Enzyme immunoassays | | | | |
| RIDASCREEN® IFX Monitoring | Enzyme immunoassay for the quantitative detection of infliximab (IFX) and its biosimilars | Serum/plasma | 96 | G09041 |
| RIDASCREEN® Anti-IFX Antibodies | Enzyme immunoassay for the quantitative detection of antibodies to infliximab (IFX) and its biosimilars | Serum/plasma | 96 | G09042 |
| RIDASCREEN® ADM Monitoring | Enzyme immunoassay for the quantitative detection of adalimumab (ADM) | Serum/plasma | 96 | G09043 |
| RIDASCREEN® Anti-ADM Antibodies | Enzyme immunoassay for the quantitative detection of antibodies to adalimumab (ADM) | Serum/plasma | 96 | G09044 |
| RIDASCREEN® VDZ Monitoring | Enzyme immunoassay for the quantitative detection of vedolizumab (VDZ) | Serum/plasma | 96 | G09045 |
| New RIDASCREEN® GLM Monitoring | Enzyme immunoassay for the quantitative detection of golimumab (GLM) | Serum/plasma | 96 | G09047 |
| Rapid test | | | | |
| RIDA®QUICK IFX Monitoring | Immunochromatographic lateral flow assay for the quantitative determination of infliximab (IFX) and its biosimilars | Serum/plasma | 25 | GN3041 |
| New RIDA®QUICK ADM Monitoring | Immunochromatographic lateral flow assay for the quantitative determination of adalimumab (ADM) | Serum/plasma | 25 | GN3043 |

Accessories

| | | | | |
|--|---|--|----|--------------|
| Accessories for RIDASCREEN® | | | | |
| RIDA®TUBE | For collection and preparation of stool samples, <ul style="list-style-type: none"> • Unfilled; to use after internal validation | | 50 | GZ3013 |
| RIDA®TUBE Haemoglobin | For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with RIDASCREEN® Haemoglobin G09030 and RIDASCREEN® Haemo-/Haptoglobin Complex G09031 | | 50 | GZ3012 |
| RIDA®TUBE Calprotectin | For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with RIDASCREEN® Calprotectin G09036 | | 50 | GZ3016 |
| Stool Preparation Set | For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with Pancreatic Elastase ELISA G09038 and G09040 | | 45 | GZ3008 |
| RIDA®QUICK SCAN II - IVD SET | Lateral flow reader (CE-IVD,) and 2D barcode scanner <ul style="list-style-type: none"> • for read out of GN3041 and GN3043 | | | ZRQS2-KD-SET |
| RIDA®QUICK IFX Monitoring Control Set | Positive controls <ul style="list-style-type: none"> • accessory for GN3041 | | | GP3041 |
| New RIDA®QUICK ADM Monitoring Control Set | Positive controls <ul style="list-style-type: none"> • accessory for GN3043 | | | GP3043 |