



PRODUCT INFORMATION

Product Name:

PROCEEDx Influenza B and Adenovirus Positive Sample

Catalogue Number:

VP-14-03

Product Description:

This product contains inactivated influenza B and adenovirus in storage buffer.

SPECIFICATIONS

Strains:

Influenza B/Hong Kong/5/72. ICTV 00.046.0.04.001

Human adenovirus 2, strain Adenoid 6, ICTV 00.001.0.01.010.00.102.006

Storage:

Store this preparation at 2 °C to 8 °C.

Shelf Life:

72 weeks from manufacturing date.

Volume:

1 mL

USE

Instructions For Use:

Use PROCEEDx product as per sample handling instructions specified by device or instrument manufacturer.

Intended Use:

The product is to be used **for immunoassay** instrument (device) verification/validation, manufacturing QMS, educational and R&D purposes.

NOT FOR USE IN DIAGNOSTIC PROCEDURES. FOR RESEARCH USE ONLY.

ASSISTANCE

If you have any questions regarding the production, testing or use of this antigen, please send them together with any relevant data to Microbix Technical Services, at customer.service@microbix.com or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

HAZARDS


We are aware of no specific biological hazards associated with this product. The preparation has been inactivated and should contain no infectious material however it should be handled as if potentially infectious. Generally accepted good laboratory practices appropriate to biological reagents should be employed. The buffer solution contains 0.02% sodium azide (CAS No. 26628-22-8) as a preservative.




SPECIAL HANDLING

No special handling required for this product.

QUALITY CONTROL INFORMATION

| Test Kit | Target | Result |
|--|-----------------------------|---|
| Binax Now Influenza A & B Catalogue # 416-000 | Influenza A & B | Influenza B Positive Influenza A Negative |
| BD Directigen Flu A & B Catalogue # 256050 | Influenza A & B | Influenza B Positive Influenza A Negative |
| Quick Vue Flu A & B Catalogue # 20183 | Influenza A & B | Influenza B Positive Influenza A Negative) |
| Xpect Flu A & B Catalogue # : R24600 | Influenza A & B | Influenza B Positive Influenza A Negative |
| Binax NOW RSV Catalogue # 430-022 | Respiratory Syncytial Virus | RSV Negative |
| Combi-Strip (Rota and Adeno) Catalogue # C-1004 | Rotavirus & Adenovirus | Adenovirus Positive Rotavirus Negative |
| SAS Adeno Test Catalogue # 067020 | Adenovirus | Adenovirus Positive |
| SAS Rota Test Catalogue # 069020 | Rotavirus | Rotavirus Negative |



| | |
|---|------------------------|
| REF | Catalog Number |
| LOT | Lot Number |
|  | Temperature Limitation |
|  | Expiration Date |
|  | Biological Risk |

