



## PRODUCT INFORMATION

**Product Name:**

PROCEEDx Parainfluenza 3 Positive Sample

**Catalogue Number:**

VP-10-M1

**Product Description:**

This product contains inactivated Parainfluenza 3 virus in storage buffer.

## SPECIFICATIONS

**Strains:**

C243

**Storage:**

Store this preparation at 2 °C to 8 °C.

**Shelf Life:**

24 weeks from manufacturing date.

**Volume:**

1 mL

## USE

**Instructions For Use:**

Use PROCEEDx product as per sample handling instructions specified by device or instrument manufacturer.

**Intended Use:**

The product is to be used **for nucleic acid test** instrument (device) verification/validation, manufacturing QMS, educational and R&D purposes.

**NOT FOR USE IN DIAGNOSTIC PROCEDURES. FOR RESEARCH USE ONLY.**

### ASSISTANCE

If you have any questions regarding the production, testing or use of this antigen, please send them together with any relevant data to Microbix Technical Services, at [customer.service@microbix.com](mailto:customer.service@microbix.com) or fax 905-361-8911. Your complete satisfaction with the performance of this product is important to us.

### HAZARDS




We are aware of no specific biological hazards associated with this product. The preparation has been inactivated and should contain no infectious material however it should be handled as if potentially infectious. Generally accepted good laboratory practices appropriate to biological reagents should be employed. The buffer solution contains 0.02% sodium azide (CAS No. 26628-22-8) as a preservative.

### SPECIAL HANDLING

No special handling required for this product.

# QUALITY CONTROL INFORMATION

qPCR Test Kit	Analyte	Result
Nx TAG Respiratory Pathogen Panel Luminex Molecular Diagnostics	Parainfluenza 3	Detected
	Influenza A (H1, 2009 H1N1, H3)	Not Detected
	Influenza B	Not Detected
	Respiratory Syncytial Virus A	Not Detected
	Respiratory Syncytial Virus B	Not Detected
	Parainfluenza 1	Not Detected
	Parainfluenza 2	Not Detected
	Parainfluenza 4	Not Detected
	Coronavirur (229E, NL63, OC43, HKU)	Not Detected
	Metapneumovirus	Not Detected
	Rhinovirus/Enterovirus	Not Detected
	Adenovirus	Not Detected
	Bocavirus	Not Detected
	Chlamydia pneumoniae	Not Detected
Legionella pneumophilla	Not Detected	

REF	Catalog Number
LOT	Lot Number
	Temperature Limitation
	Expiration Date
	Biological Risk

