

Infosheet

RIDA®SEEK Software product and platform compatibility

		Real-time PCR instruments				
		Roche LightCycler® 480 II	Bio-Rad CFX96™	RIDA®CYCLER	ABI 7500 Fast	RotorGeneQ Cycler
Gastrointestinal infections						
Viruses						
PG1315	RIDA®GENE Viral Stool Panel I	●	—	—	—	—
PG1325	RIDA®GENE Viral Stool Panel II	●	●	●	◆	◆
PG1335	RIDA®GENE Viral Stool Panel III	●	●	●	◆	◆
PG1405	RIDA®GENE Norovirus	●	●	●	◆	◆
PG1415	RIDA®GENE Norovirus I & II	●	●	●	◆	◆
PG4705	RIDA®GENE Enterovirus	◆	◆	◆	◆	◆
Bacteria						
PG2205	RIDA®GENE EHEC/EPEC	●	●	◆ July 22	◆	◆
PG2215	RIDA®GENE EAEC	●	●	◆ July 22	◆	◆
PG2225	RIDA®GENE ETEC/EIEC	●	◆ Aug 22	◆ July 22	◆	◆
PG2285	RIDA®GENE E. coli Stool Panel I	●	●	◆ July 22	◆	◆
PG2305	RIDA®GENE Helicobacter pylori	◆	◆	◆	◆	◆
PG2405	RIDA®GENE Bacterial Stool Panel	●	●	◆ July 22	◆	◆
PG2415	RIDA®GENE Bacterial Stool Panel I	◆	—	—	—	—
Parasites						
PG1715	RIDA®GENE Parasitic Stool Panel I	●	—	—	—	—
PG1725	RIDA®GENE Parasitic Stool Panel II	●	◆ July 22	◆ June 22	◆	◆
PG1745	RIDA®GENE Dientamoeba fragilis	◆	◆	◆	◆	◆
Hospital acquired infections						
Viruses & bacteria						
PG0605	RIDA®GENE MRSA	●	◆	◆	◆	◆
PG0705	RIDA®GENE Hospital Stool Panel	●	◆	◆	◆	◆
PG0835	RIDA®GENE Clostridium difficile	●	●	◆ June 22	◆	◆
PG2705	RIDA®GENE CAP Bac	◆ July 22	◆ Aug 22	◆	◆	◆
Respiratory infections						
Viruses						
PG0505	RIDA®GENE Flu	●	●	◆	◆	◆
PG0545	RIDA®GENE Flu & RSV	●	◆	◆	◆	◆
PG5805	RIDA®GENE Parainfluenza	●	◆	◆	◆	◆
PG5905	RIDA®GENE RSV & hMPV	◆	◆	◆	●	◆
PG6815	RIDA®GENE SARS-CoV-2	●	●	●	●	◆
PG6825	RIDA®GENE Flu & SARS-CoV-2	●	●	●	●	◆
Bacteria						
PG2505	RIDA®GENE Bordetella	●	●	◆	◆	◆
Fungi						
PG1905	RIDA®GENE Pneumocystis jirovecii (qual)	◆	◆	◆	◆	◆
Human genetics						
PY0205	RIDA®GENE HLA-B27	◆ Aug 22	◆	◆	◆	◆
PY1205	RIDA®GENE Factor II	◆ Aug 22	◆	◆	◆	◆
PY1210	RIDA®GENE Factor V	◆ Aug 22	◆	◆	◆	◆
PY4215	RIDA®GENE Lac Intol	◆ Aug 22	◆	◆	◆	◆

● = available ◆ = in process