

RIDA®QUICK SCAN II

Easy and quick quantitative lateral flow measurements



03/2022



Simple operation:
Intuitive colored touchscreen



Imaging technique:
Wavelength independence and image documentation of test lines



Reliable measurement:
Internal calibration and control strip with certificate

More information:



<https://r-bio/15>

Our new LTF reader brings your measurements to a more precise and accurate level

Benefits



Efficient

Methods and test assay specifications via QR code technology



Unique

Developed for R-Biopharm's Clinical Diagnostics



Documentary advantage

Results as PDF, TXT and XML files



Practical

Firmware update with USB Stick



Compact

Small and portable instrument

Product overview

Product	Description	Art. No.
RIDA®QUICK SCAN II – IVD Set	Lateral flow reader and 2D barcode scanner	ZRQS2-KD-SET
RIDA®QUICK SCAN II	Lateral flow reader	ZRQS2-KD
Honeywell Xenon 1900	2D barcode scanner for ZRQS2-KD-SET	ZBS
RIDA®QUICK SCAN II control strip clinic	Control strip for ZRQS2-KD and ZRQS2-KD-SET	ZRQS2-KS-KD

Technical data

Size (LxBxH)	229 x 192 x 95 mm
Weight	1.470 g
Display	Resistive colored touchscreen
Measuring time	< 23 sec
Memory capacity	Up to 1400 measurement
Connectivity	USB, LAN
Inter-instrument CV	< +/- 3 %
Ambient conditions	
Temperature	5 - 40 °C

RIDA®QUICK SCAN II – IVD Set



- 1. Device
- 2. Control strip
- 3. Power adapter
- 4. Barcode scanner
- 5. Adapter
- 6. Certificate
- 7. User manual