


RIDA[®]SEEK

Your evaluation and documentation assistant

Benefits and functions - compact



Contents

What is RIDA [®] SEEK?	p. 02
Functions and benefits	p. 03
Get your result in 5 steps	p. 09
Analysis confirmation	p. 10
Your contact	p. 12

What is RIDA®SEEK?

RIDA®SEEK is an automated interpretation software which allows a qualitative interpretation of the results from raw laboratory data.

Benefits and functions of RIDA®SEEK:

- Automated interpretation of results
- Quality management made easy e.g. documentation of runs, assay lots and much more
- Communication with the laboratory information system (LIS)
- Operates with most common cyclers



Automated interpretation

Doubtful results are highlighted

- ▶ Proposal of result (positive/negative)
- ▶ Declaration of Ct-values
- ▶ Validity of controls



Your benefits:

Time savings through automated result interpretation as no detailed curve examination – especially for multiplex assays – is required.

Simplification of work processes by proposing a result.

Standardized interpretation of results because the algorithm was trained specifically for RIDA®GENE assays in combination with the cycler.

Your benefits:

Automatically documents all results in a comprehensive report. Thereby, saves time and avoids manual transmission errors.

Extensive quality management with an archive module and the possibility to trace runs, assay lots, controls etc. as well as a quality module to check the progress of controls.

Archive with information about runs, assay lots, controls and more

Quality management

Quality control module with progress of controls

To implement a connection of the software to your LIS please contact your LIS contact person.

Communicates with LIS



Your benefits:

Simplified data management and **error prevention** compared to manual data entry and transmission due to LIS-connection.

Time savings through standardized results reports.

Your benefits:

Time savings compared to manual documentation and **error prevention** during data entry and transmission. Work processes are simplified fundamentally.

Results are summarized in a report.

Documents everything

The scope of reports is configurable.

Intuitive user interface



Your benefits:

Easy training of new employees as well as **efficient use of diagnostic equipment** due to a clear user interface when using different cyclers.

Your benefits:

Data security – sensitive data will be stored on a local server, no cloud connection.

All data will be stored on a local server.

Data security



GET YOUR RESULT IN 5 STEPS.

Step 1: You select the test run in the software and check the plate set-up.



Step 2: You choose lot and (if required) color compensation.




Step 3: You confirm the analysis. If necessary you can check the curve and verify the decision.



Step 4: You approve the analysis.



Step 5: You receive a final report and export the results – directly into the laboratory information system!



Excerpt from the RIDA®SEEK software:
- Step 3: Confirmation of the analysis.

If necessary you can check the curve and verify the decision.

r-biopharm®



Questions? Please contact us!

- Automated interpretation of results
- Quality management made easy
e.g. documentation of runs, assay lots and much more
- Communication with the laboratory information system (LIS)
- Operates with most common cyclers



R-Biopharm

An der Neuen Bergstraße 17 · 64297 Darmstadt

www.r-biopharm.com · mdx@r-biopharm.de