

# Molecular diagnostics of gastrointestinal infections

Product portfolio



**Comprehensive:**  
GI diagnostics for all relevant pathogens



**Individual:**  
The suitable solution for every testing strategy



**Efficient:**  
Harmonized processing for virus, bacteria  
and parasite detection


More information:



<https://r-b.io/gastroen>

# Focus on GI diagnostics: Solve complex requirements efficiently

Gastrointestinal infections pose a significant global health threat. GI infections are the fourth leading cause of death worldwide in children under 5 years of age. The diagnosis of gastrointestinal infections poses a variety of challenges for laboratories, particularly in terms of time and resource efficiency.



## Workflow optimization using PCR



### Save time:

Results within hours instead of days enable quick and appropriate treatment decisions and support the fight against antibiotic resistance.



### Broad pathogen coverage:

A test can cover a wide range of possible pathogens.

a) Small multiplex assays: These assays test for a limited number of pathogens and are ideal for specific clinical questions.

b) High-multiplex assays: These assays can test for 20 or more pathogens simultaneously, which is particularly useful in unclear clinical situations where multiple pathogens must be considered.



### Reduced workload:

Fewer manual steps and handling errors.



### Precise diagnoses:

Higher sensitivity and specificity compared to traditional methods.



### Harmonized processing:

Test requirements for bacteria, viruses, and parasites with a single workflow.


# What are your laboratory's requirements?

## Fully automated processing and flexible combination of assays

### RIDA®UNITY

#### From extraction to sample analysis – all in one place

- Option of a syndromic approach with up to 6 assays
- Sample processing tailored to requirements
- 1-click evaluation of results
- Medium to high sample throughput
- IVDR certified system




## High-multiplex approach for comprehensive pathogen coverage

### RIDA®Xplore & RIDA®Jump

#### Detection using Multiplex Tandem PCR

- Efficient processing through parallel detection of viruses, bacteria, parasites; up to 20 pathogens in one run
- Semi-automated processing
- Can be combined with existing laboratory equipment
- Medium to high sample throughput




## Flexible processing of individual samples

### Vivalytic



#### Fully automated Point-of-Care system

- Cartridge: ready-to-use
- Easy handling
- Connection to LIS possible
- Cross-indication portfolio




## Different sample requirements in low to medium throughput

### RIDA®GENE

#### Direct detection using real-time PCR

- Efficient combination of all RIDA®GENE assays regardless of RNA or DNA detection
- Possibility of syndromic testing
- 20 thaw/freeze cycles
- Manual processing with existing laboratory equipment



# RIDA®UNITY

RIDA®UNITY assays are based on real-time PCR technology and are optimized for fully automated processing on the RIDA®UNITY system.

This system solution offers the possibility of processing 96 extractions and 192 reactions per run. Depending on the sample volume, two to three runs or up to 384 reactions can be processed in one working day (8 hours).



## Convincing due to

### Flexibility

- Up to six different assays can be combined in one run
- Samples can be processed according to specific requirements
- Modular processing
- Integration of external cyclers (BioRad)

### Efficiency

- Fully automated
- Reuse of reagents and eluates
- Matrix-independent extraction


### Safety

- Complete traceability
- Controlled loading
- Complete process controls
- IVDR certified system

### Speed

- Intuitive software handling
- Automatic evaluation software
- Ready-to-use reagents
- Time to result: 2.5 - 4.5 hours

## Simple, fully automated workflow



Products	Description	Art. No.
<b>Viruses</b>		
RIDA®UNITY Norovirus I & II	Norovirus GI Norovirus GII	UN1415
RIDA®UNITY Viral Stool Panel II	Rotavirus Adenovirus Astrovirus	UN1325
RIDA®UNITY Viral Stool Panel IV*	Norovirus Enterovirus Sapovirus	UN1345
RIDA®UNITY Hospital Stool Panel*	Norovirus Rotavirus <i>C. difficile</i>	UN0705
<b>Bacteria</b>		
RIDA®UNITY <i>C. difficile</i> **	<i>C. difficile</i> <i>C. difficile</i> toxin A/B	UN0835
RIDA®UNITY Bacterial Stool Panel	<i>Campylobacter</i> spp. <i>Salmonella</i> spp. <i>Yersinia enterocolitica</i>	UN2405
RIDA®UNITY EHEC/EPEC	EHEC EPEC STEC EIEC/ <i>Shigella</i> spp.	UN2205
RIDA®UNITY <i>E. coli</i> Stool Panel I*	Shiga toxin 1 Shiga toxin 2 EPEC	UN2285
RIDA®UNITY <i>E. coli</i> Stool Panel III*	ETEC EIEC/ <i>Shigella</i> spp. EAEC	UN2265
<b>Parasites</b>		
RIDA®UNITY Parasitic Stool Panel II	<i>Cryptosporidium</i> spp. <i>Entamoeba histolytica</i> <i>Giardia lamblia</i>	UN1725

\* Coming soon

\*\* Fully automated processing on the RIDA®UNITY, as well as manual processing with MagNAPure96 + LightCycler®480 II possible



**Research use only**


RIDA®UNITY Shiga Toxin 2 Typing RUO	<i>stx2a/c/d</i> <i>stx2f</i>	UN2245RUO
RIDA®UNITY Parasitic Stool Panel III RUO	<i>Cyclospora cayetanensis</i> <i>Blastocystis hominis</i> <i>Dientamoeba fragilis</i>	UN1755RUO



**Your advantages**

- > 98 % of GI pathogens covered by the RIDA®UNITY portfolio
- Reliable EHEC diagnostics with all relevant Shiga toxin subtypes (1a, 1c, 1d, 2a-g)
- Complete coverage of relevant *C. difficile* infection markers

**Seasonal pathogen distribution**



RV = Rotavirus, AdV = Adenovirus F40/41, AsV = Astrovirus, V = Vibrio, P = *Plesiomonas*, YE = *Yersinia enterocolitica*; \* Worldwide data

# RIDA®Xplore/ RIDA®Jump

The TandemPlex® panels are based on our patented multiplex tandem PCR (MT-PCR) technology. This is based on a two-step PCR, in which the gene fragments to be amplified are preamplified in a few copies before a further PCR with nested primers is performed. Processing is carried out on the RIDA®Xplore or RIDA®Jump devices for medium to high sample throughput.



## Convincing due to

### High-level multiplexing

- Syndromic detection of pathogen groups
- Up to 30 targets in one panel



### Simplicity

- Low hands-on time
- Ready-to-use reagents
- Intuitive evaluation software



### Safety

- Pathogen detection via real-time amplification curve and melting curve
- Low risk of contamination




### Modular design

- Semi-automatic
- Can be combined with existing laboratory equipment
- Optimal utilization due to separately operable systems



## Optimal capacity utilization – workflow



# RIDA®PrepX Stool Kit (Art. No. 90240)

## Rapid diagnosis of diarrhea samples

The RIDA®PrepX Stool Kit enables fast extraction to optimize the workflow of liquid stool samples. Without time-consuming DNA/RNA purification, samples can be processed directly and evaluated with high sensitivity using state-of-the-art PCR technology – for maximum efficiency and reliable results in everyday laboratory work.



Validated with	Panels
• RIDA®Xplore	• RIDA®Plex Faecal Pathogenes M 16-well
• RIDA®Jump	• RIDA®Plex Faecal Pathogenes A 16-well
	• RIDA®Plex Faecal Bacteria and Parasites 12-well
	• RIDA®Plex Enteric Viruses 8-Well

### Convincing due to

#### Speed


- 20-minute protocol
- Minimal manual effort

#### Reliability

- Suitable for bacteria, viruses, and parasites
- Efficient binding of PCR inhibitors

#### Cost savings

### Flexible workflow combinations



\* Conversion of plate

# RIDA®Xplore/RIDA®Jump Product portfolio

## Syndromic testing approach

### RIDA®Plex Faecal Pathogens M 16-well

Art. No. 25039



#### Bacteria

*Salmonella* spp.  
*Shigella* spp. and EIEC  
*Campylobacter jejuni*  
*Campylobacter coli*  
*Clostridioides difficile* toxin A & B  
Shiga toxin 1  
Shiga toxin 2  
*Yersinia enterocolitica*  
*Yersinia pseudotuberculosis*  
*Aeromonas* spp.

#### Parasites

*Giardia lamblia*  
*Cryptosporidium parvum*  
and *C. hominis*  
*Entamoeba histolytica*

#### Viruses

Norovirus genotype I  
Norovirus genotype II  
Rotavirus A  
Adenovirus types F and G  
Astrovirus  
Sapovirus

### RIDA®Plex Faecal Pathogens A 16-well

Art. No. 25031



#### Bacteria

*Salmonella* spp.  
*Shigella* spp. and EIEC  
*Campylobacter jejuni*  
*Campylobacter coli*  
*Clostridioides difficile* toxin A & B  
*Yersinia enterocolitica*  
*Yersinia pseudotuberculosis*  
*Aeromonas* spp.

#### Parasites

*Giardia lamblia*  
*Cryptosporidium parvum*  
and *C. hominis*  
*Entamoeba histolytica*  
*Blastocystis hominis*  
*Dientamoeba fragilis*

#### Viruses

Norovirus genotype I  
Norovirus genotype II  
Rotavirus A  
Adenovirus types F and G

#### Your advantages



- Universal solution for detecting viruses, bacteria, and parasites
- Covers 98 % of the most important GI parameters
- Option for extraction-free processing

#### Your advantages



- Universal solution for virus, bacteria, and parasite detection with a focus on comprehensive parasite diagnostics
- Option for extraction-free processing

## Specific detection of individual pathogen groups

### RIDA®Plex Enteric Viruses 8-well

Art. No. 85037



Norovirus genotype I  
Norovirus genotype II  
Rotavirus A  
Adenovirus types F and G

Astrovirus  
Sapovirus  
Enterovirus

### RIDA®Plex Parasites 8-well

Art. No. 85021



*Giardia lamblia*  
*Cryptosporidium parvum*  
and *C. hominis*  
*Entamoeba histolytica*

*Cyclospora cayentanensis* (RUO)  
*Blastocystis hominis* type 1  
*Blastocystis hominis* type 3  
*Dientamoeba fragilis*

### RIDA®Plex Faecal Bacteria and Parasites 12-well

Art. No. 85041



#### Bacteria

*Salmonella* spp.  
*Shigella* spp. and EIEC  
*Campylobacter jejuni*  
*Campylobacter coli*  
*Clostridioides difficile* toxin A & B  
*Escherichia coli* O157  
Shiga toxin 1  
Shiga toxin 2  
*Yersinia enterocolitica*  
*Yersinia pseudotuberculosis*

#### Parasites

*Giardia lamblia*  
*Cryptosporidium parvum*  
and *C. hominis*  
*Entamoeba histolytica*

#### Your advantages



- Targeted pathogen group-specific processing
- Option for extraction-free processing (exceptions: item no. 85021)

- More efficient processing: Detection of bacteria and parasites in a single workflow to replace two methods in the laboratory: cultivation and EIA/microscopy

# Vivalytic

Vivalytic enables simple and rapid PCR diagnostics with a cross-indication portfolio (gastrointestinal, respiratory, dermatological, CSF infections, STI, UTI, and HAI) both in the laboratory and at the Point of Care, where fast and reliable results are needed for therapeutic decisions, e.g., in hospitals, emergency rooms, or doctors' offices.



## Convincing due to

### Simplicity

- Fully automated processing
- Cartridge system
- 4 steps to results

### Speed



- Testing of individual samples
- Option for early finish

Product	Description	Art. No.
<b>Viren, Bakterien</b>		
Vivalytic Norovirus	Norovirus GI/GII	V1025
Vivalytic <i>C. difficile</i>	<i>C. difficile</i> Toxin A/B	V1019
Vivalytic Noro-, Rotavirus, <i>C. diff.</i>	Norovirus GI /GII Rotavirus Type A <i>C. difficile</i> Toxin A/B	V1026

### Your advantages

- Focus on the 3 most common GI pathogens (~ 50 %)
- Time to result < 1 hour

## Simple workflow



# RIDA®GENE

The RIDA®GENE kits are based on real-time PCR technology and represent a flexible manual solution that can be combined with various existing laboratory devices. The broad portfolio enables direct qualitative detection of gastrointestinal, respiratory, sexually transmitted, and hospital-acquired infections.



## Convincing due to

### Flexibility

- Combination of viruses, bacteria, and parasites in a single run
- Compatibility with various extraction and detection systems



### Requirement accuracy

- Single-source or multiplex approach possible
- Stable up to 20 thaw/freeze cycles
- Low to medium sample throughput




### User friendly

- Uniform test execution
- Ready-to-use reagents
- Internal control, positive control, negative control included
- Automatic evaluation with RIDA®SEEK



## Benefit from the combined processing of RIDA®GENE DNA and RNA assays



Product	Description	Art. No.
<b>Viruses</b>		
RIDA®GENE Viral Stool Panel I	Norovirus GI/GII Rotavirus Adenovirus Astrovirus	PG1315
RIDA®GENE Viral Stool Panel II	Rotavirus Adenovirus Astrovirus	PG1325
RIDA®GENE Viral Stool Panel III	Norovirus GI/GII Rotavirus Adenovirus	PG1335
RIDA®GENE Norovirus	Norovirus GI/GII	PG1405
RIDA®GENE Norovirus I & II	Norovirus GI Norovirus GII	PG1415
RIDA®GENE Sapovirus	Sapovirus	PG1605
RIDA®GENE Enterovirus	Enterovirus	PG4705
RIDA®GENE Hospital Stool Panel	Norovirus GI/GII Rotavirus <i>C. difficile</i> Toxin A/B	PG0705
<b>Bacteria</b>		
RIDA®GENE EHEC/EPEC	EHEC EPEC STEC EIEC/ <i>Shigella</i> spp.	PG2205
RIDA®GENE EAEC	EAEC	PG2215
RIDA®GENE ETEC/EIEC	ETEC EIEC/ <i>Shigella</i> spp.	PG2225
RIDA®GENE E. coli Stool Panel I	Shiga Toxin 1 Shiga Toxin 2 EPEC	PG2285
RIDA®GENE Bacterial Stool Panel	<i>Campylobacter</i> spp. <i>Salmonella</i> spp. <i>Yersinia enterocolitica</i>	PG2405
RIDA®GENE Bacterial Stool Panel I	<i>Campylobacter</i> spp. <i>Salmonella</i> spp. STEC EIEC/ <i>Shigella</i> spp.	PG2415
RIDA®GENE Clostridium difficile	<i>C. difficile</i> <i>C. difficile</i> Toxin A/B	PG0835
<b>Parasites</b>		
RIDA®GENE Parasitic Stool Panel I	<i>Cryptosporidium</i> spp. <i>Entamoeba histolytica</i> <i>Giardia lamblia</i> <i>Dientamoeba fragilis</i>	PG1715
RIDA®GENE Parasitic Stool Panel II	<i>Cryptosporidium</i> spp. <i>Entamoeba histolytica</i> <i>Giardia lamblia</i>	PG1725
<b>Research use only</b>		
RIDA®GENE Parasitic Stool Panel IV RUO	<i>Giardia intestinalis</i> <i>Entamoeba histolytica</i> <i>Cryptosporidium</i> spp. <i>Cyclospora cayetanensis</i>	PG1765RUO
RIDA®GENE Bacterial Stool Panel III RUO	<i>Yersinia enterocolitica</i> <i>Aeromonas</i> spp. <i>Vibrio cholerae</i> <i>Plesiomonas shigelloides</i>	PG2435RUO
RIDA®GENE Bacterial Stool Panel V RUO	<i>Campylobacter</i> spp. <i>Salmonella</i> spp. <i>Yersinia enterocolitica</i> EIEC/ <i>Shigella</i> spp.	PG2455RUO



## Advantages for your *E. coli* diagnostics



- Comprehensive PCR portfolio of human pathogens *E. coli*
- Reliable EHEC diagnostics
- Differentiation between Shiga toxin 1 and 2
- Detection of all relevant Shiga toxin subtypes (1a, 1c, 1d, 2a - g)
- Validated for stool samples and culture

	<i>stx1</i>	<i>stx2</i>	<i>eae</i>	<i>ipaH</i>	<i>aat/aggR</i>	LT	ST
PG2205	STEC						
			EHEC				
			EPEC				
				EIEC/ <i>Shigella</i> spp.			
PG2215					EAEC		
PG2225				EIEC/ <i>Shigella</i> spp.			ETEC
PG2285		STEC					
			EHEC				
			EPEC				

# Overview of pathogens

	RIDA®UNITY	RIDA®GENE	RIDA®Xplore/RIDA®Jump	Vivalytic
<b>Viruses</b>				
Adenovirus	●	●	●	
Astrovirus	●	●	●	
Enterovirus	*	●	●	
Norovirus	●	●	●	●
Rotavirus	●	●	●	●
Sapovirus	*	●	●	
<b>Bacteria</b>				
<i>Aeromonas</i> spp.			●	
<i>Campylobacter</i> spp.	●	●	●	
<i>C. difficile</i>	●	●	●	●
<i>E. coli</i> O175			●	
EHEC	●	●	RUO	
EPEC	●	●	RUO	
ETEC	*	●		
EIEC	●	●	●	
EAEC	*	●		
<i>Helicobacter pylori</i>	(RUO)	●		
<i>Salmonella</i> spp.	●	●	●	
<i>Shigella</i> spp.	●	●	●	
STEC	●	●	●	
<i>Yersinia enterocolitica</i>	●	●	●	
<i>Yersinia pseudotuberculosis</i>	*		●	
<b>Parasites</b>				
<i>Blastocystis hominis</i>	RUO	●	●	
<i>Cryptosporidium</i> spp.	●	●	●	
<i>Dientamoeba fragilis</i>	RUO	●	●	
<i>Cyclospora cayetanensis</i>	RUO		RUO	
<i>Entamoeba histolytica</i>	●	●	●	
<i>Giardia lamblia</i>	●	●	●	

\* In development  
(RUO) = Research Use Only, for manual processing only