

Diagnostics of enteropathogenic infections by RIDASCREEN® EasyKit ELISA

Convenient – fast – reliable

- Ready-to-use reagents
- Minimum time expenditure thanks to simple preanalytics
- Reagent vials optimized for direct use in fully automated instruments




RIDASCREEN® EasyKit ELISA sets new standards for economical and flexible diagnostics

- Greater precision when manual processed
- Easy analysis using software – RIDA®SOFT Win.net
- Perfectly adapted for use in fully automated instruments (DSX®, Agility®)

Dynex' Agility® fully automated instrument avoids cross-contamination by use of disposable tips. Both instruments can also be connected quickly and easily to the laboratory information system (LIMS).

Dynex Agility®



EasyKit reagent sets can be used directly in the device-specific mounts of the fully automated machine Agility®. A 2D barcode with lot-specific test information ensures safe management of the reagents.

During test processing, the lids of the reagents are kept in corresponding stackable trays to avoid mix-ups.





Proven simplified preanalytics using universal sample diluent and uniform sample dilution for all parameters.

Dynex DSX®



A rack specially adapted to the DSX® and the EasyKit reagents with an attached tray for the reagent lids ensures safe and contamination-free processing.

- 30 spaces for reagents
- 24 spaces for controls
- 2 spaces for sample buffers

Product Overview – RIDASCREEN® EasyKit ELISA

ELISA

Product	Art. No.
Virus	
RIDASCREEN® Adenovirus	C1001
RIDASCREEN® Astrovirus	C1301
RIDASCREEN® Norovirus	C1401
RIDASCREEN® Rotavirus	C0901
Bacteria & Bacteria toxins	
RIDASCREEN® Campylobacter	C2401
RIDASCREEN® Clostridium difficile GDH	C0701
RIDASCREEN® Clostridium difficile Toxin A/B	C0801
RIDASCREEN® Clostridium perfringens Enterotoxin	C0601
RIDASCREEN® Helicobacter pylori	C2301
RIDASCREEN® Verotoxin	C2201
Parasites	
RIDASCREEN® Cryptosporidium	C1201
RIDASCREEN® Cryptosporidium/Giardia Combi	C1121
RIDASCREEN® Entamoeba	C1701
RIDASCREEN® Giardia	C1101

Instruments

DSX®	Z65400
Agility®	ZAGILITY

Controls for RIDASCREEN®*

Product	Art. No.
Virus	
Adenovirus	CRP1004
Astrovirus	CRP1304
Norovirus	CRP1404
Rotavirus	CRP0904
Bacteria	
Campylobacter	CRP2404
Clostridium difficile GDH	CRP0704
Helicobacter pylori	CRP2304
Bacteria toxins	
Clostridium difficile Toxin A/B	CRP0804
Clostridium perfringens Enterotoxin	CRP0604
Verotoxin	CRP2204
Parasites	
Cryptosporidium	CRP1204
Cryptosporidium/Giardia Combi	CRP1124
Entamoeba	CRP1704
Giardia	CRP1104

* All controls are supplied as a set of 1.8 ml positive controls and 1.8 negative controls.

R-Biopharm Contacts:

Product Manager

Infectious Diseases/Antigen Detection:

Helmut Leidinger
 Phone: +49 (0) 61 51 - 81 02-26
 Fax: +49 (0) 61 51 - 81 02-40
 E-mail: h.leidinger@r-biopharm.de

Clinical Sales International:

Phone: +49 (0) 61 51 - 81 02-0
 Fax: +49 (0) 61 51 - 81 02-688
 E-mail: clinical.sales@r-biopharm.de