



Your competent partner in gastroenterology

R-Biopharm AG

ELISA tests with diagnostic focus on:

Colon cancer prevention

RIDASCREEN® Haemoglobin
RIDASCREEN® Haemo-/Haptoglobin Complex

Inflammatory bowel disease

RIDASCREEN® Calprotectin
RIDASCREEN® α_1 -Antitrypsin
RIDASCREEN® sIgA

Therapeutic drug monitoring

RIDASCREEN® IFX Monitoring
RIDASCREEN® Anti-IFX Antibodies
RIDASCREEN® ADM Monitoring
RIDASCREEN® Anti-ADM Antibodies
RIDASCREEN® VDZ Monitoring
RIDA®QUICK IFX Monitoring

Pancreatic diagnostics

Pancreatic Elastase ELISA
Pankrin® ELISA

Haemoglobin and Haemo-/Haptoglobin Complex
Markers for colon cancer screening. A combined testing of both parameters increases the probability to detect adenomas.

Calprotectin

Indicator of gastrointestinal inflammation, used for differentiation of IBD and IBS and for treatment monitoring.

α_1 -Antitrypsin and sIgA

Indicators for the integrity of the gut wall, relevant for monitoring and diagnosing a Leaky Gut Syndrome.

Therapeutic Drug Monitoring (TDM)

TDM of Infliximab (IFX) and its biosimilars as well as Adalimumab (ADM) and Vedolizumab (VDZ) supports the therapy optimization and allows a treatment cost reduction at the same time.

Pancreatic elastase and pankrin

Markers for the detection of pancreas insufficiency and pancreatitis.

Accessories:

Stool collection tubes

Easy, rapid and clean collection of stool samples.

R-Biopharms Gastroenterology product portfolio is dedicated to the advanced diagnosis of gut-related diseases. The ELISA tests are suited for single plate or medium-high-throughput analysis on automated readers.

ELISA

Product	Description	Matrix	Tests	Art. No.
Enzyme immunoassays				
RIDASCREEN® Haemoglobin	Enzyme immunoassay for the detection of haemoglobin	Stool	96	G09030
RIDASCREEN® Haemo-/Haptoglobin Complex	Enzyme immunoassay for the detection of haemoglobin/haptoglobin complex	Stool	96	G09031
RIDASCREEN® α ₁ -Antitrypsin	Enzyme immunoassay for the detection of α ₁ -Antitrypsin	Stool	96	G09034
RIDASCREEN® sIgA	Enzyme immunoassay for the detection of secretory IgA	Stool	96	G09035
RIDASCREEN® Calprotectin	Enzyme immunoassay for the detection of calprotectin	Stool	96	G09036
Pancreatic Elastase ELISA	Enzyme immunoassay for the detection of pancreatic elastase	Stool	96	G09038
Pancreatic Elastase ELISA (SK15)	Enzyme immunoassay for the detection of pancreatic elastase; additional standard (SK15)	Stool	96	G09040
Pankrin® ELISA	Enzyme immunoassay for the detection of pancreatic elastase and other pancreatic enzymes	Serum	96	G09039

Therapeutic Drug Monitoring (TDM)

Enzyme immunoassays				
RIDASCREEN® IFX Monitoring	Enzyme immunoassay for the detection of infliximab and its biosimilars	Serum/plasma	96	G09041
RIDASCREEN® Anti-IFX Antibodies	Enzyme immunoassay for the detection of antibodies to infliximab and its biosimilars	Serum/plasma	96	G09042
RIDASCREEN® ADM Monitoring	Enzyme immunoassay for the detection of adalimumab	Serum/plasma	96	G09043
RIDASCREEN® Anti-ADM Antibodies	Enzyme immunoassay for the detection of antibodies to adalimumab	Serum/plasma	96	G09044
New RIDASCREEN® VDZ Monitoring	Enzyme immunoassay for the detection of vedolizumab	Serum/plasma	96	G09045
Rapid test				
New RIDA®QUICK IFX Monitoring	Immunochromatographic lateral flow assay for the quantification of infliximab (IFX, Remicade®) and its biosimilars Remsima® and Inflectra®	Serum/plasma	25	GN3041

Accessories

Accessories for RIDASCREEN®				
New RIDA®TUBE	For collection and preparation of stool samples, • Unfilled; to use after internal validation		50	GZ3013
RIDA®TUBE Haemoglobin	For collection and preparation of stool samples, • Only use with RIDASCREEN® Haemoglobin G09030 and RIDASCREEN® Haemo-/Haptoglobin Complex G09031		50	GZ3012
New RIDA®TUBE Calprotectin	For collection and preparation of stool samples, • Only use with RIDASCREEN® Calprotectin G09036		50	GZ3016
Stool Preparation Set	For collection and preparation of stool samples, • Only use with Pancreatic Elastase ELISA G09038 and G09040		45	GZ3008
RIDA®QUICK SCAN II	Lateral flow reader – for read out of GN3041			ZRQS2-KD
New RIDA®QUICK IFX Monitoring Control Set	Positive controls – accessory for GN3041			GP3041

R-Biopharm contacts:

Clinical Sales International:
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-688
E-mail: clinical.sales@r-biopharm.de

Order Department:
Phone: +49 (0) 61 51 - 81 02-0
Fax: +49 (0) 61 51 - 81 02-20
E-mail: orders@r-biopharm.de