



## RIDA®SEEK software – product and platform compatability

		Real-time PCR instruments			
		Roche LightCycler® 480 II	ABI 7500 Fast	Bio-Rad CFX96™	Agilent Techn. Mx3005P
<b>Gastrointestinal infections</b>					
<b>Virus</b>					
PG1315	RIDA®GENE Viral Stool Panel I	•			•
PG1325	RIDA®GENE Viral Stool Panel II	•	•	•	•
PG1405	RIDA®GENE Norovirus	•			•
PG1415	RIDA®GENE Norovirus I & II	•	•		•
PG1605	RIDA®GENE Sapovirus	•			
<b>Bacteria</b>					
PG2205	RIDA®GENE EHEC/EPEC	•	•		•
PG2215	RIDA®GENE EAEC	•			
PG2225	RIDA®GENE ETEC/EIEC	•			
PG2285	RIDA®GENE E. coli Stool Panel I	•			•
PG2405	RIDA®GENE Bacterial Stool Panel	•	•	•	•
<b>Hospital acquired infections</b>					
<b>Virus &amp; Bacteria</b>					
PG0605	RIDA®GENE MRSA	•	•	•	•
PG0705	RIDA®GENE Hospital Stool Panel	•			•
PG0825	RIDA®GENE CD Toxin A/B	•			•
PG0835	RIDA®GENE Clostridium difficile	•	•		•
<b>Respiratory infections</b>					
<b>Virus</b>					
PG0505	RIDA®GENE Flu	•	•	•	•
PG0545	RIDA®GENE Flu & RSV	•		•	•
PG5805	RIDA®GENE Parainfluenza	•	•	•	
<b>Bacteria</b>					
PG2505	RIDA®GENE Bordetella	•		•	•
PG4305	RIDA®GENE Mycoplasma pneumoniae	•			
PG8005	RIDA®GENE Legionella	•			
<b>Sexually transmitted infections</b>					
<b>Bacteria</b>					
PG4945	RIDA®GENE STI Mycoplasma Panel	•			
<b>Parasites</b>					
PG4975	RIDA®GENE Trichomonas vaginalis	•			