



Your competent partner in gastroenterology

R-Biopharm AG

ELISA tests with diagnostic focus on:

Colon cancer prevention

RIDASCREEN® Haemoglobin
RIDASCREEN® Haemo-/Haptoglobin Complex

Inflammatory bowel disease

RIDASCREEN® Calprotectin
RIDASCREEN® α_1 -Antitrypsin
RIDASCREEN® sIgA

Pancreatic diagnostics

Pancreatic Elastase ELISA
Pankrin® ELISA

Therapeutic drug monitoring

RIDASCREEN® IFX Monitoring
RIDASCREEN® Anti-IFX Antibodies
RIDASCREEN® ADM Monitoring
RIDASCREEN® Anti-ADM Antibodies
RIDASCREEN® VDZ Monitoring
RIDASCREEN® GLM Monitoring
RIDA®QUICK IFX Monitoring
RIDA®QUICK ADM Monitoring

Haemoglobin and Haemo-/Haptoglobin Complex
Markers for colon cancer screening. A combined testing of both parameters increases the probability to detect adenomas.

Calprotectin

Indicator of gastrointestinal inflammation, used for differentiation of IBD and IBS and for treatment monitoring.

α_1 -Antitrypsin and sIgA

Indicators for the integrity of the gut wall, relevant for monitoring and diagnosing a Leaky Gut Syndrome.

Therapeutic Drug Monitoring (TDM)

TDM of Infliximab (IFX) and its biosimilars as well as Adalimumab (ADM), Vedolizumab (VDZ) and Golimumab (GLM) supports the therapy optimization and allows a treatment cost reduction at the same time.

Pancreatic elastase and pankrin

Markers for the detection of pancreas insufficiency and pancreatitis.

Accessories:

Stool collection tubes

Allow an easy, rapid and clean collection of stool samples.

R-Biopharm's Gastroenterology product portfolio is dedicated to the advanced diagnosis of gut-related diseases. The ELISAs are suited for single plate or medium-high-throughput analysis on automated ELISA systems.

ELISA

Product	Description	Matrix	Tests	Art. No.
Enzyme immunoassays				
RIDASCREEN® Haemoglobin	Enzyme immunoassay for the quantitative detection of haemoglobin	Stool	96	G09030
RIDASCREEN® Haemo-/Haptoglobin Complex	Enzyme immunoassay for the quantitative detection of haemoglobin/haptoglobin complex	Stool	96	G09031
RIDASCREEN® α ₁ -Antitrypsin	Enzyme immunoassay for the quantitative detection of α ₁ -Antitrypsin	Stool	96	G09034
RIDASCREEN® sIgA	Enzyme immunoassay for the quantitative detection of secretory IgA	Stool	96	G09035
RIDASCREEN® Calprotectin	Enzyme immunoassay for the quantitative detection of calprotectin	Stool	96	G09036
Pancreatic Elastase ELISA	Enzyme immunoassay for the quantitative detection of pancreatic elastase	Stool	96	G09038
Pancreatic Elastase ELISA (SK15)	Enzyme immunoassay for the quantitative detection of pancreatic elastase; additional standard (SK15)	Stool	96	G09040
Pankrin® ELISA	Enzyme immunoassay for the quantitative detection of pancreatic elastase and other pancreatic enzymes	Serum	96	G09039

Therapeutic Drug Monitoring (TDM)

Enzyme immunoassays				
RIDASCREEN® IFX Monitoring	Enzyme immunoassay for the quantitative detection of infliximab (IFX) and its biosimilars	Serum/plasma	96	G09041
RIDASCREEN® Anti-IFX Antibodies	Enzyme immunoassay for the quantitative detection of antibodies to infliximab (IFX) and its biosimilars	Serum/plasma	96	G09042
RIDASCREEN® ADM Monitoring	Enzyme immunoassay for the quantitative detection of adalimumab (ADM)	Serum/plasma	96	G09043
RIDASCREEN® Anti-ADM Antibodies	Enzyme immunoassay for the quantitative detection of antibodies to adalimumab (ADM)	Serum/plasma	96	G09044
RIDASCREEN® VDZ Monitoring	Enzyme immunoassay for the quantitative detection of vedolizumab (VDZ)	Serum/plasma	96	G09045
New RIDASCREEN® GLM Monitoring	Enzyme immunoassay for the quantitative detection of golimumab (GLM)	Serum/plasma	96	G09047
Rapid test				
RIDA®QUICK IFX Monitoring	Immunochromatographic lateral flow assay for the quantitative determination of infliximab (IFX) and its biosimilars	Serum/plasma	25	GN3041
New RIDA®QUICK ADM Monitoring	Immunochromatographic lateral flow assay for the quantitative determination of adalimumab (ADM)	Serum/plasma	25	GN3043

Accessories

Accessories for RIDASCREEN®				
RIDA®TUBE	For collection and preparation of stool samples, <ul style="list-style-type: none"> • Unfilled; to use after internal validation 		50	GZ3013
RIDA®TUBE Haemoglobin	For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with RIDASCREEN® Haemoglobin G09030 and RIDASCREEN® Haemo-/Haptoglobin Complex G09031 		50	GZ3012
RIDA®TUBE Calprotectin	For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with RIDASCREEN® Calprotectin G09036 		50	GZ3016
Stool Preparation Set	For collection and preparation of stool samples, <ul style="list-style-type: none"> • Only use with Pancreatic Elastase ELISA G09038 and G09040 		45	GZ3008
RIDA®QUICK SCAN II - IVD SET	Lateral flow reader (CE-IVD,) and 2D barcode scanner <ul style="list-style-type: none"> • for read out of GN3041 and GN3043 			ZRQS2-KD-SET
RIDA®QUICK IFX Monitoring Control Set	Positive controls <ul style="list-style-type: none"> • accessory for GN3041 			GP3041
New RIDA®QUICK ADM Monitoring Control Set	Positive controls <ul style="list-style-type: none"> • accessory for GN3043 			GP3043